aHealth smart Watch

aWatch

INDEX

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aHealth smart Watch

PRODUCT DESCRIPTION

aHealth watch uses real-time photoelectric detection technology for measuring blood pressure, electrocardiogram, heart rate, fatigue, emotional level, daily activities level, sleep quality and other healthy related physiological data. Once aHealth watch found the user's physiological data is abnormal, user can press and hold SOS emergency button for promptly notifying families or doctor in case of emergency. aHealth watch can minimize and reduce risk of serious health problems, it is an advantageous health pre-caution device for elders who have blood pressure or heart rate issues.

aHealth watch has GPS location tracking feature, and it will be very easy and convenience to find the watch location by using the aHealth application which is free. App can be download from both stores iOS and Google Play. You may also view the watch GPS tracking history, it will be very helpful for family members to track the previous route of the members with Alzheimer's disease or amnesia.

aHealth watch can measure health related physiological data at any time or anywhere. And you may read testing result immediately on screen of the watch.


@Health smart Watch

The data can be upload to the cloud server for tracking of your health 24/7. Your health data in the cloud server will be helpful for sending reports to your doctor or professional health care providers. Users will have more precisely and objectively health care maintenance with big data analyzing.

aHealth watch is a wearable health data measuring device that provides timely access to you and your families an advantage health protection plan. A very simple device for those patients who experience high risk health issues.


@Health smart Watch

1-Power on/off 2-Electrodes 3-SOS button

@Health smart Watch


Note: Both electrodes must be held simultaneously to measure blood pressure.

4-Photoelectric sensor 5-Magnetic charging port

@Health smart Watch

CHARGING AND TURNING THE AHEALTH WATCH ON/OFF

Connect the watch to the charger as shown in illustration. When the watch has been successfully charged, the charger will glow green. Press and hold the power button to turn the watch on.

APPLICATION INSTALLATION AND SYNCING THE AHEALTH WATCH

1. Installation: On your mobile device, the App Store search for the “aHealth” app, or use the built in QR code on the watch to download and install the application.

2. Syncing: Register and login to the aHealth application, then choose “Device Management” . Sync your aHealth watch to the app by scanning the QR code on the watch.


Note: Due to the watch build-in GSM module may not support some countries, when using without 2G network ability (Japan, Taiwan), the watch and application can be synced with a Bluetooth connection.

3. Bluetooth connection: switch on your mobile phone Bluetooth, and search the watch “aHealth watch” for pairing.

@Health smart Watch

POWER ON

Press and hold top right button to turn watch on.

STANDBY/SLEEP MODE

Press right button. To access the menu, swipe left.

MOVEMENT


Calorie output is calculated using the total number of steps, height and weight. Access the steps menu to set your height, weight and targeted goals.

HEART RATE


Click heart rate icon to begin measuring. Watch must be worn tightly on the wrist for best results.

@Health smart Watch

BLOOD PRESSURE

Click blood pressure icon and hold both electrodes as shown on illustration. Keep the electrodes contacted with skin firmly and relax both arm for best result. The watch should be worn on the left wrist

@Health smart Watch


NOTE: BLOOD PRESSURE MEASURING POSTURE

(1) Watch should be worn above left hand wrist joint, and watch right side close to ulna styloid process.

(2) Watch should be fastened around wrist, do not leave gaps, but no effect to normal circulation of the blood.

(3) Relax the left hand, put the left forearm on a desk, keep right thumb and forefinger inside moist, then please put right thumb and forefinger inside for contacting with two metal electrodes on the watch side, making, the right ring finger and little finger gently to the left wrist, and put right hand wrist and forearm on the desktop and relax.


(4) Relief your shoulder, maintain a comfortable sitting, then start blood pressure measuring, please keep quiet and don't move during whole measurement process, even breathing, keep the gesture as shown on illustration

@Health smart Watch

ELECTRO CARDIOGRAM

Electrodes with fingers as shown on illustration. The screen will start counting down 120 - 0, the ECG report will display on screen. Select OK/SAVE to upload to cloud server.

Please follow below illustration indicated posture to wear aHealth watch for measuring the ECG (if there has light leakage between watch and wrist skin, the measurement results will be inaccurate).

@Health smart Watch


SLEEP MONITORING

Sleep patterns are created based on the time shown on watch. You may also manually input your own sleep schedules.

SOS EMERGENCY CALL

Press and hold SOS button for 3 seconds to send alert to your family members and your doctor for getting immediately help in case of emergency.

@Health smart Watch

GPS TRACKING

Switch on GPS mode on the watch and start tracking your location on the smart phone application.

WEATHER


The weather menu displays the current date, time as well as current and future weather conditions.

@Health smart Watch

CALENDAR

View current date and calendar viewing.

ALARM CLOCK


Click the alarm clock icon to enter your alarm information. When an alarm goes off, the watch vibrating and ringing to notify you.

@Health smart Watch

SET UP

(1) GPS Positioning: You may turn off GPS location settings, however, if you turn the GPS off you cannot be located in emergencies;

(2) Bluetooth: Open Bluetooth settings to pair the aHealth Watch with a sound source.

(3) Preference Settings:

Screen Brightness: The aHealth Watch will automatically adjust to surrounding light conditions.

Date and Time: Enter the menu to set the date and time as well as change the watch display.

Language Settings: Enter the menu to switch between languages.

@Health smart Watch

Clock: In this menu, you can customize the look of the clock when in standby mode, and manage clock settings.

Display and Brightness Adjustment: The display automatically adjusts, but can be adjusted by hand.

(4) **Syncing the Watch:** Sync the aHealth watch to the aHealth application using the two-dimensional code in this menu.

(5) **Restore Factory Settings:** Click this menu to reset your aHealth Watch and erase all stored data.

(6) **A:** Display device information.

SUPPORT AND AFTER SALE SERVICE

Still have questions? We've got answers and we're standing by. Give us a shout, we're here to help.

FCC Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

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