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support@trihear.com  
https://trihear.com  
Trihear Technology Co., Ltd.  
507, Building 4, Pingshan Private Enterprise  
Science and Technology Park, No.65 Lishan Road,  
Taoyuan Street, Nanshan District, Shenzhen, China

FCC ID: 2A26H-THN18M

Trihear SmartMic

User Manual  
Model Name: TH-N18M

**Preface**

Thank you for purchasing your new Trihear SmartMic. SmartMic is a Bluetooth hearing-assist microphone designed to help people who use hearing aids hear speech more clearly in everyday situations, such as conversations, TV audio, or phone audio when connected via a 3.5mm cable. This user manual will guide you through everything you need to know to get started with your SmartMic. If you have any problems setting up or using your Trihear ClearTV, please contact us by email at: support@trihear.com or scan the QR code to log in our website: [trihear.com](https://trihear.com)

**Package Contents**

Smart Mic, Type-C charging cable, 3.5mm audio cable, User Manual

**Product Overview**

3.5mm audio port, Pairing button / Volume down, Power button / Volume up, Indicator light, Smart Mic (built-in mic), Type-C charging port, Back clip

1

**Product Specifications**

|                  |                                        |
|------------------|----------------------------------------|
| Frequency range  | 2402-2480GHz                           |
| Connection range | 10m (33ft)                             |
| Battery life     | 8 hours (typical scenario: within 3 m) |
| Battery type     | Lithium-ion battery                    |
| Charging time    | 2 hours                                |
| Dimensions       | 45*26*18mm                             |
| Net weight       | 18g                                    |

Up to 10m (33ft)

There is no line-of-sight requirement between the hearing aid and Smart Mic. However, environmental interference such as walls or furniture may reduce the connection distance.

2

**How to Use**

**Charging**

Charge the device using the included Type-C charging cable. While charging, the green indicator flashes; when fully charged, the green indicator stays on. Charging time is about 2 hours. (Power adapter not included. Please use a standard 5V/1A power adapter. We recommend fully charging the device before first use.)

**Power On / Off**

(1) Power on: Press and hold " " for 2 seconds. The device turns on and the green indicator lights up.

(2) Power off: Press and hold " " for 2 seconds. The device turns off; the green indicator flashes quickly and then goes out.

3

Indicator light

Power off:

Press and hold 2 seconds

Smart Mic can be used normally while charging. To extend service life, it's recommended to turn the device off when it is not in use.

**Connection Methods**

(1) Use the 3.5 mm audio cable to connect to a set-top box or mobile phone.

4

(2) Clip the device to the other person's collar or chest.

5

**2. Reconnection**

If a pairing record exists, turn on Smart Mic, then turn on the hearing aid within 10 minutes; it will automatically connect.

If Smart Mic does not connect to any hearing aid or other device within 10 minutes, it will automatically power off.

**3. Force Pairing**

Press and hold the Pairing button for 3 seconds. The blue indicator flashes rapidly, indicating pairing mode.

**4. Clear Bluetooth Pairing Records**

With Smart Mic powered on, press and hold the Pairing button for 10 seconds. The blue indicator stays on for 2 seconds and then flashes rapidly, indicating the pairing record has been cleared and the device has re-entered pairing mode.

**Volume Adjustment**

After Smart Mic is connected to the hearing aid, single-press the Power button to increase the volume, and single-press the Pairing button to decrease the volume.

6

**Troubleshooting**

**1) Cannot connect to the hearing aid**  
Cause: Smart Mic is not in pairing mode; the hearing aid is not powered on.  
Solution: Press and hold the Pairing button for 3 seconds to enter pairing mode; turn on the hearing aid.

**2) No sound from 3.5mm audio input**  
Cause: The 3.5mm audio cable is not securely connected.  
Solution: Unplug and replug the 3.5mm audio cable, and make sure it is fully inserted.

**3) Low volume**  
Cause: Volume set too low / mic opening blocked.  
Solution: Increase volume using the buttons; clean the mic opening to ensure it is unobstructed.

**4) Cannot charge**  
Cause: A non-included charging cable is used.  
Solution: Use the included charging cable and connect it as described in this manual.

**5) Buttons unresponsive**  
Cause: Device has frozen.  
Solution: Press and hold the Power button while plugging in the charging cable for 2 seconds, then unplug the cable to recover.

7

**6) No operation / cannot power on**  
Cause: Battery is depleted.  
Solution: Charge the device for 10 minutes and try again.

**7) Disconnects after connection**  
Cause: Out of range, or walls/obstacles causing interference.  
Solution: Keep the distance within 10m; remove obstacles between the devices and avoid signal interference.

**CE Statement**

Use the Trihear SmartMic in the environment with the temperature between -10°C and 40°C, otherwise, it may damage your Product Name. It can be operating under 2000m.

Hereby, (Trihear Technology Co., Ltd.) declare that the radio equipment type Trihear SmartMic is compliance with Directive 2014/53/EU.

The full text of the EU declaration of conformity is available at the following internet address: [www.trihear.com](https://www.trihear.com)

CE

8

**FCC Statement**

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:  
(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC RF Exposure Statement:  
This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment can be used in portable exposure conditions.

9

**Warranty Policy**

Warranty Policy  
Purchasing Date:  
Model No:  
Serial no:  
Instruction: Thank you for choosing Trihear products, please note that for the main units, we offer 1 year warranty time from your purchasing date. For other accessories in the product package, 6 months Warranty time be offered. Any man-made damaged parts are not covered by the warranty service.  
When you are asking for warranty service, please take the product together with the purchase invoice or receipt with you, and get after sell service from where you buy or local designated after-sales service center. Any more question or demands you can contact us directly by email: [support@trihear.com](mailto:support@trihear.com)

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