





Size:248*159*64mm
250克铜版纸（双铜） 1200灰板

产品使用说明书

超声波洁牙仪

·感谢您购买我们的产品。

·使用前请务必仔细阅读此说明书并妥善保管，以备将来参考。

·本产品仅用于家庭个人牙齿清洁，不可用于任何医疗用途。

产品介绍

感谢选择本产品作为您的口腔健康助理。请按下图所示产品清单检查相关部件是否齐备。

平头洁牙头 X2

尖头洁牙头 X2

充电器

底座

刮舌器

探针

口镜

小扳手

保护罩

充电器底座

刮舌器

探针

口镜

小扳手

产品使用

警告

开启洁牙仪时,切勿将金属清洁头对准眼睛等人体敏感部位以及电源、灯具等物。产品开启时,请勿用皮肤触碰金属清洁头以免受伤。

注意:首次使用前请将产品充满电,因为出厂时为了安全考虑,主机并未充电。充电时,产品不能使用。

拆装清洁头

1. 安装清洁头

将清洁头对准底座上的卡槽,将清洁头插入到底,直到听到“咔哒”声为止。

2. 拆卸清洁头

握住底座上的卡扣,将清洁头向上提起,直到听到“咔哒”声为止。

切换模式

长按底座上的“模式”按钮,直到听到“滴滴”声为止。

短按底座上的“模式”按钮,直到听到“滴滴”声为止。

长按底座上的“模式”按钮,直到听到“滴滴”声为止。

清洁牙齿

1. 清洁牙齿时,将清洁头对准牙齿,轻轻移动清洁头,直到听到“滴滴”声为止。

2. 清洁牙齿时,将清洁头对准牙齿,轻轻移动清洁头,直到听到“滴滴”声为止。

3. 清洁牙齿时,将清洁头对准牙齿,轻轻移动清洁头,直到听到“滴滴”声为止。

产品充电

将底座充电器侧面的USB接口通过USB充电线连接至电脑或充电器。当底座充电器上的指示灯亮起时,即可开始充电。

充电时,请勿将产品放入水中或任何液体中。

充满电后,充电器指示灯将熄灭。

产品参数

产品名称: 超声波洁牙仪

额定功率: 3W

防水等级: IPX6

产品规格: 3档调节

额定电压: 5V

无线充电电流: 1A

电容量: 800mAh

产品尺寸: 208.5*28*28mm

重要注意事项

危险

1. 如果电源插头或电源线损坏,请勿使用。

2. 请勿将本产品放入水中或任何液体中。

3. 请勿将本产品放入火中或任何热源附近。

4. 请勿将本产品放入微波炉中。

5. 请勿将本产品放入洗碗机中。

6. 请勿将本产品放入洗衣机中。

7. 请勿将本产品放入烘干机中。

8. 请勿将本产品放入烘干机中。

9. 请勿将本产品放入烘干机中。

10. 请勿将本产品放入烘干机中。

11. 请勿将本产品放入烘干机中。

12. 请勿将本产品放入烘干机中。

警告

1. 请勿将产品放入水中或任何液体中。

2. 请勿将产品放入火中或任何热源附近。

3. 请勿将产品放入微波炉中。

4. 请勿将产品放入洗碗机中。

5. 请勿将产品放入洗衣机中。

6. 请勿将产品放入烘干机中。

7. 请勿将产品放入烘干机中。

8. 请勿将产品放入烘干机中。

9. 请勿将产品放入烘干机中。

10. 请勿将产品放入烘干机中。

11. 请勿将产品放入烘干机中。

12. 请勿将产品放入烘干机中。

注意

1. 请勿将产品放入水中或任何液体中。

2. 请勿将产品放入火中或任何热源附近。

3. 请勿将产品放入微波炉中。

4. 请勿将产品放入洗碗机中。

5. 请勿将产品放入洗衣机中。

6. 请勿将产品放入烘干机中。

7. 请勿将产品放入烘干机中。

8. 请勿将产品放入烘干机中。

9. 请勿将产品放入烘干机中。

10. 请勿将产品放入烘干机中。

11. 请勿将产品放入烘干机中。

12. 请勿将产品放入烘干机中。

故障排除

故障现象

故障原因

解决方法

清洁头无法正常工作

清洁头未正确安装或损坏

重新安装清洁头或更换清洁头

充电器无法充电

充电器未正确连接或损坏

重新连接充电器或更换充电器

产品无法启动

产品未正确充电或损坏

重新充电产品或更换产品

产品无法关闭

产品未正确关闭或损坏

重新关闭产品或更换产品

产品无法调节

产品未正确调节或损坏

重新调节产品或更换产品

产品无法防水

产品未正确防水或损坏

重新防水产品或更换产品

产品无法使用

产品未正确使用或损坏

重新使用产品或更换产品

Instruction Manual

Ultrasonic therapeutic apparatus

· Thank you for purchasing our products.

· Be sure to read this manual carefully before use and keep it safe for future reference.

· This product can only be used for personal dental cleaning at home, not available for any medical purpose.

Product description

Thank you for choosing this product as your oral health care assistant. Please check the product list shown below to check whether the relevant parts are available.

Flat head cleaning head x 2

Tip cleaning head x 2

Charging Cable

Protective cover

Scaler hand

Charging base

Tongue scraper

Probe

Mouth mirror

Small wrench

Protective cover

Scaler hand

Charging base

Tongue scraper

Probe

Mouth mirror

Small wrench

Product use

Warning

When turning on the scaler, do not point the metal cleaning head at sensitive parts of the human body such as eyes, as well as power supplies, lamps and other equipment. When the product is turned on, please touch the metal cleaning head with your skin to avoid injury. Note: Please fully charge the product before using it for the first time, because the main unit is not fully charged when it leaves the factory for safety reasons. While charging, the product cannot be used.

Disassembling the cleaning head

1. Install the cleaning head

Align the cleaning head with the interface on the top of the transducer, use a small wrench to rotate the top of the transducer clockwise to replace or remove.

2. Disassemble the cleaning head

Align the cleaning head with the interface on the top of the transducer, use a small wrench to rotate the top of the transducer counter-clockwise to replace or remove.

Switch mode

Long press and hold the Power Mode switch button to gear and the work indicator will be always on in the gear mode.

After starting, the default mode is soft mode and continuous cleaning mode is available. You can switch the mode by pressing the Power Mode switch button. The product will enter the soft mode or the continuous mode according to the mode you select.

After power on, long press the Power Mode switch button to enter the gear mode. The product will enter the gear mode and the work indicator will be always on in the gear mode.

Clean the teeth

1. When cleaning the calcification, align the scaler hand with the teeth. Do not use the scaler hand to clean the teeth directly.

2. When scraping the surface of calcification, use the scaler hand to scrape the surface of calcification. Do not use the scaler hand to scrape the surface of calcification directly.

3. When using the scaler hand to clean the teeth, use the scaler hand to clean the teeth. Do not use the scaler hand to clean the teeth directly.

4. When using the scaler hand to clean the teeth, use the scaler hand to clean the teeth. Do not use the scaler hand to clean the teeth directly.

5. When using the scaler hand to clean the teeth, use the scaler hand to clean the teeth. Do not use the scaler hand to clean the teeth directly.

6. When using the scaler hand to clean the teeth, use the scaler hand to clean the teeth. Do not use the scaler hand to clean the teeth directly.

Product charging

Connect the USB port on the side of the base charger to the power adapter via USB charging cable, then place the scaler handle into the base to charge.

When fully charged, the charging indicator is always on.

When fully charged, the charging indicator is always on.

Charging Precautions

1. Do not immerse the charging port in water or any liquid.

2. Do not charge, use and place the product in a high temperature place.

3. It is recommended to use a "CCC" certified 5V DC 1A DC output power adapter.

4. This product cannot be used during charging.

Product parameters

Product name: Ultrasonic therapeutic apparatus

Ultrasonic tooth scaler Rated voltage: 5V

Wireless charging current: 1A

Battery capacity: 800mAh

Product size: 208.5*28*28mm

Rated power: 3W

Waterproof rating: IPX6

Product gear: 3-speed

Important Notes

Danger

1. Do not use if the power plug or power cord is damaged.

2. When turning on the scaler, do not point the metal head at sensitive parts of the human body such as eyes, and damaged equipment such as power supplies and lamps.

3. Do not touch the skin with the metal cleaning head for a long time to avoid injury.

4. The product contains a built-in rechargeable battery, please do not heat in fire or charge, use or place it at high temperature.

5. The product is not suitable for people under the age of 18, please keep this product away from children.

6. Others who do not have full capacity for civil conduct or lack relevant experience and knowledge must be under the supervision and guidance of the supervisor use, clean and maintain the product.

Warning

1. Do not immerse the product in water or any liquid. If it falls into the water when connected to the power supply, please unplug the power plug first and then take out the product.

2. Do not press the metal cleaning head hard on the gums or teeth to avoid damage.

3. If you need to charge, please confirm that the power plug and the interface are completely dry, and then connect the power.

4. Do not modify, disassemble or repair this product by yourself.

5. Before connecting the power supply, please check whether the local power supply voltage is consistent with the voltage marked on the product.

6. Do not use the product in a bath or bath.

7. Do not wrap the power cord around the product when it is powered on.

8. Ultrasonic scaling is contraindicated in patients with cardiac pacemakers to avoid symptoms such as dizziness and arrhythmia caused by the interference of electromagnetic radiation.

Note

When using this product for the first time, there may be soreness, itching or slight bleeding. This is because the gums are not used to the efficient cleaning method. It is recommended to use the gentle mode of the product to clean first. This phenomenon will generally stop within 1-2 weeks. If the above phenomenon improves significantly over 2 weeks, stop using this part and consult a dentist in time.

Removal of fault

Troubleshooting Summary

Trouble causes

Cause Solutions

Can not power on

The battery is exhausted or too low in charging status.

Charge timely, using the charging cable.

This cleaning head has a weak vibration

The plug is not properly connected to the power supply.

Insert properly, remove foreign objects.

The handle cannot be connected to the charging base

There is a foreign object in the socket at the bottom of the handle.

Remove foreign objects.

No display of charging

The plug is not properly connected to the power supply.

Check if the power supply is connected.

Will resume to normal after charging

FCC Warning Statement: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC Radiation Exposure Statement

The antennas used for this transmitter must be installed to provide a separation distance of at least 20 cm from all persons and must not be co-located for operating in conjunction with any other antenna or transmitter.

Size:100*80mm
128克的铜版纸