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


VISUAL COLD HOT COMPRESS BLACKHEAD
REMOVAL APPLIANCE

Manual

Product Parts

Product Introduction

The product is a multi-functional blackhead removal appliance with cold and hot compress, visualization.

A: Using vacuum suction technology of the product to stretch and massage your skin gently to help suck out the blackheads, dirt deep in the pore, which makes your skin smoother and shinier.

B: The product is equipped with an HD camera with 300,000 Pixel and a 20x magnifier, which can make you see the impurities such as the blackheads, grease in the pores. Watching while you suck out the blackheads, and clean up the pores without leaving a blind spot.

C: Adopting cold and hot compress to improve skin elasticity and shape firm and plump skin.

Product Parameters

Model: ZSX-H300-P
Product size: 212.0X44.0mm
Net weight: 155g
Charging voltage: 5.0V
Battery capacity: 1000 Mah
Rated power: 6.5W
Pixel: 300,000 HD/ 20-times magnifier
Battery life: 2H
Hot compress temperature: 42°C-45°C
Cold compress temperature: <10°C
Color: Black and white
Suction power: <65 Kpa

Operating Instructions

Large round hole suction head

This suction head with a powerful force that is not suitable for fragile, sensitive skin or the skin around your eyes. It is used to remove blackheads, grease, dirt, cosmetic residues on the skin. It needs to be moved slowly on your skin when you operate the product and please do not use it in the same area for a long time.

Small round hole suction head

This suction head with low force, which is suitable for tender, sensitive skin, and the weak parts of the corner of the nose.

Oval suction head

This suction head is suitable for removing blackheads on cheeks, corners of eyes, and wings of the nose. Use the vacuum absorption function of the product to absorb facial skin, lift the skin, and then quickly remove the product. Repeated operations will promote the blood circulation of the skin, increase skin elasticity, and reduce wrinkles and fine lines.

Microcrystalline head

The head can gently exfoliate aging skin, leaving the skin delicate and smooth.

How to use

1: Press the switch/shift key for 3s, the product enters into the standby state, and the display indicator lights up.

A. Power on

1. Press the switch/shift key once to enter the first gear (low gear) of the blackhead suction mode, and the indicator on the display will light up.
2. Press the switch/shift key three times to enter the 3rd gear (high and medium gear) of the blackhead suction mode, and the indicator on the display will light up.
3. Press the switch/shift key three times to enter the 3rd gear (high and medium gear) of the blackhead suction mode, and the indicator on the display will light up.
4. Press the switch/shift key four times to enter the 4th gear (high gear) of the blackhead suction mode, and the indicator on the display will light up.
5. Press the switch/shift key five times to return to the standby state.

B. Visualization

1. Press the camera key for once, and the indicator will light up.
2. Open the App on your mobile phone/enter into the setting bar//turn on WiFi /Use the corresponding network connection.
3. Exit the Settings to return to the desktop, click the Play button to enter the Visual window to see the picture.
4. Click "record" to save blackhead removal Video transmitted by the product camera in the phone.
5. After clicking "photograph", the images transmitted by the product camera can be saved in the mobile phone in the format of a photo.
6. Then press the camera switch button to turn off the visual function.

C. Power-on State

1. Press the cold and hot compress key for once to enter into the hot compress mode, and the sun indicator will light up.
2. Press the cold and hot compress key twice to enter into the cold compress mode, and the snow indicator will light up.

1: Press the switch/shift key and hold it for 3s to power off.

Precautions :

1. Before using, please dry up your hands to prevent water or other liquids from entering the products from the vents or buttons, thereby damaging the product.
2. Please do not use exfoliating creams and other exfoliating products at the same time.
3. Do not press hard on the skin with the suction head for a long time, which will create a strong pulling force, thus damage your skin.
4. Because of the exfoliating function of the microcrystalline heads, it is recommended to use it once a week for dry skin and twice a week for oily skin, and should not exceed 5 minutes for each use. It depends on the skin.
5. Please do not use the product continuously in the same position, or it will make the skin red and injured. If this happens, please stop using it immediately.
6. This product is not recommended to be used by people who lack five senses, people who have a mental disease, and inexperienced people. Especially not for children.
7. The product is not suitable for people who get cosmetic surgery, such as eyelids surgery/ face lift, nose reshaping, or other cosmetic treatments.
8. The product is not suitable for pregnant women and patients with heart disease.

Operation of blackhead removal

① Face hot compress

Using the hot compress function of the product, or put a warm towel to your face for 3 to 5 minutes to help open the pores of the skin.

③ Suction level choice

Choose a proper suction level for your skin, can apply a little lotion before using (lubrication function)

⑤ Clean up your face

Wash up your face with water after sucking out the blackheads.

⑥ Face care

Put on some toner or a mask, or use the cold compress of the product to help close up the pores.

② Apply the blackhead derived liquid

Drip the liquid onto the makeup cotton and then put it on where the blackheads are for 3-5 minutes. You can use the hot compress function of the product at the same time

④ Adsorption

Put the product on the part where there are blackheads, and then slowly to suck them out (1-2 minutes each time)

Visual function operations

1. Scan the QR code on WeChat to download the corresponding APP

2. Scan the QR code to download the corresponding APP


QR code for ISO QR code for Android system

3. Press the camera key for once, and the indicator will light up.

4. Camera connection

Remarks:
The speed of drawing pictures will be different according to the different mobile phone systems, and pictures can be generated within 3-5s normally.

1. Open the app and click the setting button in the lower left corner to connect to Wi-Fi, and then click the connection button with Sheng Long WAH successfully.

2. After the WIFI connection is successful, return to the main interface and click the drawing button in the middle of the screen to draw pictures successfully.

These suction heads can be removed and rinsed with clean water, but do not make the body of the product directly will water.

Precautions

To avoid any accidents this product is not recommended to be used by the following people:

1. Those who are being treated by a doctor and have abnormalities in their body.
2. Patients with malignant tumors.
3. Pregnant women, puerpera, children.
4. Infectious patients.
5. People with skin sensation impairment or skin abnormalities.
6. Acutely ill patients.

Warranty Card

Product name	Visual cold hot compress blackhead removal appliance	Model	ZSX-H300-P
Buyer Information	Name Address Telephone	Province City/ County District	
Dealer information			
Fault description			
Please fill in the content required on the warranty card in detail when repairs are required.			This product is guaranteed for 1 year

Certification

Certification _____

Quality inspection _____

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.