RingConn

English	1
Deutsch	4
Français	7
Español	10
日本語	14
Italiano	17
简体中文	20
繁體中文	23
Regulatory	26

Please scan the QR code above to
obtain more language versions.

Thank you for choosing RingConn Ring. We look forward to accompanying you through every moment of your smart and healthy daily life. This manual contains information about the RingConn Ring and its charging accessories, including technical details and safety instructions. Please read it carefully and use the product safely to enjoy the best experience.

Important Safety Instructions

This product contains electronic devices. To ensure your safety, please carefully read the following safety instructions before you use it:

1. Do not disassemble the product casing.
2. Do not replace the built-in battery.
3. Do not use corrosive cleaning agents to clean the product.
4. Do not charge the product when it is damp or not completely dry.
5. Do not expose the product to direct sunlight for extended periods.
6. Do not use the product in saunas or steam rooms.
7. Do not place the product in extreme temperature or humidity environments.
8. Do not continue using the product if it is deformed or internally damaged.
9. Do not put the product in a fire as the battery may explode.
10. Please use RingConn's standard charging case to charge the ring, and use a certified power adapter to charge the charging case.
11. If the product becomes excessively hot during charging, stop charging immediately.
12. This product contains a built-in lithium battery. Due to the nature of lithium batteries, prolonged periods of non-use may lead to capacity loss, swelling, aging, circuit board corrosion, and other safety hazards or damage. It is recommended to charge the battery to approximately 70% before storing the product. Keep the ring in the charging case to maintain its battery level, and periodically recharge both the ring and the charging case.
13. The red and green light-emitting diodes (LEDs) in the product emit light within the visible color spectrum. This makes it safer for both the wearer and others nearby.
14. The charging time depends on the initial battery level of the ring and the ambient temperature.
15. Please charge the battery at room temperature. The charging temperature ranges from a minimum of 0°C/32°F to a maximum of 41.5°C/106.7°F.
16. The RingConn Ring communicates with the RingConn app using encrypted Bluetooth communications.

Package Contains

RingConn Gen 3
Smart Ring

Smart Charging Case

Type-C Cable

Instruction Manual,
Quick Guide

- This information is for reference only and does not constitute any form of commitment. Please view the actual product for exact details (including but not limited to color, size, etc.).
- This information is subject to change without prior notice. In case of any inconsistency between this information and the product description on the official website, that on the official website shall prevail.

Product Information

Product Name	RingConn Gen 3 Smart Ring
Product Model	RCP-01-6 RCP-01-7 RCP01-8 RCP-01-9 RCP-01-10 RCP-01-11 RCP-01-12 RCP-01-13 RCP-01-14 RCP-01-15
Battery Capacity Of The Ring	3.78V,13mAh(#6,#7),15mAh(#8~#11),16.5mAh(#12~#15)
Charging Method/Time	Magnetic Charging/approx. 90 minutes for full charge
Transmission Mode	Bluetooth 5.0
Compatibility	iOS14.0, Android 8.0, Bluetooth 4.0 and above
Offline Storage	10 Days
Material	Metallic titanium
Ring Weight (net)	2-3g
Working Temperature	0°C-40°C/32°F-104°F
Protection Grade	IP68

Product Name	Smart Charging Case
Battery Capacity Of The Charging Case	500mAh
Charging Method/Time	Type-C Cable & rated 5W/approx. 120 minutes for full charge
Charging Case Weight(net)	80g
Bluetooth LE Frequency Range	2400~2483.5 MHz
Maximum RF Output Power	0.00050 W

Medical Disclaimer

• Warning:

This product is not a medical device. RingConn and its associated applications are not to be used for the diagnosis, treatment, cure, or prevention of any disease or medical condition. Please consult a medical professional before making any changes to your exercise, sleep, diet, or birth control plans.

• Important Safety Instructions:

1. If you have any pre-existing health conditions that may be affected by the use of this product, please consult your doctor.
2. Ring surface materials are safe for skin contact (ISO 10993 standard).
This product is designed to record your heart rate, blood oxygen, sleep, and other related data. Prolonged use of the product may cause skin irritation. If you have any redness, swelling, itching, or other allergic symptoms, please consult your doctor and seek medical advice promptly.
3. Please keep the ring at the right level of looseness so that it is not too tight to prevent poor blood flow.
4. Do not use or view the product while driving or in situations that require concentration.
5. This product is not a toy and is not intended for use by minors. It contains small parts that may pose a choking hazard if swallowed. Please keep it out of the reach of minors.
6. The performance of this product may be affected by external factors. Please do not place the ring in an environment with excessive heat, cold, or outdoor bright light. Therefore, the performance data in this information is for reference only.

7. Do not pierce the ring or the battery in the center as this could risk the product short circuiting.
8. This product may pose risks to users with certain health conditions. If you experience any of the following, stop using the product and consult a physician:

Epilepsy or sensitivity to light;	Poor blood circulation or abrasions;
Heart disease or other medical conditions;	If you are taking any photosensitizing medication.

Usage And Maintenance

1. The size of your finger changes at different times of the day, and sometimes it may be difficult to remove the ring. If you have difficulty removing the ring, please try the following methods:
 - Wet your finger with mild soapy water or cold water and slowly turn it to remove the ring;
 - Raise your hand above your chest until your blood pressure drops and then try to remove it;
 - In case of emergency and if you feel uncomfortable and are unable to remove the ring on your own, please seek medical advice immediately.
2. Remove the ring to avoid wear and tear when performing strength training, holding a shovel or other heavy tool, or holding objects made of metal, ceramic, or stone.
3. Please do not put the ring into an ultrasonic cleaner to avoid damaging some components.
4. This product is waterproof for everyday use, and can be worn in the rain, while bathing, in swimming pools, and during shallow-water activities like snorkeling. Do not use it in high-temperature and high-humidity environments such as hot springs or saunas (steam rooms). The product has an IP68 waterproof rating (10 ATM/100 meters/328 feet).
5. Avoid applying excessive skincare products on your hands when using the product, as it may cause surface corrosion.
6. Regularly clean the ring with water and mild soap, and wipe it with a soft cloth, paying special attention to the parts that come into contact with your skin.
7. Occasionally remove the ring to give your skin a chance to relax. This is also a good opportunity to charge the ring.
8. The ring may get scratched or scratch soft materials it comes into close contact with, such as skin, gold/silver jewelry, or electronics with soft coatings. Please handle with care.
9. To protect your battery and maintain its service life, please place the ring in its charging case when you stop using it, and charge it and the charging case regularly. It's recommended that you charge the battery at least once every 15 days.
10. To ensure continuous and complete data recording, it is recommended that you sync the data from the RingConn Ring to the RingConn App daily. The RingConn Ring can store data for up to 10 days, though the actual storage duration may vary depending on the ring model and usage frequency. Timely synchronization can prevent data loss.
11. The charging accessories are not waterproof and are for indoor use only. Do not use or store them in damp environments such as bathrooms.

Warranty

RingConn LLC provides a 1-year warranty for your RingConn Smart Ring.

To register your product for warranty please visit:

<https://ringconn.com/return-refund-policy/>

For a warranty claim, we will require:

a copy of the purchase receipt from an authorized retailer.

For general support or inquiries please contact us at

cs@ringconn.com

If your local law in effect at the time of purchase requires a warranty period longer than

one (1) year, this warranty shall be extended to the extent required by such law.

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<https://ringconn.com/pages/euro-compliance>

English

- Hereby, Shenzhen Ninenovo Technology Limited declares that the radio equipment type Smart Health Monitoring Ring is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address:
<https://ringconn.com/pages/euro-compliance>
- The apparatus described above meets the essential requirements of the Electromagnetic Compatibility Directive (2014/30/EU) Electromagnetic compatibility.
- RF exposure safety
This product is a radio transmitter and receiver. It's designed not to exceed the emission limits for exposure to radio frequency (RF) energy set by the Federal Communications Commission. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
- Policy & Privacy Policy Terms: Our Terms of Use (<https://ringconn.com/policies/terms-of-service>) and Privacy Policy (<https://ringconn.com/pages/privacy-policy>) will apply when using RingConn's products and services.
- RingConn charger complies with EN/IEC 62368-1 international safety standard for ICT and AV equipment.

• **FCC Warning**

FCC ID: 2A854JZRCP-01

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with FCC radiation exposure limits set forth for an uncontrolled environment.

• **IC Warning**

IC ID: 31941-JZRCP01

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

The device has been evaluated to meet general RF exposure requirement.

The device can be used in portable exposure condition without restriction.

• **Warning:**

This product has passed the testing under Proposition 65 of California and complies with the relevant safety standards. Changes in the product's usage scenarios may increase the risk of chemical exposure. Before using the product, please be sure to carefully read the product instruction manual and avoid being in abnormal usage scenarios to ensure your health and safety.

The product should be disposed of separately from household trash. The product and the internal battery should be disposed of in accordance with local laws and regulations. Recycling this product helps to preserve natural resources and prevent potential negative impacts on the environment and human health caused by improper waste disposal.

- This product is compliant with Directive 2012/19/EU. Please dispose of used electrical and electronic equipment (WEEE) according to national laws and regulations along with electrical waste. It should not be mixed with general household waste.

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