

SKU: K5102

Item Name	Smart Pill Organizer
Weight	0.33lb
Product Size	4*3.5*1.14in
Material	Food-grade materials PP, ABS
Color	White
Compartments	4
Input	5V/1A
Battery Capacity	1000mAh
Package Includes	1*Smart Pill Organizer 1*USB-C Charging Cable 1*User Guide

The diagram shows the Pillbox in two states. The top part shows the Pillbox open, revealing the internal compartments. Labels with red dashed lines point to the 'indicator light' (a small circular light on the lid), the 'Number of pills' (a small display on the lid), and the 'Pill Compartments' (the internal storage areas). The bottom part shows the Pillbox closed. Labels with red dashed lines point to the 'USB' port (a small port on the side) and the 'ON/OFF' button (a small circular button on the top).

1. Display of number of pills taken;
2. Medication reminder;
3. Storage of medication records;
4. Memory function: the alarm clock set before the shutdown will be restored by default when the kit is turned off and turned on again;
5. Reminding method: sound effect+light reminder+vibration;

1. Short press ON/OFF to start the machine, long press the ON/OFF button for 3s to shut down the machine, press the ON/OFF button three times to enter the DUT mode, and then press the ON/OFF button three times to exit the DUT mode;
2. Display of the number of pills taken: within 30 min after the time point of taking the medicine, only the number of pills taken in the set medicine bin will be displayed after opening the medicine box, and other medicine bins will not be displayed. During this period, the number of pills taken in this period will period will not be displayed if the medicine bin is opened once
3. The cover open: Only records when the phone is connected to the kit will be stored;
4. Bluetooth connect: find the Bluetooth K5J02 Pill Dispenser on your mobile phone and connect it;
5. APP control: Download ComfyTemp App from the Google Play or App Store and register an account(Supports IOS 12 or Android 5.0+). Adjust Bluetooth to "on" status in the "Settings" menu of your mobile phone. Connect the pill organizer in the app and enjoy the app control

6. APP control: Download Comfytemp App from the Google Play or App Store and register an account(Supports IOS 12 or Andro -id 5.0). Adjust Bluetooth to "on" status in the "Settings" menu of your mobile phone. Connect the pill organizer in the app and enjoy the app control.

 Danger	When the surface of the pill dispenser is dirty, wipe it with a The battery of this product is a special rechargeable battery which cannot be used in other products. Please do not charge the battery that has been taken out.
--	--


Danger

- Do not throw into the fire, make it overheated
- Do not pierce with sharp objects, gravitational impact, or disassemble or modify
- Do not connect battery anode and cathode with metal
- Do not store it with necklaces, hairpins, etc.
- Do not charge, use and place in high temperature places such as near the stove or under the scorching sun to prevent heat generation, fire, rupture, etc.
- Please keep the removed battery out of reach of infants and children to avoid swallow by mistake
- When the battery fluid leaks, please do not touch it with bare hands, and handle it as follows:
 - When battery fluid enters the eyes, there is a danger of blindness. Do not rub your eyes directly with your hands, rinse your eyes with clean water immediately and consult a doctor
 - When battery fluid gets on your body or clothing, it may cause skin inflammation or injury. Please rinse with clean water. If you still feel unwell, consult a doctor immediately

Battery Removal Steps:

1. Use professional tools to pry open the shell on the switch side;
2. Use a screwdriver to unscrew the fixing screws of the internal bracket;
3. Unscrew the screws on the circuit board;
4. Use scissors to cut the battery connection leads;
5. Remove the battery

1. Inside with a lithium battery, please keep it away from the high-temperature environment.
2. Please keep this product out of reach of children.
3. As a personal product, please do not mix it with others.
4. Please store in a well-ventilated place, keep away from the environment with strong sunlight and corrosive gas and liquid erosion.
5. Please place the product in a safe location, avoid placing heavy objects on the product to prevent damage.

We offer all customers
5 years warranty,
365 days free replacement,
90 days free return
and 24 hours quick customer service.



24Hrs Quick
Customer Support


www.comfytemp.com

Product Name: Smart Pill Organizer
Brand: Comfytemp
Model: K5102

UK	REP
-----------	------------

Company Name: UKCA SERVICE LTD
Address: 21 Ellesmere Avenue,Worsley,
Manchester,UK,M28 0AL
E-mail: ukca.service@gmail.com

EC	REP
-----------	------------

Company Name: TopToll GmbH
Address: Einsteinstr.2,Neuss 41464,
Germany US: +49-188-DE339651736
E-mail: toppollde@gmail.com

Manufacturer: Shenzhen Yixiaiwanglu Technology Co.,Ltd
Company Name: Xinxiaohendeshizhangshiyuaxiangongsi
Address: Youngseng Second Industrial Zone, Fukang community,
Longhua, 4168M, 4F, 38Bld, 1980 Technology&Cultural Industry
Park,Shenzhen GuangDongShen 518000 CN

FCC Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference it receives, including interference that may cause undesired operation.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

Interference has been tested and found to comply with the limits for a Class B digital device pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference with radio communications.

However, this equipment may cause interference to other equipment. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

RF Exposure Information

The device has been evaluated to meet general RF exposure requirement. This equipment should be installed and operated with minimum distance 20mm between the radio and your body.

157g铜版纸+双面印刷过光油 7.4*10.4cm 5折6页