

CONNECT WITH A REMOTE CONTROLLER

Please follow the steps before first-time usage.

Connection

1



Long press the  button for 5 seconds **on the massager**.

2





While flickering **on the massager**, press the  button on the remote controller once to connect.

Disconnection




Press the  button for 5 seconds **on the massager** and press the  button **on the remote controller**.

※ It is highly recommended to disconnect the massager and the remote controller after use to conserve the battery, as it will stay connected and continuously consume it.



Guide for the Gel Pads



It is recommended change gel pads after 20-25 times usage.

Guide for Mute



Press the Mute button once when the massager is on; Press one more time to disable the muting function.

Guide for Connection



A remote controller can connect with multiple massagers. (Please follow the steps on the left to connect more massagers.)

GUIDE FOR OPERATING

Mode Switch

Power OFF
(hold for 3 sec)

Decrease Intensity
(Low)

<Massager>

※ Each time pressing the mode switch button, the level will be 0; please press + to increase the level.

Power ON (press 2 times)

Increase Intensity (High)

Power OFF

Power ON


Mode Switch

Increase Intensity

Mute

<Remote controller>

Warning



Before taking off the gel pads from your body, turn off the power first.

分割线


分割线

分割线

分割线

背面


PLEASE Do NOT Use THE PRODUCT UNDER THE FOLLOWING CONDITIONS FOR YOUR SAFETY.



PLEASE Do NOT Use THE PRODUCT IN A WET ENVIRONMENT.





PLEASE Do NOT Use THE PRODUCT FOR FOLLOWING PEOPLE.



kure.fit

Kure.fit
Quick Manual



⚠ Please read the manual carefully and follow the instruction.

背面

分割线

分割线

分割线

分割线

FCC Caution.

(1)§ 15.19 Labelling requirements.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

§ 15.21 Changes or modification warning

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

§ 15.105 Information to the user.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.

- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

- Consult the dealer or an experienced radio/TV technician for help.
