ACRS-III WRITER

User Manual

21AE3-0013

※ 품질인증 : ISO 9001:2000 / KC인증 번호 : MSIP-CRM-JVM- JV-RFRM10

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
21AE3-0013

1. Outline

-EEPROM or RFID type WRITER that can save medicine information to CANISTER

< ACRS III WRITER >

This Product Support RFID.

2. Applied to equipment : ATDPS

3. System Requirements

- Frequency : 13.56MHz
- Protocol : ISO15693 (TI Tag-it)

4. How to use the Cartridge


- ① The ACRS III WRITER and PC are connected


- ② Run the "PAM-Pro" Program.


- ③ Click on Medicine Setting.

④ Select "Com Port" and Press "Connect".

⑤ When the RFID Tag is recognized, related information is displayed in the message window. Also, "Read" and "Write" button can be used to store medicine information on the canister

Specification

Specification	ACRS III WRITER
Size (WxDxHmm)	75mm X 130mm x 85mm
Weight (g)	630g
Power	5V
Communication	USB
Temperature	-20℃ ~ +50℃
Frequency	13.56MHz

Certified Company Name	JVM CO., Ltd
Equipment designation	RFID/NFC Radio device (13.56MHz band)
Model Name	ACRS III WRITER
Manufacture Date	Separately marked
Manufacturer / Manufacture country	JVM CO., Ltd / Republic Korea
Service center number	1588-7587

- This device may cause radio interference during use and may cause harmful interference from other device.

WARNING & SAFETY INFORMATION

Safety Symbols

Warning

- When encountering this symbol in the manual, you must follow these recommendations to avoid irreparable damage to your car, system or connected devices or to avoid accidents with injuries or death.

IC Information

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference

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FCC Information

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1.This device may not cause harmful interference, and
- 2.This device must accept any interference received, including interference that may cause undesired operation.

FCC notification to users

This equipment has been tested and found to comply with the limits for a CLASS B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference, the user is encouraged to try to correct the interference by consulting with a dealer or an experienced technician for technical assistance.

Any changes or modifications to the equipment not expressly approved by the party responsible for compliance could void user's authority to operate the equipment.