

CosBEAUTY

Model:CB-S006



aiDirect Pro
Smart Skin Analyzer Refined Version
Instruction Manual



Touch A New Beautiful Life

Contents


Packing List	2
Product Components	3
Safety Precautions	4
Instruction Manual	6
Charging	12
Maintenance	13
FAQs	14
Specs	17
Disclaimer	19
Warranty	20

Packing List

Smart Skin Analyzer Refined Version×1
AC adapter×1
Power cord×1
Pouch×1



Product Components



Safety Precautions

Please read this safety precautions carefully and follow the instruction manual to operate. It is able to ensure the operation of product is correct and safe, to prevent users or other people from personal injuries or property damage. Please make sure to follow it.



- Please do not use the device under any instruction other than this instruction manual, otherwise it may cause breakdown or accident.
- Long-time use or incorrect use of charging device may result in power cord breakage. Therefore, if any following situations occurs, please do not use the device any more: 1. Power cord is twisted and unable to restore; 2. Power cord or power connector is broken, out of shape, cracked.
- Please do not wrap the cord around the power connector.
- Please do not damage the power cord, do not bend, twist, stretch, squeeze, heat, tie, clamp or process it, so as not to cause fire, electric shock and other dangers.
- Do not place the device in fire.



Attention

- Do not immerse the device in water or wash it.
- Do not use it on infants or children, keep it beyond children's reach.
- Do not shoot at eyes, it may cause vision damage.



Warning

- Do not touch or operate with wet hands.
- While not charging, please pull out the adapter plug from socket.
- Please do not use or store the adapter in the environment with high humidity such as bathroom.
- Do not disassemble, repair, transform the device.
- Do not toss or throw the device, the impact may cause damage.
- Please use the original power adapter ONLY.
- Do not use the original adapter on any other products.
- Please stop using the device immediately if any damage shown.

User Manual

Preparation

- Please download App: CosBeauty

Scan to
Download

Supported system :
iOS/Android




- Please make sure the device is fully charged.
- Before testing, please wipe the surface of testing window to make sure it is dry and clean to reduce the effect on testing accuracy.
- Before testing, please make sure the device is connected to the App.
- Before testing, please make sure your skin is clean(cleanse your skin if necessary) and test after 20mins to ensure the accuracy.
- Please refer to App Store or Google Play for upgrading.

Steps

iPhone device connecting steps

- 1 Press the power button, light on means turning on the device successfully, Wi-Fi indicator is lighting on and blinking means that the Wi-Fi is ready to connect.



- 2 Go to <setting> in your smart phone, find Wi-Fi named CBS006_XXXXXX. When connected successfully, white Wi-Fi indicator light will turn into solid light. When connection failed, white Wi-Fi indicator light will remain flashing.



7

- 3 Open CosBeauty App, select one detection dimension, prepare for testing.
- 4 Wipe to clean the testing window, touch the prepared skin gently, wait until the preview image is clear, click <Take Photo> button to capture the picture.

Attention: Do not move the device from skin before the image is captured completely.



- 5 While shooting, you can click button to switch different lights in order to see your skin problem.


Note: 3 different lights are normal lights, which is not harmful to skin at all.

- 6 After shooting, APP will analyze the data that has been collected and offer you the customized solution.

8


- 7 After testing, you can view your testing history and the variation trend of your skin.
- 8 While shooting, you can click preview image to switch different lights in order to see your skin problems.

Note: In case of special circumstances that the device cannot be shut off, please hold the power button for 10 seconds.



Android device connecting steps

- 1 Press the power button, light on means turning on the device successfully, Wi-Fi indicator is lighting on and blinking means that the Wi-Fi is ready to connect.



- 2 Open CosBeauty App, click [aiDirect Pro Smart Skin Analyzer], select any function. App will be connected to Wi-Fi named CBS006_XXXXXX automatically.
 - When connected successfully, white Wi-Fi indicator light will turn into solid light.
 - When connection failed, white Wi-Fi indicator light will remain flashing. App will try again to connect.
 - If the Wi-Fi cannot connect automatically, please go to [setting] and connect Wi-Fi named CBS006_XXXXXX manually.
- 3 When connected successfully, select one detection dimension, prepare for testing.



- 4 Wipe clean the testing window, touch the prepared skin gently, wait until the preview image is clear, click <Take Photo> button to capture the picture.
Attention: Do not move the device from skin before the image is captured completely.
- 5 While shooting, you can click button to switch different lights in order to see your skin problems.
Note: 3 different lights are normal lights, which is not harmful to skin at all.
- 6 After shooting, APP will analyze the data that has been collected and offer you the customized solution.
- 7 After testing, you can view your testing history and the variation trend of your skin.
- 8 After testing, hold the power button for 2.5 secs to shut off the device. When you finishing using, wipe to clean the testing window with soft cloth or tissue to keep it clean, put the cap on, place it into the pouch and store it at a dry cool place.
Note: In case of special circumstances that the device cannot be shut off, please hold the power button for 10 seconds.

Charging

Battery specs

1000mAh Li-ion, charging time 180mins,
working time 110mins.

Power indicator


Low battery: Red light flashing, and indicate
in CosBeauty App.

Charging way

Connect power cord and adapter
(5V/1A) with device to charge. During
using, charging and after fully charged,
indicator light will show red or white
according to the battery power.

Charging indicator

Charging: red solid light.
Fully Charged: red light turns into white
solid light.



Maintenance

- Before and after use, please wipe clean the testing window.
- Please put the cap on after use to avoid dust and dirt.

Storage

- Store it in the pouch when it's not be used, keep it in cool dry place, protect it from light. Avoid contacting with hard objects, prevent from damage. Do not store it in wet or humid environment, otherwise it may affect the product function.
- Power off while not using. Continuously charging may cause battery decrease.
- If you not used it for a long time, you must release the battery power completely every 6 months, then recharge until fully charged.
- As aged lithium batteries can not be repaired, when the device is scrapped, please do not throw away the battery freely. Return it to the recycle station to avoid environment pollution.

FAQs

Please pay attention:

- Gently place the testing window against skin and do not push too hard. Take picture when you can see skin texture clearly, otherwise the accuracy may be affected.
- Sunscreen test can recognize most sunscreen and SPF, there are very few chemical sunblock could not be detected.

How do I know if the picture is clear?

- When you preview the image on the screen, you can see clearly the skin texture, equal light distribution, no blur on any part, thus it is a clear picture. Do not press too hard on the skin, it may cause image deformation.

Is the test result accurate if wearing make-up?

- Testing with make-up will affect the results of complexion and smoothness, we highly suggest you to test without wearing any make-up, therefore, it reflects real skin condition precisely.

How to observe skin when you move device on skin?

- We suggest to hold device with testing window below testing probes direction. When you move the device up and down on skin, the image shown on screen moves in the same direction; when you move the device left and right on skin, you will see the image on screen moves in opposite directions.

Why is it failed to power on?

- Please check if the battery is fully charged, if it's not used for a long time, it may be low battery that caused unable to power on, please fully charge the battery before use.

Unable to connect device

- It may because of low battery, please fully charge the device first.

Why can't it be shut off?

- In case of special circumstances that the device cannot be shut off by holding power button for 2.5 secs, please hold the power button for 10 seconds.

Why is it not able to be charged?

- Please check if the groove on the device is contacting the metal parts on the cord perfectly. If not, please plug the metal part into the groove again.

- If the device is not used for long time, the battery may enter the state of protection, it takes 15mins charging to activate it, after 15mins, the charging light will be on indicating the device is charging. If you still cannot charge your device, please contact CosBeauty after sale service.

Network swiched in Android system

- In AP(Access point) mode, Android system does not support Wi-Fi and carrier network in the same time. when Android phone connected to aiDirect pro device, it is not connected to network, for your convenience, we offered automatic switching function, due to system limitation, the connection speed is rather slow, we suggest to finish all the testing first, then you can check other articles and activities in the APP.
- If the APP is not connected for long time, you can manually connect to "CBS006_XXXXXX" WiFi in setting.

Upgrade

- Please refer to CosBeauty official website for product hardware upgrade information: <http://www.cos-beauty.com>
App is upgrading all the time, please subject to App guidance if any changes appear.

Specs

Dimension	105.50mm×37.00mm×58.50mm
Net weight	100g
Battery	1000mAh lithium battery
Working time	110mins
Charging time	180mins
Rated power	2W
Working voltage	3.7~4.3V
Rated input	5V ---1A
Working temperature	0~45℃

Storage temperature	-10℃~60℃
Working humidity	≤90%RH
Material	ABS, Metal
Available color	White, Pink, Blue
Picture spec	2 megapixel, 100times zoom
Supported system	iOS 7.0 or above, Android 4.3 or above

• This product uses lithium battery.

If the country you stayed have specific regulation, please put the used battery in regulate area for nuified disposal.

Disclaimer

CosBeauty reserves the final explanation of CosBeauty aiDirect Pro Smart Skin Analyzer Refined Version, any third party applies device related information without CosBeauty's authorization is not allowed, we will take legal action. In any case, CosBeauty Technology Co.,Ltd is not responsible for any consequences caused by improper use of consumers.

CosBEAUTY

Warranty card

Name : _____ (Mr./Ms.)

Tel : _____

Address : _____

Purchase date : _____

Serial no. : _____

Model no. : _____

Purchase store : _____

Retailer tel. : _____

After-sale service hot - line : 400 9600 630

Signature _____

Warranty Coverage

1. Warranty period:

12 months beginning on the date of purchasing.

2. For the duration of the Warranty Period, CosBeauty warrants the Product against defects in materials and workmanship arising from Normal Use of the Product.

3. During the warranty period, if any of the followings occurs, certain maintenance costs will be charged:

- Product damage due to fire, flood, abnormal voltage or other natural disasters.
- The malfunction and damage caused by incorrect use, the unauthorized repairs and alteration.
- Not able to present warranty card and purchase receipt.
- The user information on warranty card is not recorded clearly or correctly.
- Product malfunction caused by the fall and improper transportation after purchase.
- Product damage caused by disassembling and repairing by user without being checked by maintenance staff from authorized store.
- The warranty card is missing product model no., serial no., purchase date or retailer information, or texts have been changed.

4. Whether or not during the warranty period, if any of the followings occurs, no maintenance provided:

- Purchase for business instead of individual consumption, it is not covered by warranty terms.
- Color fading, wearing, consuming under proper use is not covered by warranty terms.
- If product model no. and serial no. filled in warranty card is not consistent with official product model no. and serial no., it is not covered by warranty terms.

5. Notice:

- This warranty card is limited valid in the mainland of the People's Republic of China [not valid in Hong Kong, Macao, Taiwan].
- This warranty card cannot be reissued in any case, please keep it properly and present it when repairing.
- If you have any question regarding after-sale service, please contact us or the retailer.

FCC Caution

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.