

说明书1

正



Instructions for iite S1 Toothbrush
www.iite.co

Prohibitions!

- Do not modify, dismantle, or repair the S1 toothbrush (hereafter referred as "The Product") on your own.
- Do not soak the charging cradle in water or any other liquid.
- After cleaning the Product, make sure the plug of the charging cradle has been dried completely before connecting the power.
- The Product has a built-in rechargeable battery, so do not throw it into fire, heat it, or charge, use, or discard in high-temperature environment.
- Do not knock or drop the Product and do not use it any more if it has been damaged due to collision or drop.
- Keep the product far from electromagnetic fields, such as induction cooker when it is being charged.

Cautions!

- Check for visual damages each time before using the product. Do not use it if it has been damaged because of drop or any other reason.
- Before connecting the product to power, check whether the built-in charge is equal to the rated voltage of product.
- If the charging cradle has undergone any damage, do not use it.
- Do not use the charging cradle outdoors or close to high-temperature surfaces.
- If any warning sign is found on the product including brush head and handle, do not use it any more.
- Do not use the product when bathing or showering.



If any operation was made on your mouth or gum in the past two months, contact your dentist for advice before using this product.
If your teeth are found bleeding upon using this Product or the bleeding lasts over one week, contact your dentist for help.
If a heart pacemaker or other device has been implanted into your body, contact your surgeon or the device manufacturer before using the Product. Do not use a toothbrush that is not recommended by the manufacturer.
This Product can be used only to clean the teeth, gum and tongue. Do not use it for other purposes. If you feel uncomfortable or painful, stop using and go to see a doctor for help.
If you use the toothpaste containing hydrogen peroxide, sodium bicarbonate, or bicarbonates locally, "whitening toothpaste", clean the brush head completely with soap and clear water after brushing teeth every time.

Notes!

- Do not place any metal object, for example, coin or clip on the charging cradle.
- Do not leave the power lines around this adapter tightly for storage.
- Unplug the power the Product off if it is not be used for a long time, preventing fire or other hazard.
- Do not share one toothbrush with other family members or other people.

Note: Although your teeth are in good condition, bleeding or bleeding may happen at the beginning when you use this Product, because your gum has not been accustomed to the stimulus for the first time. Once the bleeding or bleeding stops for one or two weeks, otherwise your gum may have something wrong, and you should stop using the Product and consult a dentist for help.

Components



How to Use

- Perform the following steps before the first use:
- Place the charging cradle at a safe, flat, dry, and shadowy location to avoid falling into water or being in high-temperature environment.
 - Connect the power plug.
 - Insert the handle vertically into the charging cradle with buttons facing you. If the charging indicator goes on, the Product is being charged.
 - Notes:
 - Charge the Product for 24 hours before the first use.
 - It is recommended to use the Product in temperatures ranging from 5°C to 35°C, and moisture ranging from 30% to 85% RH.
 - Charge the Product under the following circumstances:
 - Before the first use.
 - Have not been used for at least three months.
 - The toothbrush vibrates less strongly than usual when used.
 - To avoid being affected by low battery during usage, place the Product on the charging cradle after each use. Do not worry about overcharging the battery or reducing the battery life.
 - After the increase of charging time, the Product can be used for fewer times after each charging. If the Product is fully charged and out of power after being used only once, the battery can no longer work.
 - LED indicator status corresponding to different battery levels:


State of Green LED	Battery level		
	Over 90%	20%-91%	Less than 1%
	Off	Flash slowly (every 1s)	Flash rapidly (every 0.2s)

5.5cm

反

For daily use, perform the following steps:

1. Cover the brush head on the shaft by lip. (See Figure A).




2. Wash and squeeze some toothpaste onto the brush head.

3. Place the brush on the gumline in a small circle. (See Figure B).



- Move the brush along the teeth slowly and tenderly and do not use excessive force.
- Cleaning teeth: Hold the handle tenderly with fingers, put the brush head into mouth, switch the brush on, keep brushing in a moderate angle, and click on the connecting teeth and gum. Rotate the brush head slightly on a fixed location for several seconds and then change to another location and rotate again. Perform the same operations on every part of teeth without any omission, including the internal surface and the chewing surface, as well as the rear of back teeth. Move the brush along the gum line, no press or rub. This is the whole process of brushing teeth.
- Clean the brush head and handle. Wash the brush head and handle with clear water and knock the brush head tenderly to remove food debris. (See Figure C).
- Put the brush head as shown in Figure D and wipe the head to make it dry. (See Figure D).



Charging

- This Product is charged in inductive manner. Put the brush on the standard-configured charging cradle that has been connected to power and then the charging begins.
- Charging indicator in the process of charging: the LED indicator flashes slowly and keeps lighting when the battery has been fully charged. The charging takes 24 hours for the first time and 7 hours each time thereafter.
- When the charging begins, the brushing mode is shut off, the **On/Off** function button is disabled (can not be pressed), and the motor stops rotating.
- During charging, if the handle leaves the charging cradle, the Product will start standby mode (the LED is off and the motor does not rotate).
- When the Product is switched on and the handle is placed on the charging cradle, the motor stops rotating and the charging mode begins.

Bluetooth

State of blue LED	Bluetooth status		
	Red light on (every 1s)	Blue light on (every 1s)	Blue light on (every 1s)
	Off	Flash rapidly (every 1s)	Keep lighting

Cleaning

Charging cradle

- Wipe the charger with soft cloth.
- Wipe the metal part of the plug with dry cloth every 6 months to avoid poor contact caused by accumulated dust.
- Brush head:
 - Take the brush head off the handle.
 - Wash the brush head and brush with clear water.
 - Take the brush head off the handle and wash the connection part at least once a week.
 - Handle:
 - Take the brush head off the handle.
 - Wipe the outer surface of the handle with soft and wet cloth and put back on the charging cradle after it has got dry.

Product Information

Name: iite S1 led S1 Super Toothbrush
Model: iite-S1M
Rated voltage: 220~240V AC
Rated frequency: 50/60Hz
Rated input power: 2W
Product standards: GB4706.1-2005, GB4706.15-2008
Origin: China

Battery life

Around three years, if the number of times that the battery can be used after being fully charged has reduced significantly, the battery life has possibly come to an end. (The battery life varies according to usage and storage conditions.)

Troubleshooting

Faults	Possible causes	Corrective measure
The indicator fails to light during charging	Product is not fully charged or connection between the engine	After being in poor charging, the indicator may neither flash nor keep on. After it is used normally, it will start charging in a few minutes. If after a certain long duration the indicator continues to light, contact the store where you buy the product or the maintenance center authorized by the repair service.
The product battery can only stand for several minutes even after being fully charged.	Insufficient charging time	Charge for over 10 hours.
The brush vibrator is abnormally weak.	The battery life comes to an end.	If repair is needed, contact the store where you buy the product or the maintenance center authorized by the repair service.
The brush does not work and the motor does not turn after being fully charged has been pressed.	A new motor has been installed on the brush.	Place the brush along with the handle.
The brush does not work and the motor does not turn after being fully charged has been pressed.	The battery life has come to an end.	Contact the store where you buy the product or the maintenance center authorized by the repair service.

Warranty

Dear iite toothbrush users:
Thank for using our iite rechargeable toothbrush. In order to provide the optimal service for you, to ensure satisfactory service, please read the following terms and keep the warranty valid only in PRC and Hong Kong, Mexico and Taiwan are not included in the scope of the warranty.

Warranty terms

- This warranty proves that our company fully produces the product.
- Liabilities for the products within the warranty period take effect from the purchasing date.
- (No repair is needed when the product is purchased).
- Within the warranty period, the products can be repaired for free at sites designated by the site upon issue of the warranty and proof showing the purchasing date.
- Pay special attention to the following:
 - When the warranty period is not released or missed.
 - The warranty is valid only in PRC and Hong Kong, Mexico and Taiwan are not applicable.
 - Some products have repair services for products subject to fault or failure caused by normal use from the purchasing date within the following period:
 - One year from the purchasing date (Products for commercial use are excluded).
 - Under the following circumstances, the repair service can only be provided upon payment even in the warranty period:
 - The product is beyond the warranty period, that is, one year.
 - No warranty or other proof showing purchasing date can be issued.
 - The product has been caused by force majeure (Such as fire, earthquake, typhoon, or war).
 - Some visual defects, such as scratches and stains appear during usage.
 - The product has been caused by unreasonable damage, or incorrect repair.
 - Consumables including LED and other accessories are beyond the repair service scope.
 - The product has been caused by other uses, such as commercial use, instead of home use, in the PRC (except Hong Kong, Mexico, and Taiwan).
 - Customer service center telephone: 4006610604 (beijing) Official website: <http://www.iite.co>


说明书2

正



iite S1 quick start guide
www.iite.cc

Components




Protection case for travelling

Brush head

General type (standard configuration)

Body



On/off Shift button

Blue indicator (blue)


Charging indicator (green)

Handle

Charging cradle

How to use


Step 1: Put the motor-driven brush on the charging cradle for 24 hours before the first use.




24 hours

Note: This product applies only to children above 4 years old.

Step 2: Cover the brush head on the shaft tightly.




Step 3: Wet the brush head and squeeze some toothpaste onto the head.




反

Step 4: Place the brush head on the line connecting teeth and the gum.




45°

Step 5: Switch on and begin to clean your teeth. Make sure every part is cleaned.




Step 6: To avoid sprinkling of the toothpaste, switch off and take the brush out of your mouth.




开关

Step 7: Wash the brush head with flowing clear water after brushing teeth, wave the brush head to make it dry, and wipe the handle with soft cloth.



Step 8: Put the toothbrush back on the charging cradle.



Note: If the brush has not been used for a long time, unplug and clean the plug and then store it in a shadowy and dry location with no direct sunlight.

Download the iite smart toothbrush App as follows:

a. Swipe the following QR code.



b. Visit www.iite.cc and download the App.

c. Search for iite in Apple's App Store, SnapRes, or 91 Mobile Assistant and download the App.

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.