Model # 5330887

TAYLOR®

Bluetooth weight scale

INSTRUCTIONS FOR USE & WARRANTY INFORMATION

CUSTOMER SERVICE INFORMATION PLEASE READ!

For questions regarding this product, please contact us at:

TAYLOR division of Lifetime Brands, Inc.
2220 Entrada Del Sol, Suite A | Las Cruces, New Mexico 88001 USA
1-866-843-3905 | www.taylorusa.com
Email: taylorusa@lifetimebrands.com

We value your thoughts on our product and services. Please post an on-line review to help us continue to meet your needs and expectations.

Please contact us before returning to retailer.

Preparation Before Use

INSERT THE BATTERIES

1. Use 3 AAA (1.5V) batteries (included). Remove any plastic wrap from the batteries before proceeding.
2. Open battery compartment cover on bottom of scale.
3. Place the new batteries into the battery compartment according to the polarity illustration in the compartment.
4. Replace the battery compartment cover.

This scale can be used without the app. How to use as a weigh only scale:

WEIGHT UNIT SETTING

- Your scale is set at the factory to weigh in pounds (lb).
- You may select kilograms units of measurement by pressing the button on the bottom of the scale while the scale is on (press lightly on the platform to turn the scale on).

SET THE SCALE IN THE RIGHT POSITION

- Always use the scale on a flat and hard floor surface. Do not use on carpets.

START-UP THE SCALE

- Press firmly once on the scale platform to start-up the scale. The display shows "0.0" and then turns off. The scale is ready for use.


Note: If you move the scale at all, please repeat the steps above for best accuracy.

Press on platform to turn on scale

Allow the scale to turn off

Step evenly on platform

Weight Display

- The scale will turn off automatically.

SETTING UP YOUR BLUETOOTH® SMART SCALE USING BLUETOOTH® 4.0 TECHNOLOGY

WeighSync™ application transmits your results to an app on a mobile device using Bluetooth® 4.0 technology.

BLUETOOTH® WEIGHT SCALE REQUIREMENTS

The Bluetooth weight scale is designed to be used with the following mobile devices:


- **iOS**
iPhone – using iOS 13.0 and newer
- **Android**
9.0 (API 28) and newer
- **Compatible with devices running Bluetooth® 4.0 technology and above**

Scan the QR code to download the WeighSync™ app from the App Store (iOS) or Google Play (Android). The QR code will have the most up to date version.

Quick Overview of Using the App

1. Your Weight Tracking Smart Scale will look for the mobile device automatically. Please make sure the Bluetooth setting is "ON" and your device is within 3-30 feet of the scale.
2. The Privacy Agreement in the app must be accepted to use the app with the scale.
3. Edit the user profile to set up your personal data, including activity levels and goal weight.
4. Proceed to the next step to store your user profile.
5. If you would like to change personal data, including activity levels and goal weight, simply edit the user profile again.
6. An Onboarding/Tutorial option will appear. Review or skip to move to the next step.
7. Add device to the app by stepping on the scale and pressing 'New' in the upper right hand corner in the app.
8. Select your user slot.
9. To take your weight reading, in the app tap on the scale icon and then step directly on to the platform. Your weight will appear on the scale. Daily calorie estimates, BMI, and weight difference compared to goal weight will appear on the app.

Due to continuous improvements the app may have some modifications.






Summary of Getting Started with app set up

1. Choose the type of device
2. Download app
3. Open the app
4. Allow WeighSync™ to access Bluetooth

Note: A verified email cannot be changed once verified.

CREATE YOUR PROFILE

1. Register a new account or log in if you have an existing account.
2. Enter your email address to use it as your login credentials. Create a username and password for your profile.
3. You must accept the privacy policy to use this app.
4. Enter your personal profile here.
 - Birthday
 - Gender
 - Height
 - Select Language (English or Spanish)


5. Select Activity level
Sedentary : Your day involves mostly sitting and minimal physical activity. You have limited exercise routines.

Limited Activity:
While you're not completely inactive, your daily routine mostly involves light movements and occasional walks.

Moderate Activity:
You engage in regular physical activities like brisk walking or light workouts, maintaining a balanced routine.

Very Active: Your lifestyle includes frequent exercise sessions and active pursuits, contributing to your overall well-being.

Extremely Active:
You lead a highly active lifestyle, involving intense workout, physical challenges, and a strong dedication to fitness.




6. Enter your target weight.
Use the slider or tap on box to type in your target weight.

ONBOARDING/TUTORIAL

This will walk you through the application on and where to find data within the app. Press 'Skip' to bypass this step.

PAIR YOUR DEVICE


1. See below for the device page option. This is where you add or remove the scale from your personal device.
2. To pair your device, make sure your phone Bluetooth is on, the display on the scale is lit up, and the Bluetooth light is blinking on the scale. If the scale is not connecting to the app, click on cancel and try again. If the scale continues to not connect after 3 tries, please contact technical support.


USING YOUR SCALE AFTER PAIRING

1. After you successfully paired a device, choose a memory slot to store your data to use it offline when app is not open.


Up to 50 sets of data can be stored on the scale before the oldest data is recycled.

2. Press on the scale icon to begin taking your weight.

3. After a successful measurement, the app will calculate your BMI, as well as estimate Cal-Max, if you enter the necessary personal data. Press "Take a Weight" to weigh again.

What is the Cal-Max™ daily calorie intake estimator?


The Cal-Max™ function estimates the number of calories required based on your body composition and user entered personal data. This tool can be used as a guide when setting calorie goals during weight loss and exercise programs. This estimate shows on the app as Kcal (Basal Metabolic Rate).

Sharing your data with Apple Health or Health Connect

Go into settings, check the box next to 'Health' to enable syncing to either Apple Health or Health Connect. The Apple Health or Health Connect app will then ask you to grant access and update your health data.

Verify Your Account

Verifying your account allows you to sync data to the cloud.



Go to settings tab and select My Profile

Tap here to verify your account. A verification email will be sent. A notification that the email has been sent will appear at the bottom of the screen. If you do not see this, you may need to tap verify again to prompt the email.

Tap on the verification link in the email to to verify your account.

Your email


Once successfully verified this page will show.

"Check your spam folder if the email does not appear"

Records


The Records page is where your weight and BMI records are stored, and you can access your history.

- Daily-Shows the average weight for each day in a week.
- Weekly-Shows the average weight for each week in a month.
- Monthly-Shows the average weight for each month in a year.
- Detailed information for each measurement shows below the graph as a list.
- Press < or > to change the date range view.
- Press Today and it will send you back to the current date.

App Settings


- My Profile - To update your personal profile, as well as your password and account verification.
- Sync Data -Toggle on and off to enable data upload and download function **(only available for registered and verified account)**.
- FAQ -Frequently asked questions.
- Reset Usage Tutorial -View the user tutorial again.
- Feedback -To send feedback and questions to technical support.
- Privacy -View the Privacy Policy and User Agreement.
- Sign Out -To sign out of current account.

Battery Warnings!

- **Always purchase the correct size and grade of battery most suitable for the intended use.**
- **NOTE: Please recycle or dispose of batteries per local regulations.**
- **WARNING: Batteries may pose a choking hazard. As with all small items, do not let children handle batteries. If swallowed, seek medical attention immediately.**
- **PRECAUTION: Do not dispose of batteries in fire. Batteries may explode or leak.**
- **Remove the batteries if the product will not be used for a long period of time. Remove used batteries promptly.**
- **Clean the battery contacts and also those of the device prior to battery installation.**
- **Non-rechargeable batteries are not to be recharged.**
- **Do not force discharge, recharge, disassemble, heat or incinerate. Doing so may result in injury due to venting, leakage or explosion resulting in chemical burns.**
- **Always replace all batteries at the same time; do not mix old and new batteries. Do not mix alkaline, standard (carbon-zinc) or rechargeable (ni-cad, ni-mh, etc.) batteries.**
- **When batteries of different brand or type are used together, or new and old batteries are used together, some batteries may be over-discharged due to a difference of voltage or capacity. This can result in venting, leakage, and rupture and may cause personal injury.**

Error Icons

This scale has a low battery indicator. It will display “Lo” or readings will grow dim or become irregular if the batteries need replacement.

“Err” indicates an overload warning. The maximum capacity of this scale is 400 lbs / 180 kgs. Should “Err” appear on the display during a weighing, this indicates that the maximum weighing capacity of the scale has been exceeded. Remove the weight immediately; otherwise, permanent damage may occur.

Precautions

- Clean exposed parts with a soft, slightly damp cloth. To remove stains, use a mild soap. Never use detergents, excess water, treated cloths, harsh cleaning agents, or sprays. Do not attempt to lubricate or open scale casing as this will void your warranty.
- Although your scale is designed to be maintenance-free, it contains sensitive electronic parts. Avoid rough treatment. Do not drop, kick, or jump on it. Treat it with care to ensure the best performance.
- Do not store scale where you store cleaning chemicals. The vapors from some household products may affect the electronic components of your scale. Do not store scale on its side.
- This scale is an extremely sensitive weighing device. In order to prevent run down of the batteries, do not store anything on the scale.

User Tips

1. Place your scale on a hard, flat surface, do not use on uneven tile, flooring or carpeting. Try to keep your scale in the same location or position as much as possible.
2. Step on the scale with one foot applying enough pressure to turn the scale on - then remove your foot.
3. The digital readout should show dashes or zeros.
4. The scale will automatically turn off.
5. Now your scale is ready for you to take a weight reading. Step on the scale evenly with both feet. **NOTE:** If you've weighed yourself and a “C” is displayed after the zero, weigh yourself again for a more accurate reading. Your weight cannot be locked-in if you do not stand still. If the scale displays “Err”, you must step-off and weigh yourself again.
6. If you pick up, move or bump your scale, please repeat the steps above for best accuracy.
 - Weigh yourself without clothing whenever possible. Clothes can vary in weight.
 - Move the scale to a location where it won't be bumped. If your scale is bumped while being activated or while in use, an inaccurate reading may occur.
 - For the most consistent readings, weigh yourself at the same time of day. Allow for weight fluctuations due to food or drink in your system.

One (1) Year Limited Warranty

This scale is warranted against defects in materials of workmanship (excluding batteries) for one (1) year for the original purchaser from date of retail purchase. It does not cover damages or wear resulting from accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair. Do not return to retailer. For product support & warranty information: Please contact us with proof of purchase (purchase receipt) at:

TAYLOR
2220 Entrada Del Sol, Suite A
Las Cruces, New Mexico 88001 USA
Phone: 1-866-843-3905
Email: taylorusa@lifetimebrands.com

There are no express warranties except as listed above. This warranty gives you specific legal rights, and you may have other rights which vary from state to state. For additional product information, or warranty information outside the USA, please contact us through www.taylorusa.com.

©2024 Taylor, a division of Lifetime Brands, Inc. All rights reserved. Oak Brook, IL, 60523

Made to our exact specifications in China.
Not legal for trade.

FCC ID:2AJ5E-5330887

Scan Here for WeighSync™
User Agreement & Privacy
Statement

REF# 5330887

CP5330887 25.25 WC

FCC INFORMATION

Party issuing Supplier's Declaration of Conformity & Responsible party – U.S. Contact information:

Taylor division of Lifetime Brands Inc.
2311 W. 22nd Street, Suite 200
Oak Brook, IL 60523

Equipment:

Product Name: Bluetooth weight scale

Model No.: 5330887

Trade Name: Taylor

We declare under our sole responsibility that the above referenced product has been tested according to FCC requirements and this device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Caution: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on. The user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter.

Radiation Exposure Statement

The device has been evaluated to meet general RF exposure requirement in portable exposure condition without restriction.

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© 2024 Taylor división de Lifetime Brands Inc. Todos los derechos reservados.
Taylor división de Lifetime Brands Inc.