

DURING MEASUREMENT

- Place your feet in the same position every time you weigh yourself. Stand still during measurement.
- Weigh yourself without clothing whenever possible. Clothes can vary in weight.
- Place the scale where it won't be bumped while being activated or while in use, as an inaccurate reading may occur.
- For the most consistent readings, weigh yourself at the same time of day, for the best results early in the morning, as weight can fluctuate due to food or drink in your system.
- User Type Modes** – Standard/Athlete/Baby/Pet can be selected in the Precision Hub App. See the Helpful Hints guide for app instructions.

FCC INFORMATION

Party issuing Supplier's Declaration of Conformity & Responsible party – U.S. Contact information:

Taylor division of Lifetime Brands Inc.
2311 W. 22nd Street, Suite 200
Oak Brook, IL 60523

We declare under our sole responsibility that the above-referenced product has been tested according to FCC requirements, and this device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: [1] this device may not cause harmful interference, and [2] this device must accept any interference received, including interference that may cause undesired operation.

Equipment:

Product Name: Bluetooth Body fat scale

Model No.: 5341063

Trade Name: Taylor

Caution: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on. The user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC RF Exposure Warning Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

This equipment can be used in portable exposure conditions.

Operation is subject to the following two conditions:

- [1] This device may not cause interference.
- [2] This device must accept any interference, including interference that may cause undesired operation of the device.

Not for diagnostic use. Not a medical device.

ISED Radio Frequency Exposure Statement:

This equipment complies with ISED exposure limits set forth for an uncontrolled environment.

This equipment can be used in portable exposure conditions.

This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter.

Radiation Exposure Statement

The device has been evaluated to meet general RF exposure requirement in portable exposure condition without restriction.



1-866-843-3905 (8-5pm MST, M-F)

www.taylorusa.com

taylorusa@lifetimebrands.com

©2025 Taylor, a division of Lifetime Brands, Inc. All rights reserved. Oak Brook, IL, 60523

Made in China.



Not legal for trade.

IB5341063
CP5341063 12.25 WC

Model# 5341063

TAYLOR

Bluetooth Body fat scale

with Body Fat, Body Water, Muscle Mass Estimation, Weight Tracking, KCAL (BMR) and BMI

INSTRUCTIONS FOR USE


IMPORTANT SAFETY INFORMATION

The TAYLOR Smart Body Composition Scale is intended for home use for weight measurement and estimation of body composition metrics for informational purposes only. The scale and the associated Precision Hub™ app are designed to be used as tools for tracking changes in body composition over time to support a healthy lifestyle.

This product is not intended to diagnose, treat, mitigate, or prevent any disease or health condition.

The measurements and data provided by this device are estimates and should not be used as a basis for medical diagnosis or treatment. The information provided by this scale is not a substitute for professional medical advice, diagnosis, or treatment. Users should always consult with a qualified healthcare professional for any health concerns or before making any decisions related to their health. This includes, but is not limited to, starting a new diet or exercise program. Do not disregard professional medical advice or delay in seeking it because of something you have read or interpreted from the measurements of this device or the associated application.


Scale Setup & Calibration



- Insert 3 AAA (1.5V) batteries (included).
- Place the scale on a hard, flat surface and toe tap to turn the scale on and to calibrate.
- Scan QR code and download the free Precision Hub app.
- Set up your user profile within the Precision Hub app.
- Non-app user should manually create the user profile (P1-P8) from the scale before measuring.



Precision Hub app





- Scale is ready to weigh.

Bluetooth Body fat scale


This scale is designed to be used with the following mobile devices:
iOS: iPhone using iOS 15.6 and newer
Android: 14.0 API level 34 and newer

Getting to know your scale

Front & back view of scale




LED readout description



1. Body Water % Reading
2. Muscle Mass % Reading
3. Basal Metabolic Rate (kcal)
4. BMI Readout
5. Weight Tracker Indicator (as compared to previously saved measurement)
6. Body Fat % Reading
7. Weight Readout in Kilograms (kg)
8. Weight Readout in Pounds (lb)
9. Height Displayed in Centimeters (cm)
10. Height Displayed in Feet/Inches
11. User Number
12. Age
13. Athletes or non-athletes
14. Gender
15. Current Weight (C) or Previous Weight (P)
16. ▲ Forward Button
17. SET Button
18. ▼ Back Button


Proper Foot Placement

Stand barefoot on the scale, ensuring your feet make full contact with the metal electrodes.



Interpreting body composition readings

If your body fat percentage, for example, is 25.8% and you need some help understanding these results, please refer to the app for more detailed information.



Care & Maintenance

1. The platform will be slippery when it is wet. Keep it dry! However, do not expose to direct sunlight, extreme hot/cold temperatures, or humidity.
2. This scale is designed to read weights between 11 and 440 lbs / 5 and 200 kgs. To prevent damage, do not place more than 440 lbs / 200 kgs on it.
3. Treat the scale with care as it is a precise instrument. Always keep your scale in a horizontal position. Do not store the scale on its side.
4. Clean the scale with a damp cloth and avoid water getting inside the scale.
5. DO NOT use chemical/abrasive cleaning agents.
6. Always keep your scale in a horizontal position.
7. If your scale fails to turn on, check whether batteries are installed or battery power is low. Insert or replace with new batteries.
8. Do not use for commercial purposes and follow all instructions provided in this manual.

DETAILED INFORMATION ON HOW TO INSERT THE BATTERIES

1. Insert the 3 'AAA' (1.5V) batteries (included).
2. Remove any plastic wrap from the batteries before installing.
3. Open the battery cover on the back of the scale unit.
4. Place the batteries into the battery compartment according to the "+" and "-" polarity symbols in the battery compartment.
5. Replace the battery compartment cover and screws if applicable.
6. If you do not intend to use the scale for a prolonged period of time, it is advisable to remove the batteries before storing.
7. Do not mix old and new batteries.

Battery Warnings!

- Always purchase the correct size and grade of battery most suitable for the intended use.
- NOTE: Please recycle or dispose of batteries per local regulations.
- WARNING: Batteries may pose a choking hazard. As with all small items, do not let children handle batteries. If swallowed, seek medical attention immediately.
- PRECAUTION: Do not dispose of batteries in a fire. Batteries may explode or leak.
- Remove the batteries if the product will not be used for a long period of time. Remove used batteries promptly.
- Clean the battery contacts and also those of the device prior to battery installation.
- Non-rechargeable batteries are not to be recharged.
- Do not force discharge, recharge, disassemble, heat above 55°C, or incinerate. Doing so may result in injury due to venting, leakage, or explosion, resulting in chemical burns."
- Always replace all batteries at the same time; do not mix old and new batteries. Do not mix alkaline, standard (carbon-zinc), or rechargeable (Ni-Cad, Ni-MH, etc.) batteries.
- When batteries of different brands or types are used together, or new and old batteries are used together, some batteries may be over-discharged due to a difference of voltage or capacity. This can result in venting, leakage, and rupture, and may cause personal injury.

Troubleshooting Common Issues

1. If the scale does not turn on or no weight is displayed?

Check if the scale is powered up and started up. If not, please refer to the Section "Scale Setup & Calibration".

2. If the scale is inaccurate or the weight varies?

Place your feet in the center of the platform and stand still. Allow the weight to lock in before looking down at the display. Flooring may be a reason for fluctuation. Use the scale on the hardest floor available for accuracy. Do not move the scale unnecessarily.

3. I have tried all corrective actions, but still can't solve the problem ...

Please contact the Taylor Customer Service Department at 1-866-843-3905 (from 8:00 am to 5:00 pm, Mountain Standard Time, Monday through Friday) or by email to taylorusa@lifetimebrands.com.

4. Why do I get a different body fat reading when I use a different brand of body fat scale?

Different body composition scales take estimations around different parts of the body and use different mathematical algorithms to estimate the percentage of body fat. The best advice is not to make comparisons from one device to another, but to use the same device each time to monitor any change.

5. Profile Selection When Weights Are Similar

- If two users have weights within ± 6.6 lb, the scale will display both profiles (for example, P1 P2). Press the Up arrow to select P1, or the Down arrow to select P2.
- If the two users are a scale-only user and an app user with similar weights (for example, P1 and A1), the display will show P1 A1. Press the Up arrow to select P1, or the Down arrow to select A1.

6. Why isn't the app syncing my offline measurements?

The app will display a "Synchronize Data" button only if your specific scale requires manual offline measurement transfers. Tap the "Synchronize Data" button and wait for your offline measurements to upload. Some models sync automatically - in those cases, you won't see the button at all.

7. Issues with device connection?

- The scale can only be connected to one cell phone and app at the same time. If the Bluetooth icon is shown on the scale, it may have been connected to another cell phone.
- If the scale is already connected to another cell phone, you can close the app or temporarily turn off the Bluetooth on that cell phone. If the Bluetooth icon disappears, you can try to connect to the scale again.

8. Why isn't my body fat percentage being measured when I connect with the app?

- You are not standing on the scale with bare feet. Wearing shoes or socks can cause the electrodes to fail to measure body data.
- The scale is not placed on a flat, hard surface. Placing the scale on a carpeted or raised surface will result in deviations in the data.
- Weighing without standing still. After weighing, wait until the weight stabilizes or the display on the scale flashes 3 times before the weighing is complete.
- The soles of the feet are not stepping on the electrodes. The soles of your feet need to touch the electrodes, and stand still to measure.
- Incorrectly setting personal information. When personal information is filled out incorrectly, usually the age is too young or the height is abnormal, which is beyond the measurement range of the body fat scale, the body fat percentage will not be measured, only the weight will. Please check if the basic personal information is incorrect. For example, the user's height was originally 150cm, but 50cm was selected by mistake, resulting in the measurement range being exceeded. After modifying 150cm, the measurement is normal.
- Out of age range. Users younger than 10 years old or older than 99 years old can only measure weight and BMI data, and will not be able to measure body fat percentage and other body data.

WARNING INDICATIONS!

Lo

Low Battery Indication:

The Battery is running low. Please replace it with a new battery.

Err

Over-Load warning:

Remove the weight immediately, or damage can be caused.

Err

Contact Error:


It can occur during body fat, body water, and muscle mass estimates if not standing still. Remember bare feet and maximum contact between feet. Go to Page 3 for reference.

If Error messages continue after following the correct actions above, remove the batteries and then reinsert after 1 minute.

How to Delete a Saved Profile

In the event that two profiles conflict (within +- 6.6 lbs of each other) it is possible to delete one of the profiles by doing the following.

1. Step on the scale. Your weight will appear.
2. The scale will try to pull up a user profile. If there are two or more conflicting profiles, the screen will display the conflicting numbers.
3. Press the "▲" button to select **P1** or press the "▼" button to select **P2**.
4. Press and hold the "SET" button for 3 seconds, and the display will show "dEL".
5. The next time the scale is turned on, all information pertaining to the deleted user will be cleared.



TO ENTER DELETION MODE (FOR SINGLE OR ALL PROFILES):

1. Press and hold the "Unit" button for 5 seconds. The screen will flash "ALL" and then show "dEL".

TO DELETE A SINGLE PROFILE:

1. From the "dEL" screen, use the "▲" or "▼" arrow to select the profile you wish to delete.
2. Press "SET" to choose that profile.
3. The screen will flash "YES" for confirmation. Press "SET" again to confirm the deletion.

TO DELETE ALL PROFILES:

1. From the "dEL" screen, press "SET" while the screen displays "ALL".
2. The screen will flash "YES" for confirmation. Press "SET" again to confirm the deletion.

SAFETY WARNING

- The warning signs and icons shown appear in order for you to use this product safely and correctly, as well as to prevent product damage, risk, and injury to you and others.

PRECAUTION NOTICE

- Do not use the scale on people who have body implants such as a pacemaker, artificial limbs, contraceptive devices, or metal plates or screws. It may cause the device to malfunction or produce an inaccurate result. When in doubt, consult your physician.
- Do not disassemble the scale as incorrect handling may cause injury.
- Do not step on the scale when your body or feet are wet, especially after bathing or showering, to prevent slipping.

IMPORTANT NOTICE TO USERS

- This product is intended for humans ages 18+.
- Make sure to use only the type of battery stated (see Section "Scale Setup & Calibration").
- The "Athlete" mode is applicable only to people 18+ years of age.
- Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or over-hydration due to such factors as alcohol consumption, menstruation, illness, intense exercise, etc.
- Do not use on pregnant women. The result is inaccurate and the effects on the fetus are unknown.
- For body fat and body water estimates, always estimate in bare feet.

Product Specifications

1. Bioelectrical Impedance Analysis (BIA) technology for body fat and total body water percentage, muscle mass, calorie intake estimations, body mass index (BMI) calculation, and weight tracking
2. High precision STRAIN GAUGE technology for weight measurement
3. Athlete (for ages 18-80 years) or Normal Mode selection
4. 8-user memories
5. 3-button operation
6. Auto-on and auto-off functions
7. Power-saving LED readout
8. Low battery indicator: "Lo"
9. Overload indicator: "Err"
10. Capacity: 440 lb or 200 kg
11. Graduation: 0.2lb or 0.1 kg
12. Body fat graduation: 0.1% (in the range from 5% - 55%)
13. Body water graduation: 0.1% (in the range from 5% - 85%)
14. Muscle mass graduation: 0.1% (in the range from 5% - 85%)
15. Age range from 18 to 80 years
16. Height range from 3'3" to 7'2.5" (100 to 220 cm)
17. Body fat range: 5% to 55%
18. Total body water range: 5% to 85%
19. Muscle mass range: 5% - 85%
20. Operates with 3 AAA(1.5V) batteries (included)
21. Product dimension: 11.8"×11.8"×0.91" (Approx.)
22. Gift box dimension: 13.31" x 13.15"×2.0" (Approx.)
23. Accuracy of weight measurement: ±1.1 lb (11 lb ~ 143.3 lb); ±1.8 lb (143.3 lb ~297.6 lb); ±2.26 (297.6 lb ~ 440 lb)
24. Output power for Body Fat Analyzer: <300uA

One (1) Year Limited Warranty

This scale is warranted against defects in materials of workmanship (excluding batteries) for one (1) year for the original purchaser from date of retail purchase. It does not cover damages or wear resulting from accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair. Do not return to retailer. Should this scale require review (or replacement at our option) while under warranty, please contact us for help and troubleshooting. For best service, please take a photo or note the information on the product label, located on the bottom of the scale.

Phone: 1-866-843-3905 (8-5pm MST, M-F)

Email: taylorusa@lifetimebrands.com

Web: taylorusa.com



There are no express warranties except as listed above. This warranty gives you specific legal rights, and you may have other rights which vary from state to state. For additional product information, or warranty information outside the USA, please contact us through www.taylorusa.com.