SCS + SCS-B SmartStand

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[ES] Manual de instrucciones

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Contains Non-Chargeable 3V Coin Battery (CR1225)

- [DE] Enthält nicht wiederaufladbare 3V Münze Batterie (CR1225)
- [ES] Contiene no recargable de la moneda 3V batería (CR1225)
- [FR] Contient non rechargeable 3V pièce Batterie (CR1225)
- [JP] 含まれている非充電式3Vコイン電池 (CR1225)
- [KO] 포함 비 충전식 3V 코인 배터리 (CR1225)
- [PT] Contém não-exigível 3V bateria moeda (CR1225)
- [ZHS] 包含不可充电的3V硬币电池 (CR1225)
- [ZHT] 包含不可充電的3V硬幣電池 (CR1225)

Warnings

- [DE] Warnungen
- [ES] Advertencias
- [FR] Avertissements
- [JP] 警告
- [KO] 경고
- [PT] avisos
- [ZHS] 警告
- [ZHT] 警告

[EN] SmartStand SCS and SCS-B

It is essential that the operating instructions are read before the tool is operated for the first time. Always keep these operating instructions together with the tool. Ensure that the operating instructions are with the tool when it is given to other persons.

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Intended Use

[EN] This product is designed to charge electronic Rainin pipettes, confirm/report calibration state of both electronic and manual Rainin pipettes, and when used with the proprietary Rainin software, provide key asset management records including pipette use history, calibration status, and location/recent location of individual pipettes.

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




















12.2


www.mt.com/legal
 End User License Agreement


[DE] Endbenutzer-Lizenzvereinbarung
 [ES] Acuerdo de licencia de usuario final
 [FR] Contrat de licence de l'utilisateur final
 [JP] ユーザー使用許諾契約
 [KO] 사용자 사용권 계약
 [PT] Contrato de licença de usuário final
 [ZHS] 用户许可协议
 [ZHT] 用戶許可協議

13.1 Connect via Bluetooth

[DE] Anschließen über Bluetooth
 [ES] Conectar a través de Bluetooth
 [FR] Connecter via bluetooth
 [JP] Bluetooth経由で接続します
 [KO] 블루투스를 통해 연결
 [PT] Conectar via Bluetooth
 [ZHS] 通过蓝牙连接
 [ZHT] 通過藍牙連接


13.2


[DE] SmartStand über PC einrichten
[ES] Configuración del SmartStand a través del PC
[FR] Configuration de SmartStand via un PC
[JP] PC経由でSmartstandをセットアップ
[KO] PC를 통한 스마트스탠드(SmartStand) 설치
[PT] Configuração da SmartStand pelo PC
[ZHS] 通过PC设置SmartStand
[ZHT] 透過個人電腦設定SmartStand

13.3 Setup SmartStand via PC

13.4

14.1 Set Next Calibration Date

- [DE]** Das nächste Kalibrierdatum setzen
- [ES]** Conjunto siguiente fecha de calibración
- [FR]** Définir la prochaine date de calibration
- [JP]** 設定次の校正日
- [KO]** 설정 다음 교정 날짜
- [PT]** Definir próxima data de calibração
- [ZHS]** 设置下一个校准日期
- [ZHT]** 設置下一個校準日期

1

Assets

Rainin	SL-200XLS	K10012	B1243	7BF183	1/3/2016	8/9/2016	872	

2

SL-200XLS

Next Calibration 09/28/2016

Save Close

14.2 Set Next QuickCheck Date

- [DE] Das nächste QuickCheck-Datum setzen
- [ES] Conjunto siguiente fecha de QuickCheck
- [FR] Set prochaine date QuickCheck
- [JP] 次QuickCheck日付を設定します
- [KO] 다음 QuickCheck 날짜를 설정합니다
- [PT] Definir próxima data de QuickCheck
- [ZHS] 设置下一个QuickCheck日期
- [ZHT] 設置下一個QuickCheck日期

10:45 AM
9/22/2015

BG136003456

06/06/2014

06/06/2015

12/06/2015

09/04/2015

E4-20XLS
B813000199

159

✓

A screenshot of a spreadsheet interface. At the top left, there is a blue circle containing the number "1". Below it, a grey header bar contains the word "Assets" in a blue box, followed by seven empty circular icons. The main area of the spreadsheet is a grid. A horizontal blue band highlights one row, which contains the following text from left to right: "Rainin", "SL-200XLS", "K10012", "B1243", "7BF183", "1/3/2016", "8/9/2016", and "872". The date "8/9/2016" is circled with two concentric blue circles.

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















[EN] Nearby equipment of all types will emit RF energy and can be a potential source of interference signals with the RFID antennas. Move the SmartStand to another location in case interference is assumed.

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Limited Warranty and Limitation of Liability

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FCC Statements

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense. Any changes or modifications not expressly approved by Mettler-Toledo Rainin could void the user's authority to operate this equipment as well as the product's certifications and warranty. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

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MSIP-CMM-M7R-SS-SCSB1

MSIP-CRM-BGT-BT121

Applicant name/manufacturer 신청인/제조사: Mettler-Toledo Rainin, LLC

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Model Number 모델명: SCS-B

Manufacture/Country of Origin 제조사/제조국: USA

Manufacture Date 제조연월일: 2016/11

www.mt.com/rainin

For more information

[DE] Für mehr Informationen
[ES] Para más información
[FR] Pour plus d'informations
[JP] 詳細については
[KO] 자세한 내용은
[PT] Para maiores informações
[ZHS] 了解更多信息
[ZHT] 了解更多信息

Mettler-Toledo Rainin, LLC
7500 Edgewater Drive
Oakland, California 94621 USA

Subject to technical changes
© 2016 Mettler-Toledo Rainin, LLC
30343930 Rev A