
ENGLISH - USER'S MANUAL


About **SOUNDSORY®**

Soundsory® is a sound based intervention embedded in a proprietary headphone equipped with Bone Conduction to develop sensory integration and establish core foundations. The specific rhythmic music is electronically processed to stimulate body movement and highlight the Attention Mechanisms in the brain as well as all cognitive functions. Soundsory® is a 20-hour program, divided into 2 times 20 days with one 30-minut session per day. It is advisable to perform physical exercises during the session according to the model available online. The program can be repeated for a deeper integration.

Part description

Accessories

Features

- Power on/off - Play/Pause - Backwards/Forwards - Volume up/down
- Bluetooth mode (TF - SD card mode / Program mode
- Battery status on the screen and battery low signal (2 seconds of 4 beeps every 10 minutes)
- Program in TF mode - Part 1 comprising days 1 to 20. Part 2 comprising days 21 to 40.
- Track name and remaining playing time of track
- Automatic switch off without Bluetooth connection for 10 minutes

How to wear it

How to charge

Charge Indication

| Representation  | Status                                   |
|-----------------|------------------------------------------|
| 1 bar           | Low battery. Headphone should be charged |
| 4 flashing bars | Headphone is charging                    |
| 4 fixed bars    | Headphone is fully charged               |

Troubleshooting

| Issue             | Possible Problem                      | Resolution                                                                     |
|-------------------|---------------------------------------|--------------------------------------------------------------------------------|
| Does not turn on  | Check if the battery is fully charged | Charge the battery off                                                         |
| Do not hear sound | Check the volume status               | Switch the device on and on again to reset the default volume                  |
| Sound is too loud | Check the volume level                | Decrease the volume by using volume -- button or switch the device off and on. |

Specifications

|                    |                                                                 |
|--------------------|-----------------------------------------------------------------|
| Model              | TS02 Bluetooth headphone                                        |
| Speaker            | Diameter 40mm Impedance 32 Ohm Sensibility 115db+-3db SPL @1KHz |
| Response frequency | 20-20KHz                                                        |
| Microphone         | Microphone diameter 4.5*2.2 mm Impedance 2.2 K                  |
| Bone Conduction    | Dual suspension BC. Sensibility 100 +/- 3db                     |
| Wireless           | Bluetooth 5 + EDR spec                                          |
| Operation          | 10 meters (33ft)                                                |
| Battery            | Polymer lithium battery 3.7V/1000mAh - 4 hours charging time    |
| Playing time       | 10 hours on SD card                                             |

Storage and maintenance

- Store the headset in a cool and dry place. The working temperature should be 0-45°C (32-113°F). Working in cold/dry environments or highly humid areas will reduce the battery life.
- After storing for a long period, please charge the headset before using again.
- The headset is not water resistant. Do not put the headset into water.
- Clean the headset with a soft dry cloth.
- Do not bring the headset into contact with any sharp object.
- Do not leave the headset in extremely high or low temperatures.

**Note:**  
Design and specifications are subject to change without notice. For more information on the product, please visit [www.soundsory.com](http://www.soundsory.com)

FCC Warning

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. The equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Charges or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Relatión Expresar Statement The device has been evaluated to meet general RF exposure requirement in portable equipment condition without restriction.

Where is the serial code?

Precaution of use, limitations and disclaimer

Soundsory® volume is set by default at a comfortable level for a majority of people. Feel free to adjust the volume for your comfort but do not turn it up too high. Children must use Soundsory® under constant adult supervision and instruction. Do not use Soundsory® more than 30 minutes a day. If you feel pain or discomfort, immediately stop using it. Resume using Soundsory® once you feel well again. Soundsory® is an electronic device, use it carefully. Soundsory® is a pedagogical device. It is not approved by Federal Drug Administration. It cannot be considered as a medical device.

How to contact us?

SOUND FOR LIFE LTD  
Unit 507, 5/F,  
Chinachem Golden Plaza, 777 Mody Road,  
Tsim Sha Tsui East, Kowloon, Hong Kong

**Note:**  
If you need to return a product under warranty, please send a detailed explanation of the problem to [info@soundsory.com](mailto:info@soundsory.com) and wait for the instructions.

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