





Connected Oral Care System

User Guide

Better, At Home, Preventive Oral Care Habits
For Happier, Healthier Smiles

Your New Truthbrush® System



2

Getting Started

Here's what you'll need:

1. A Wi-Fi connection. Have your network name and password ready.
IMPORTANT: Your Truthbrush® system will function properly only when connected to a 2.4GHz Wi-Fi signal.
2. A mobile device (phone, tablet) connected to the same Wi-Fi. You'll use this for set up.
3. The Truthbrush® app. You can get it here. Download it to the device that you will use to finish setup and to track your performance.



Scan to download app



3

The Three Step Setup

- 1 Create your Truthbrush® account. If you are new to Truthbrush®, open the app on your device and follow the wizard to create an account during setup. Got a Truthbrush® account already? Skip to Step 2.
- 2 Set up your Truthbrush® Hub. This is the device that will connect to your Wi-Fi. Got a Truthbrush® Hub already? Skip to Step 3.
- 3 Set up your Truthbrush® Tracker and take the next step toward happier, healthier smiles!

4

Your New Truthbrush® Hub

The Truthbrush® Hub connects your Truthbrush® Tracker to a Wi-Fi internet connection so that data from every brushing event is safely and securely transmitted to our servers in the cloud.

A single Truthbrush® Hub can be connected to multiple trackers. If there are 3 or more of you sharing a bathroom, you only need 1 hub. If you have a 2nd bathroom that is more than 30 feet away, we recommend you add a 2nd Truthbrush® Hub. Your tracker will connect with the hub it is closest to.

The Truthbrush® Hub is required for the Truthbrush® Tracker to function properly.


The Truthbrush® Standard Hub and the Truthbrush® Travel Hub function similarly. Follow the same set up instructions for both devices.

5

Set Up Your Truthbrush® Hub

1 Powering up your Truthbrush® Hub

- a. Choose a location for your hub (Within range of Wi-Fi and within 30 feet of where you brush your teeth).
- b. Plug in your Truthbrush® Hub (Standard or Travel version).



6

Technical Specifications

Operating Voltage	Standard Hub: 110-240v Travel Hub: 5v
Wireless Technology	Wi-Fi 2.4GHz band, 802.11b/g/n, WPA1/WPA2. Bluetooth LE
Working Range	A Truthbrush® Hub should be placed within 30 feet of your Truthbrush® Tracker(s).

Truthbrush® Hubs can be connected to multiple trackers. Software will be updated automatically.

8

Your New Truthbrush® Tracker

Your Truthbrush® Tracker is an innovative, connected oral care device that works like a fitness tracker designed to improve your oral health. Your new Truthbrush® Tracker is effective & easy!

- Track your at home, preventive oral care habits and get helpful, professional tips to improve your oral care routine.
- Use your favorite toothbrush...manual or electric; adult or child; any brand.
- No need to push any buttons. The Truthbrush® Tracker is fully automatic. It starts when picked up and stops when placed back down.

For even better results

- Share your information with your dentist.
- Add Truthbrush® Trackers for the whole family.

Which Truthbrush® Tracker do you have?



9

Set Up Your Truthbrush® Pro Tracker

1 Powering up / Inserting the battery in your Truthbrush® Pro Tracker

- a. Open your Truthbrush® Pro Tracker by placing the Truthbrush® tool (included) on the front face of the tracker aligning the tooth icon on the tracker with the cut out of the tooth icon on the tool.



- b. Unscrew or Rotate counter clockwise to open the tracker.



- c. Insert battery (CR 1620) with the + sign face up




10

- d. Close the tracker, screwing or rotating clockwise until the indicator point on the back of the tool is aligned with the indicator point on the back of the tracker helping make a watertight seal.



- e. The Truthbrush® Pro Tracker is fully automatic and will turn on with movement and turn off after each use. There is no On/Off button to be pushed.



11

