




Instructions for Use


Note: Read the instructions carefully prior to use.

BioHub Overview

Front



Back



Package Contents



KEEP THE BOX. You will use it to return the BioHub.



BioHub



USB Wall Adapter



Power Cable



Instructions for Use

Getting Started

1 Identify a Location

Place the BioHub in an area with:



A dedicated outlet

where others will not be tempted to unplug it.



Good cellular signal.

Use your phone to check your signal strength.




We suggest placing it on your bedside table.

If that's not possible, place the BioHub in a room that you visit often, like your living room.




Keep BioHub away from water

2 Plug In the BioHub



Plug the USB-C into the BioHub.

Insert the smaller end of the power cable into the BioHub.



Plug the USB-A into the wall adapter.

Insert the larger end of the power cable into the wall adapter.



Plug the wall adapter into an outlet.



Make sure the wall adapter fits snug in a standard electrical wall outlet.

3 Confirm Network Connection



Screen will turn on.

After 10 seconds, the BioHub logo will appear on the screen.



BioHub will search for a cellular network.

See our FAQs if BioHub is unable to find and make a connection.



Network found.



Network found. BioHub is ready for use.

When connected, the screen will show the current time.



Note: Keep your BioHub plugged in for the whole monitoring period. If BioHub becomes unplugged, plug it back in to re-establish radio communication.


Using the BioHub

Navigation

Use the arrows to flip between the clock and status screen.



Clock



Status

Status Icons

Bluetooth (left icon)

- This icon indicates that the BioHub's bluetooth antenna is working as expected.
- This icon indicates that the BioHub is currently connected to a device.
Note: This icon only displays occasionally.

Upload in progress (middle icon)

- This icon flashes when data is being uploaded to the BioCloud. **Note:** This icon only displays occasionally.

Cellular Signal Strength (right icon)

- These bars show the cellular signal strength of the BioHub.

Troubleshooting & FAQs

Where should I set up my BioHub?

Place your BioHub in your bedroom on your bedside table. You should be able to view and interact with the BioHub screen. If you are unable to establish a connection, try plugging in your BioHub in a frequently visited room where you have good network reception. **Note:** The BioHub is not waterproof. Please keep away from water.

How do I know my BioHub is working?

If you are able to see the clock and navigate to the status screen, your device is working as expected.

How long do I keep my BioHub for?

Keep your BioHub for your entire prescribed monitoring duration.

The display reads “Not connected. Please unplug the BioHub and plug it back in.”

What do I do?


Unplug the BioHub and plug it back into a standard wall outlet. Please watch the screen and verify that a network connection is reestablished. If this happens multiple times, you may be in an area with poor cellular network reception. Please try again in a different location.




WARNING: The BioHub is not waterproof. Please keep away from water.

Returning the BioHub

At the end of the monitoring period, use the provided, postage-paid box to return the BioHub and accessories.



Unplug the BioHub, power cord, and wall adapter.



Neatly place the items back in the box.



Seal the box and mail it via USPS:

- Outgoing mailbox
- Local USPS office
- USPS collection box

For support, please call us at

888.908.8804

or email us at

support@biointellisense.com