


Meet your toothbrush



Replaceable rotating brush head
Snaps into place, soft bristles pulse and rotate.
⚠️ Only use quip-supplied brush heads. Discontinue use if damaged.

LED power button
Power on/off with one click, hold down to cycle through 3 intensities: deep clean, daily, and gentle.


Pressure sensor
Red light with 4 pulses signals harsh brushing, helping protect teeth and gums.

Stand + base
Holds your quip in an elegant upright stance on any countertop or flat surface.

Travel case
Covers bristles and brush handle to protect on-the-go.

Charger cable
Magnetic end snaps onto the back of your brush and stand for fast, convenient charging.
⚠️ Only use quip-supplied cable. USB power adapter not included.

Get started




1. Charge up

Recycle packaging. Snap the charger cable to the back of the brush, then plug into a power adapter (not included). You can also connect the cable to the stand & mount. LED light turns solid green when fully charged.




3. Ready, set, brush

Wet bristles, add paste, then turn on brush. Gently guide bristles tooth to tooth, spending a few secs on the front, back, & chewing surfaces. **Do not scrub**, let the brush do the work.




2. Set your intensity

Press power, then hold it down to cycle through 3 intensities. Your brush will save this for next time :)



4. Clean thoroughly

Switch mouth quadrants every 30 secs when your brush pulses until the vibrations stop at 2 mins.



Track brushing, customize features, and earn rewards with the quip app from the App Store or Google Play.

Lights to know

Keep an eye on the light! The colors communicate battery life, brushing intensity, harsh pressure, and more.

Brushing intensities

- Deep clean, solid blue
- Daily, solid white
- Gentle, solid lavender

Power + battery

To check charge status, reconnect your brush to the charger cable.

- Fully charged, solid green
- Charging, breathing orange
- Low battery, solid orange
- Brushing too hard, solid red
- Error, blinking red

⚠️ **WARNING** To reduce the risk of burns, electrocution, fire or injury to persons: Use only the charger provided, and specifically an ETL or UL certified USB power supply, class II, 5V output, 5V, 5A max. Discontinue use if this product or the power supply appears damaged in any way. See FCC details at getquip.com/FCC.

Stay fresh & healthy



Refills auto-delivered

Get brush heads and more shipped every 2 months, like dentists recommend, at getquip.com/refills



Earn quip Perks

It pays to subscribe! Earn loyalty credits, exclusive discounts, and more.



Share the love

Refer friends or dentists to earn refill credits at getquip.com/referafriend

INTERIOR

Got a quip question?

We're here to help you take good habits for a spin. Get answers at getquip.com/help

Scan for a how-to video



Follow us

Pro-fresh-ional tips



Round out your routine

Dentists recommend brushing 2 mins (twice a day!), flossing daily, plus regular cleanings with a professional.



Know when to swish


Brushing right after acidic foods or bevs can hurt tooth enamel. Swish water instead, and wait 30 mins to brush (or pop gum or mints).



Don't rinse right away

After brushing, spit your paste, but don't rinse! Wait 30 mins to let the formula work its magic. Repeat AM and PM.

Ready to refill?




1. Remove

Grip handle in one hand, brush head in the other and gently pull up to remove.



2. Refresh

Use a damp cloth to scrub dried paste away from all sides of the handle and the top plastic tip.




3. Replace

Once all parts are clean and dry, slide the new brush head onto the handle until it clicks.

How to mount

Set up your magnetic mirror mount (included with select models) or shop for one at getquip.com.



1. Clean

Thoroughly clean and dry any smooth flat surface (like mirror, glass, or tile) where you'd like to attach your mount.

2. Prepare

Carefully peel off the thin plastic layer from the single-use adhesive on your mount.

3. Mount

Firmly press the mount onto its new home for 10 secs. Your brush will magnetically snap on just above the charging area.

💡 **Pro-tip!** Learn how to easily remove your mount at getquip.com/help

Meet the family

Complete your routine with more oral care faves and fresh supplies shipped to your door.



Refillable Floss String Dispenser



Anticavity Toothpaste




Cordless Water Flosser

Healthy smiles start here

quip

EXTERIOR



Healthy smiles start here

quip

ACCORDIAN FOLD