

Front

Back

Product Size

Outside diameter of lens tube 5.0mm

15mm

14mm

162.8mm

Hazardous substance

Part name	Hazardous material					
	(Pb)	(Hg)	(Cd)	(Cr(VI))	(PBB)	(PBDE)
Printed Circuit Components	X	○	○	○	○	○
Case	○	○	○	○	○	○
Plastic rubber accessory	○	○	○	○	○	○

This form is formulated in accordance with SJ/T 11364
 ○ represent the hazardous material ratio is below GB/T26572-2011
 X represent the hazardous material ratio is over GB/T26572-2011

Troubleshooting

Failure Phenomenon	Solution
No response to power on / No sign of power	Please use the matching charging cable for charging
Camera does not light up	Turn off and restart the earwax removal
App keeps crashing	Uninstall the existing APP, re-sign the QR code for installation
Ear cleaner indicator shows red light	Please use the matching charging cable for charging, the full charge time is about 60 minutes
The product successfully connected to the APP does not display the real-time image normally	Disconnect the phone from the earwax removal Wi-Fi, restart the earwax removal and reconnect the Wi-Fi
WiFi connection is unstable, sometimes disconnected in a few seconds	Restart the earwax removal and reconnect to the APP
The lens is blurred	When cleaning the lens, please wipe it carefully matching double-headed cleaning brush or alcohol swab
Wrong direction when using	Please switch the lightright ear mode at the bottom of the APP interface
Twesizer head cannot be controlled	Please switch to "Twesizer mode" at the bottom of the APP interface and click the APP Twesizer icon to control the opening and closing
If the above operation still does not solve the problem, please contact with bebird customer service. We will further provide you with practical solutions	

Attention

- When cleaning the lens, please wipe it carefully with the equipped alcohol cotton swabs.
- The Ear tip must be screened on properly and make sure the silicone tip is not loose before use, to prevent it from falling off.
- Please observe the surrounding environment before use. Do not use it in an environment that is susceptible to impact.
- Prohibit children from using this product alone to avoid inadvertent injury.
- Please do not put the device inside any form of liquid.
- The battery must be charged once a month for effective use.
- Avoid direct sunlight exposure on the accessories to prevent them from softening.
- When using, the product temperature will rise slightly (from bebird laboratory data, the maximum temperature of the whole machine does not exceed about 45 °C). Please feel free to use.

Product parameter

Brand	bebird
Model	Note5
Pixel	1000W
Lens diameter	3.5mm
Optimum focal length	12-15mm
Depth of field	10-50mm
Gyroscope	3 axis
Battery	230mAh
Battery life	About 60 minutes
Charging time	1h
Image transmission rate	20fps
Working temperature	-10-45°C
Input parameters	DC 5V 0.5A
Rod weight	30g
Working frequency	2.4Ghz
Network standard	IEEE 802.11b/g/n
Image sensor	CMOS

Warranty card

User information	
User name:	Telephone number:
Seller:	Shipping address:

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

RF Exposure Information
 The device has been evaluated to meet general RF exposure requirement.
 The device can be used in portable exposure condition, compliance with exposure requirements.

3-In-1 Multifunctional Earwax Removal