

bebird

See the best in you

Please consult the instructions before using this product

Download App

ios

Search "bebird" in App store

Android

Search "bebird" in Google Play store

or

Scan QR code to download the App

Watch Tutorial Video via this QR code

Auto Connection

1 Turn on the device and open the mobile App.

2 Click "Connect device" on the homepage.

3 After the search is successful, click "Connect".

Press and hold for 2s to power on; the blue light will flash.

Connect device

Connect device

Manual Connection

1 Click "Connect Manually" on the "Connect Device" page.

2 Go to your phone's Settings > Wi-Fi, and connect your device's Wi-Fi.

Connect device

Wi-Fi

Device Overview

Silicone Covered Tips

Camera Lens

Hose

Charging Port

Power Button

When the device is powered on, a blue light will illuminate.

During power-on and connection, the blue indicator light will flash.

Once connected, the blue light will remain lit.

While charging, the indicator light will glow red.

When charging is complete, the red indicator will turn off.

The blue indicator will flash rapidly when the battery is low.

Accessory Installation

Protective Plug :

Place the ear canal protective plug in the appropriate position on the hose to prevent it from being inserted too deeply into the ear canal.

Ear Spoon:

Unscrew the current accessory by twisting it, counterclockwise to loosen and remove it.

Attach the desired accessory by screwing it in clockwise until its tight and secure.

How to use the tweezers components?

Place tweezers compent-A onto the device body firstly, then install and tighen tweezers component-B

Gently pinch the handle of the tweezers with two fingers to close the tweezers.

Warm Reminder

Install the tweezers component correctly before use to avoid its loosening or falling off.

How to use the Hose Lens?

In most cases, you don't need to bend the hose to use it. If necessary, please note that:

Facing the scale, slowly bend the hose with your fingers to the left and right sides of the scale.

Do not bend in other directions, otherwise it may cause damage to the product.

Do not use brute force. Before using the tweezers, the bent hose needs to be straightened again.

How to use the Tacky Head or Sponge Head?

Attach the Tacky Head or Sponge Head to the hose securely, prevent loosening. Ear debris can be Stubbord Ear debris can be removed by Tacky Head. Only earwax can be absorbed by the Sponge Head.

Warm Reminder: Both types of head are disposable.

How to use the Dropper?

1 Insert the nozzle into the liquid and gently press the air bag to suck up the liquid.

2 Install the dropper on the device securely

3 With the help of a camera, liquid can be precisely dripped into the ear.

Warm Reminder:

1. Please avoid immersing the lens in liquid, otherwise it may cause the lens to become blurred and need to be wiped before continuing to use.

2. This dropper is for use with ear drops only (this product does not contain ear drops). If used with other liquids, please please seek professional advice.

Accessory List 1

Travel Box

Main Device

Inner Tray

Tacky Head

Plastic Ear Tip (hard)

Metal, Silicone Wrapped Ear Tip (Moderate)

Double-color Q-elastic Ear Tip (soft)

Smooth round-tipped metal ear spoon

Dome Ear Spoon (for practice use)

Ear Massager Tip

Sponge Head

FCC Warning

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.
Increase the separation between the equipment and receiver.
Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment The device can be used in portable exposure condition without restriction.

FCC ID: 2AU63-ULTRA

How to use the Ear/Nasal/Throat Observation tools?

Put on the ear canal speculum tightly, and slowly go into the ear.

Put on the tongue depressor tightly to prevent it from loosening.

Put on the nasal speculum tightly, and slowly go into the nose.

*Note: Please make sure the tool is securely attached. If feel uncomfortable, please do not continue.

Troubleshooting

Symptoms	Solution
Camera light off	Check if the brightness setting in the settings is set to the lowest setting.
Lens blurry	Gently wipe the lens surface with the included lense brush.
App crashes	Uninstall the existing app and reinstall.
No response on startup	Please charge it with the included charging cable.
Indicator light flashes rapidly	Battery is low, please recharge the battery in time.
The product fails to display live video after successfully connecting to the app	Check if the app permissions are disabled.
Wi-Fi connection is unstable, sometimes disconnecting within seconds	Restart the device and reconnect to the app.
Indicator light and lens light flash simultaneously	Flashing indicates a firmware upgrade to be done. Firmware upgrade will take about 50mins, please charge the device for 30 mins to avoid power cut during the upgrade. Once upgrade is done, no more flashing, and ready for normal use.
Incorrect orientation during use	Switch left/right ear mode on the bottom of the app.

If the issue persists after following the above operations, please contact our customer service, and we will provide you with further effective solutions.

Product Specifications

Brand: bebird

Name: Bebird EarSight Ultra

Usage time: Up to 80 minutes

Image sensor: CMOS

Optimum Focal Length: 12.5mm

Intelligent orientation: Gimbaled Intelligent Gyroscope

Lens grade: 4P* Slim Lens, IP67 Waterproof,

Operating Environment: -10°C to 50°C

Image Frame Rate: >= 25fps

Input Current: 5V = 0.5A

Warnings

Cleaning: Wipe the lens with a lens brush.

Assembly: Ensure all device accessories are securely attached to prevent them from falling off during use.

Operating Environment: Use the device in a suitable, private environment. Avoid using it in crowded spaces.

Battery Maintenance: The device includes a built-in rechargeable battery. If unused for extended periods, charge it once a month to maintain battery health.

Supervision for Children: Children should not use this device unsupervised to avoid accidental injury.

Sunlight Exposure: Keep the main device and accessories out of direct sunlight to prevent softening or damage.

Water Exposure: Do not submerge the device in any liquids.

Temperature Notice: During use, the device may warm slightly due to the LED lights; this is normal and is capped at approximately 45°C (113°F).

Intended Use for Tweezers: Use the tweezers only to remove earwax or other foreign objects from the ear canal.

The inside of the ear canal is relatively fragile. To avoid damaging the ear canal, please do not use this product to scrape the ear canal vigorously.

Warranty Service

Customer information

Customer Name:

Customer Telephone:

Retailer Name:

Customer Address: