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

MENU

- 1 Who Should Use Femometer
- 2 How to Use Thermometer
- 3 Alerts & Sounds
- 4 About the Display
- 5 Troubleshooting
- 6 Product Maintenance & Care
- 7 Technical Specifications
- 8 Product Functions & Features
- 9 Safety State ment
- 10 Disclaimer
- 11 Copyright Statement
- 12 Warranty
- 13 Symbol Key
- 14 Record of Purchase

2. Pair Femometer with Your Smartphone

4. Replace the Battery

4. About the Display

1. Remove the upper cover of the thermometer. 18888 will be shown on the screen. Measurement will start when the LED backlight is on and you hear a beep.	Symptoms	Possible reasons	Solutions
Lo	While measuring, the temperature is lower than 32.00°C (90.00°F).	If the device keeps showing Lo during measurement, please put the cover back, wait for one minute and measure again.	Lo
Hi	While measuring, the temperature is higher than 42.00°C (107.60°F).	Put the cover back, wait for one minute and measure again.	Hi
Err	Hardware failure	Remove the battery and wait for one minute. Put the battery back. If the message comes up again, please contact the retailer for solutions.	Err
0	Lit up and flashing	Please change the battery (CR2032).	0

6. PRODUCT MAINTENANCE & CARE

1. Battery care
 - Replace the battery when instructed via notification in the Femometer App.
 - Always dispose of used batteries according to local laws and regulations.
2. Cleaning & Disinfection
 - Clean the tip of the thermometer with a solution of 75% rubbing alcohol and 25% water before and after every use.
 - Never submerge the device in water – this can damage the device and result in erroneous data.
 - Do not use cleaning agents other than alcohol and water to clean the device as it may damage or decrease the lifetime of the product and/or present safety risks.

8. PRODUCT FUNCTIONS & FEATURES

1. Basal body temperature measurement & Bluetooth
 - The thermometer measures basal body temperature to 1% degrees Fahrenheit accuracy. Your basal body temperature is automatically uploaded to your phone via Bluetooth, which is fast and energy-efficient, and represented via temperature Charts and graphs to track and make predictions about your menstruacycle.
2. Data stored in the cloud
 - All your data is stored in the cloud, so you don't have to worry about manually uploading or losing your information.You can also sync the thermometer to different devices.
3. Share date with friends
 - Share your data easily with your partner or doctor directly from the Femometer App via Twitter or Facebook.

11. COPYRIGHT STATEMENT

– Hangzhou Bangtang Network Technology Co.,Ltd (Bangtang) reserves all rights, and retains the final right to interpret and revise this handbook and the statements contained therein. Copyright of this manual belongs to Hangzhou Bangtang Network Technology Co., Ltd (Bangtang). Without prior written permission, no company or individual is allowed to reproduce, copy or translate this many or any information contained therein.

– Hangzhou Bangtang Network Technology Co.,Ltd (Bangtang), Floor 12,92 Yueda Street, Binjiang District, Hangzhou, China

– Manufacturer/Hersteller: Joytash Healthcare Co., Ltd. No.365, Wuzhou Road, Yuhang Economic Development Zone, Hangzhou city, 311100 Zhejiang, China

12. WARRANTY

– This product is warranted from manufacturing defects for one year from the date of retail purchase. It does not cover damages or wear resulting from accident, misuse, abuse, commercial use or unauthorized adjustment and repair of product. Please direct all returns to the place of original purchase and retain your original receipt, as you may be asked to provide a copy as proof of purchase. To find the customer service menu, please visit our company website at www. bongmi.com

FCC statements:

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: The manufacturer is not responsible for any radio or TV interference caused by unauthorized modifications or changes to this equipment. Such modifications or changes could void the user's authority to operate the equipment.

The device has been evaluated to meet general RF exposure requirement, The device can be used in portable exposure condition without restriction

Federal Communication Commission (FCC) Radiation Exposure Statement

Power is so low that no RF exposure calculation is needed.

1. WHO SHOULD USE FEMOMETER

This product should be used by women to measure basal body temperature. Patient is the intended operator.

2. HOW TO USE FEMOMETER

1. Download App
 - Search "Femometer" on Apple Store or Google play to download.
 - Available on the App Store
 - GET IT ON Google Play

2. EN

3. Take Your Temperature

4. EN

5. Upload Date to the App

- Step 1: Once your temperature has been taken, your phone's Bluetooth must be on, the data can be transmitted from the thermometer to your phone.
- Step 2: The data from the thermometer will sync automatically to the Femometer App. If the data does not sync, remove the cap from the thermometer, and tap the scale ring in the Femometer App to transmit the data.
- When you remove the thermometer from the case, a single beep indicates that it is on.
- When you close the case, a single beep indicates that the thermometer is off.
- When taking your temperature, 3 consecutive beeps indicate that the measurement is complete.
- A beep lasting longer than 2 seconds indicates an error. Please close the case and reopen the thermometer until the measurement is complete.

6. EN

5. TROUBLESHOOTING

Symptoms	Possible reasons	Solutions
I never heard the end tone while the thermometer is under my tongue.	Thermometer is off. Measurement time is too short. Button battery has run out.	Close the cover and reopen it to restart the thermometer. Keep measuring for 1-3 minutes. Replace the button battery.
The data did not sync to the app after my measurement.	The Bluetooth is not open. The phone is too far to sync.	Open the Bluetooth on the phone. Put the phone close to the thermometer.
There is no sound when opening the thermometer.	Button battery has run out or the thermometer is broken.	Replace the button battery and if the thermometer still makes no sound, the thermometer is broken.

8. EN

7. TECHINCAL SPECIFICATIONS

Product name Vinca II	Measurement 99.0°F-107.0°F	Receiver Apple devices: iPhone 5 or above with version at least iOS 8.0, iPhone 4/4s excluded.
Model FM-VC-201	Accuracy ±0.09°F(95.0°F-100.4°F)	Android devices: 4.3 or above and supports Bluetooth 4.0/4.1
Size 120.2mm*24.4mm*16.4mm	Weight 24.8g (without battery)	Battery endurance 3-6 months
Weight 24.8g (without battery)	Storage conditions Temperature: -40~131°F	Operating conditions Temperature: 41°F-104°F
Transmission distance Up to 16 feet	Humidity: <93%	Humidity: <85%

10. DISCLAIMER

- This manual is written in accordance with available information and contents are subject to change without notice. Bangtang has made the best effort to guarantee the manual's accuracy and reliability. However, the company does not assume liability for loss or damage caused by omissions, inaccuracies or typographical errors in this manual.
- If there are discrepancies in this manual due to software upgrades, instructions and function dictated by the software should prevail. If improvements are made to the product or the product technology changes, please visit Bangtang's technical support website for further information.
- All symbols and logos associated with Bluetooth, including the Bluetooth word mark, figure mark, combination mark and Bluetooth Smart and Bluetooth ready marks, are all trademarks owned by the Bluetooth SIG. Any use of such marks by Bangtang is under license. Other trademarks and trade names used belong to their respective owners.

12. EN

13. SYMBOL KEY

	Device type is BF
	Serial number
	Read user manual before use
	Manufacturer
	No littering
	Recyclables

14. EN

14. PURCHASING RECORD

Model Purchased:	
Name:	
Serial Number:	
Phone Number:	
Date of Purchase:	
E-Mail Address:	

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▲ Remove the battery when not in use for a long time.

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