


SKG

Waist Massager
GS500




USER MANUAL


开机体验吧

CN


腰带正面




腰带背面



主机



指示灯(电源指示灯) 充满电后, 4个灯全亮; 白色指示灯数量代表剩余电量;
红光照射灯 当指示灯只剩余一格显示红色时, 代表电量低, 建议及时充电;



(线框图仅作参考示意参考, 外观以实物为准)

操作说明

CN

- 第一步: 收到产品, 下载SKG健康APP操控(参考APP智能操控步骤), 如果购买了遥控器, 先将遥控器配对后使用(遥控器使用方法请参照遥控器包装或说明书);
- 第二步: 长按开机键1.5秒开机, 待电极片温热后, 把按摩仪佩戴在腰部并调紧腰带 使电极片按摩头紧贴皮肤(如果皮肤干燥, 则请湿润皮肤后再佩戴);
- 第三步: 打开手机APP, 进行力度、模式、热敷档位等更多智能操控;
- 第四步: 按摩结束后, 产品自动关机。取下腰带放到安全干净的地方存放。

产品基本参数

CN


产品名称	腰部按摩仪		
型 号	S3757AA	额定充电电压	5V ㏑
电池容量	2500mAh/9.25Wh	额定输入功率	10W

充电/电源说明



CN

序号	灯状态	状态描述
1	白灯常亮	剩余电量
2	红灯常亮	低电量
3	红灯闪烁	极限电量
4	白灯循环递增亮起	充电中
5	白灯(满格)常亮	充满电

Front of belt



Back of belt



(The wireframe is for illustration purposes only and may differ in appearance from the actual product)

Step 1: After receiving the product, download the Future Wear APP (refer to the smart control procedures on APP). If you have purchased the remote control, you can also use the it to control the belt (please refer to the package insert of the remote control for how to use);

Step 2: Press and hold the On/Off button for 1.5 seconds to turn it on. After the electrode pads are warm, wear the massage belt around the waist, and then tighten the belt, so that the massage head electrode pads can be attached tightly on the skin;

Step 3: Open the mobile APP or use the remote control to control the intensity, mode and hot compress intensity level;

Step 4: After massage, the product automatically turns off. Remove the belt and store in a clean and safe place.

Product basic parameters

Product name	Waist Massager		
Model	GS500	Rated charging voltage	5V $\equiv$
Battery capacity	2500mAh/9.25Wh	Rated input power	10W

Charging instructions/power supply

No.	Indicator Status	Status Description
1	White light steady On	Remaining battery charge
2	Red light steady On	Low battery
3	Flashing red light	Extremely low battery
4	White lights turn On in cycles	Charging
5	White light steady On	Fully charged

Please support the recycling of rechargeable batteries (lithium-ion batteries)

 Danger	This product's battery is a dedicated rechargeable battery that should not be used in other products. Do not charge the battery once it has been removed. - Do not throw into the fire or expose to excessive heat. - Do not puncture with nails, subject to heavy impact, or disassemble - Do not directly connect the battery's positive and negative terminals with metal - Do not charge, use, or store in high-temperature locations such as near a stove or under direct sunlight - All of the abovementioned scenarios can result in overheating, fire, or bursting.
 Warning	Place the removed rechargeable battery out of the reach of infants to prevent accidental ingestion - If swallowed, seek medical attention immediately. - If battery leakage occurs, do not touch the liquid with bare hands. - Follow the safety instructions below: - There is the risk of blindness if battery fluid contacts the eyes. Do not rub the eye with your hands. Rinse with clean water immediately and seek medical attention. If battery liquid gets on the body or clothing, it may cause skin inflammation or injury. Rinse with clean water, and if discomfort persists, seek medical attention mediately. - Disassembly is only permitted when the product is being discarded. Failure to comply may result in fire, electric shock, injury, or other accidents. Please contact the retailer where you purchased the product for repairs (except battery replacement).

The lithium-ion battery used in this product is a valuable recyclable resource. Please remove the lithium-ion battery prior to disposal and hand over it to a professional organization capable of recycling and processing rechargeable batteries.

- The battery must be removed from the device before discarding the appliance.
- The device must be turned off before removing the battery.
- The battery should be disposed of safely.

Removal Steps:

- Open the shell using a tool.
- Unscrew the screws at the circuit board.
- Disconnect the battery terminal plug.
- Remove the battery.

Cautious

To mitigate its impacts and possible harm, do not use the product with the following medical electronic devices:

- Implanted medical electronic devices such as pacemakers.
- Life-sustaining medical electronic devices such as artificial hearts and lungs.
- Medical electronic devices such as electrocardiograms.
- Use is prohibited in flammable and explosive areas such as gas stations, high-dust workshops, and gas stations.
- This product is not equipment and does not diagnose or treat waist disorders. Patients with waist ailments should consult healthcare professionals to avoid delaying treatment.
- Do not use this product while driving, climbing, or engaging in other dangerous activities, as it may lead to accidents.
- Do not use near the heart, above the neck, head, mouth, genital areas, or on skin disorders, and avoid massaging both soles of the feet simultaneously.
- The hot compress period for this product must not exceed 15 minutes (please confirm the time), as there is a risk of low-temperature burns. Do not use while asleep or in an unconscious state, and users who are sensitive to heat should use a lower heat setting and strictly control the use time.
- Individuals with circulatory disorders such as cerebral vascular diseases, inflamed or swollen waist skin, suppuration, or waist -related ailments should not use this product, as it may worsen symptoms.
- Keep this product out of reach of infants, children, or individuals who are entirely or partially unconscious. Do not let them use it, as it may lead to accidents.
- Individuals over 70 years old are prohibited from using this product. Please select a lower setting and strictly control the usage duration since those over 65 have more delicate skin and less sensitive nerve endings.
- Gradually increase the massage setting according to personal comfort; otherwise, the discomfort from a higher setting might cause muscle strain.

Warning

- This is a personal care massage product and cannot be used to substitute medicines or medical devices.
- Children, pregnant women, and women on their menstrual period are prohibited from using this product.
- Contains a lithium battery. Do not place it in a high-temperature environment.
- Use this product under the supervision of a doctor if you have a sensitive skin, facial allergic dermatitis, or facial areas containing cosmetic materials.
- Users allergic to metal or heat must take precautions.
- This product has a heating function; please be cautious if you are insensitive to heat.

- If discomfort occurs during use, please stop using immediately.
- Do not disassemble this product yourself.
- Do not attempt to replace the battery yourself.
- If you have a waist injury, consult healthcare professionals before using it.
- The device has a heating surface, so those who are insensitive to heat must take caution while using it.
- Use within the normal operating environment: Temperature: 0°C to 40°C; Relative Humidity: <80%RH, no condensation; Atmospheric Pressure: 50Kpa to 106Kpa.
- Please choose normal transportation and storage conditions: Temperature Range: -20°C to +55°C; Relative Humidity Range: ≤95%RH, no condensation; Atmospheric Pressure: 50Kpa to 106Kpa.
- This product contains a lithium battery. Do not disassemble, modify, immerse in water, or expose to open flame, as it may result in damage or explosion.
- Do not press on the heart when using, as this may cause an accident or physical discomfort.
- Do not use this product with medical devices or with more than two of this product simultaneously, as it may cause nausea or physical discomfort.
- Do not use for purposes other than those specified in the instruction manual, as it may cause accidents.
- Do not disassemble, repair, or modify this product yourself, as it may lead to a malfunction or an electric shock accident.
- Do not use to massage parts of the body other than the waist.
- If a patient simultaneously uses high-frequency surgical equipment and this massager, it may cause burns at the massager's electrode location and may also damage this massager.
- If this product is used near shortwave or microwave therapy equipment (within 1 meter), the output of this massager may be unstable.
- Electrode pads should avoid contact with wounds and scars during massage.
- While massaging, the electrode pads must not be placed near the heart area of the chest, as this could increase the risk of cardiac fibrillation.

In the following situations, please consult with a physician before use, otherwise accidents or discomfort may occur:

- Patients currently under medical treatment and experiencing abnormal physical symptoms.
- Patients with infectious diseases.
- Patients allergic to stainless steel electrode pads.
- Patients with heart or brain nerve abnormalities.
- Patients with sensory impairments or skin abnormalities.
- Patients with a body temperature above 38 degrees Celsius (during a fever).
- Patients taking medications.
- Patients with abnormal blood pressure.
- Patients who are unable to express their own will freely.
- Patients with perception disorders due to severe peripheral circulation problems caused by conditions such as diabetes.

The following situations prohibit use:

- Patients with malignant tumors.
- Pregnant women and women in childbirth.
- Patients with acute suppurative inflammation.
- Individuals with severe cognitive impairments.
- Persons with spastic paralysis.
- Individuals with abnormal nervous stress.
- Individuals in a spasm phase following a cerebrovascular accident.
- Individuals with unclear consciousness, and limb and joint contracture deformities.
- Patients with a tendency to bleed.
- Users of cardiac pacemakers.
- Eye care
- Patients with metal implants.
- Individuals with dotting disorders.
- Near the human heart.

⚠ Note

- Please read the instruction manual thoroughly before using this product.
- If you feel physical discomfort due to abnormalities in this product, please stop using it immediately. Consult a doctor if you are experiencing unusual sensations in your body or skin, and follow medical advice.
- Do not massage with electronic devices, as it may cause malfunctions (such as a disruption in the electronic watch's timing).
- If a fault occurs and the device cannot operate, immediately disconnect the power and remove it; otherwise, it may cause overheating or short circuits, resulting in an accident. Contact customer service for after-sales consultation.
- Do not plug or unplug the charging cable with wet hands, as it may cause electric shock or accidents.
- Do not use while sleeping, as it may cause the device to malfunction and the electrode pads to move to other parts of the body, causing discomfort.
- Please wear this product correctly for the cervical spine area, placing it symmetrically on both sides to avoid accidents or discomfort.
- When disconnecting the charging cable from the device, hold the plug rather than the cable, since pulling may cause an electric shock or damage the product.
- Charge the product according to the instruction manual the first time you use it or if the battery has been drained from long-term storage to ensure it works properly.
- It is normal for the product's surface to get heated after extended use.
- When the product is working, it will moderately heat the waist area; do not massage for an extended period of time. If discomfort is felt, take a break before using again.
- Do not forcefully fold the product when using or storing it to prevent damaging its components.
- During the folding process of the product, avoid hair or skin contact to prevent pinching.

Please do not use or store this product in the following locations:

- | | |
|---|-------------------------------|
| In places exposed to direct sunlight | In places with water |
| In strong electromagnetic fields | Near fireworks |
| In hot and humid locations | In places with excessive dust |
| In locations susceptible to vibration and shock | |



Recommendation

- Prior to use, please ensure the main host buttons are working properly.
- If using the product for the first time or after the battery has been depleted from long-term storage, charge the product according to the instruction manual to ensure it works properly.
- Please dispose of this product in line with applicable city environmental protection standards to prevent environmental pollution.
- Recommended for use by individuals aged 16 and above.



Caution:

- If you experience skin irritation when using the body massager, please refrain from using it, and wait two to three days to see if your symptoms ease. If symptoms persist or worsen, please consult a physician.
- If you have eczema, allergies or asthma, you may be more likely to experience a skin irritation or allergy from using the massager.
- * The appliance it must only be supplied at safety extra low voltage corresponding to the marking on the appliance.
- * This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- * Children shall not play with the appliance.
- * Cleaning and user maintenance shall not be made by children without supervision.

This device complies with part 15 of the FCC Rules. Operation is subject to the condition that this device does not cause harmful interference (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications not expressly approved by the party responsible for compliance

could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

ISED STATEMENT:

- English: "

This device complies with Innovation, Science and Economic Development Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of the device."

- French: "

Cet appareil est conforme aux normes RSS exemptées de licence d'innovation, sciences et développement économique Canada. Le fonctionnement est soumis aux deux conditions suivantes:

- (1) l'appareil ne doit pas produire de brouillage, et
- (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement."

RF warning statement

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

L'appareil a été évalué pour répondre aux exigences générales en matière d'exposition aux RF. L'appareil peut être utilisé en condition d'exposition portable sans restriction.



This device was tested for typical body-worn operations. To comply with RF exposure requirements, a minimum separation distance of 0cm must be maintained between the user's body and the handset, including the antenna. Third-party belt-clips, holsters, and similar accessories used by this device should not contain any metallic components. Body-worn accessories that do not meet these requirements may not comply with RF exposure requirements and should be avoided. Use only the supplied or an approved antenna. This device in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU. All essential radio test suites have been carried out.

1. CAUTION : RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE. DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS

2. The device complies with RF specifications when the device used at 0mm from your body EU Compliance Statement:

SKG Health technology Co.,Ltd. Here by declares that this device is in compliance with essential requirements and other relevant provisions of the Directive 2014/53/EU. A copy of the EU Declaration of Conformity is available online at <https://www.skg.com/pages/skg-declaration-of-conformity>.

 Old electrical appliances must not be disposed of together with the residual waste, but have to be disposed of separately. The disposal at the communal collecting point via private persons is for free. The owner of old appliances is responsible to bring the appliances to these collecting points or to similar collection points. With this little personal effort, you contribute to recycle valuable raw materials and the treatment of toxic substances.



Points de collecte sur www.quefairedesdechets.fr
Privilégiez la réparation ou le don de votre appareil !



SUCCESS COURIER SL

CALLE RIO TORMES NUM.1, PLANTA 1, DERECHA,
OFICINA 3, Fuenlabrada, Madrid, 28947 Spain

successservice2@hotmail.com

OST_EU_20230214000009

