

SKG

Eye Massager
ET-EN

Manual

Product diagram

Operation guidance

• Before using this product, please adjust the elastic band to the appropriate length, and then wear it to the eye area. Properly adjust the position and tightness of the message head according to the wearing situation, and then start the host to experience the message effect.

Operation steps:

- 1) Long press the "power key" of the host for about 2 seconds, and the host starts up (with voice prompt, the prompt light is on). The first boot defaults to the vibrant mode, and the mode can be switched by short pressing the "power key". In the power-on state, press and hold the "power key" for about 2 seconds to turn off the product (with voice prompt, the indicator light turns off).
- 2) When you press the "power key" once, the mode will be switched to one way, and the four modes will be switched in a cycle (with voice prompt: "vibrant mode", "soothing mode", "sleep mode" and "eye exercises mode"). background music will be played when you are in vibrant mode, soothing mode and eye exercises mode. Under different modes, the message techniques are different, and each mode has the function of hot compress.
- 3) The host is timed for 15 minutes (± 30 seconds). When the host is powered on, it will start timing. It will automatically shut down after 15 minutes (± 30 seconds), and there will be a voice prompt with the shutdown (if it is automatically shut down in sleep mode, the shutdown prompt time will not be played).
- 4) When the mobile phone is turned on, the Bluetooth searches "S1511HB". For the first connection, you need to manually click the connection. After pairing and saving, each time you turn on the Bluetooth, it will connect automatically. After the connection is successful, the voice prompts "Bluetooth connected". After the mobile phone disconnects Bluetooth, the host voice prompts "Bluetooth disconnected". When the mobile phone connects the device with Bluetooth and plays mobile music, the music brought by the device stops and plays mobile music.
- 5) When the product is under low battery, it will be accompanied by a low power prompt, and the indicator light will change from blue light to orange light. The product will shut down automatically in case of power shortage. When the machine is shut down, the message mode before the shutdown will be recorded. After the next startup, the machine will return to the message mode before the last shutdown. When you start up for the first time after the default mode (vibrant mode).

Basic product parameters

Device Name	ET-EN	Model	S1511HB
Rated Input Voltage	5V⚡	Rated Input Power	5W
Battery Capacity	1300mAh/3.35Wh		

Note: This product is rechargeable and is suitable for charging with type-C interface. Users need to charge with a standard charger with an output of 5V⚡, 1A.

SKG Health Technologies Co., Ltd.
Website: www.skg.com
Address: 25A Floor Building 3, Donghai R & D Park No. 1005, Daxin South 1st Road, High-tech Zone, Chongqing (Tianjin) Investment Zone, Chongqing, China
If you have any questions, please feel free to send an email to support@sgk.com

Recommendations:

- Please use it in normal environment temperature 0 °C ~ 40 °C; Relative humidity: < 80% RH; no condensation; Atmospheric pressure: 100kPa ~ 109kPa.
- Please select normal transportation and storage conditions: ambient temperature range: +20 °C ~ 65 °C; Relative humidity range: < 85% RH; no condensation; Atmospheric pressure: 100kPa ~ 110kPa.
- It is recommended that you use it no more than twice a day, no more than half an hour each time, according to your personal situation.

Daily maintenance

- When the message device surface is dirty, it can be wiped with a clean soft cloth dipped in a little water. It cannot be wiped with volatile liquids such as alcohol. The product cannot be disassembled and tried to reuse.
- When not in use for a long time, please close the product to avoid dust.
- In the process of cleaning the machine, prevent water or other liquids from entering the shell to avoid short circuit during use.

Notes for Product Scrap

Please recycle the rechargeable battery (lithium ion battery)!

The battery in this product is a dedicated rechargeable battery, which is not intended for other products. Please do not charge the removed rechargeable battery.

- Do not throw it in the fire, or heat it up to high temperature.
- Do not place it with nails, impact it by gravity or decompose it for modification.
- Do not connect the positive and negative poles of the battery directly with metal.
- Do not expose it together with necklaces, bracelets, etc.
- Do not charge, use or put it in the places in high temperature such as fireproof or in the heating car.
- The situations above may cause heating, fire and leakage.

Keep the removed rechargeable battery out of the reach of infants and young children in case that it is swallowed by mistake.

- When the device is mistakenly swallowed, please consult a doctor immediately.
- When the battery fully leaks, do not touch the fluid by hand, and it can be processed in the following way:
- Do not rub your eyes with your hands when battery electrolyte enters your eyes in case of leakage. Rinse your eyes with clean water immediately and seek medical advice. When your body or clothing is contaminated by the battery electrolyte, it may cause skin irritation or injury. Please rinse with water. If you still feel discomfort, seek medical advice immediately.
- Only the unopened products can be decomposed. Otherwise, the accident of fire, electric shock, injuries may be realized. When you need to repair the battery out to replace the battery, please consult the dealer when you purchased this product.

The lithium-ion battery in this product is a valuable resource that can be recycled and reused, which is not allowed to be incinerated. Please hand it over to a professional agency that can recycle the rechargeable batteries when this product is scrapped. Before scrapping this product, remove the battery from this appliance. Before removing the battery, the appliance must be powered off. The battery should be disposed of safely, which can be taken out in the steps below:

1. With professional tools, open the upper and lower cases.
2. Remove the circuit board and the resistor.
3. Take out the battery, cut the battery wire with professional pliers, and then take out the battery.

Matters needing attention

Attention

- Before using this product, please carefully read and abide by the operation instructions in the operation manual. The company is not responsible for the consequences caused by incorrect use or non-compliance with the instructions of this product.
- As we constantly improve our products, the technical specifications are subject to change without notice.
- Do not use this product after facial makeup, which will damage your makeup and pollute the contact surface of the product. At this time, please wipe the contaminated surface with clean water.
- Do not use the eye massager while driving.
- Do not use the eye massager treatment in places with excessive high temperature and humidity.
- Do not use in bathing with shower, sauna/bath and hot springs.
- Do not use when there are symptoms of eye diseases such as glaucoma, cataract, retinal detachment and red eye disease.
- Please avoid using eye massager after drinking.
- Patients who have performed eye surgery and have heart disease, children, pregnant women and women in their critical period are forbidden to use eye massager.
- There is a lithium battery inside. Please do not put it in a high temperature environment.
- For sensitive skin, facial allergic dermatitis and plastic materials on the face, use this product under the guidance of a doctor.
- This product has heating function. Please pay attention when using it for those who are not sensitive to heat.
- If you feel uncomfortable during use, please stop using it immediately.
- Do not get massage in other parts of the body, except the eyes. It is forbidden to disassemble this product by yourself.
- The product has no waterproof function. Please do not wash it or use it in high humidity environment. If the lining is stained, you can wipe it gently with a wet towel.

FC

15.19 Labeling requirements

This device complies with part 15 of the FCC Rules. Operation is subject to the condition that this device does not cause harmful interference.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including that from any licensed radio communication system.

15.21 Extension to use

Any changes or modifications may be made only if they are approved by the party responsible for compliance with the user's authority to operate the equipment.

15.102 Information to the user

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can receive radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

CE

This device was tested for typical body-worn operations. To comply with RF exposure requirements, a minimum separation distance of 10 cm must be maintained between the user's body and the handset, including the antenna. Third-party call-, data-, headset, and similar accessories used with this device should not contain any RF components. Body-worn accessories that do not meet these requirements may not comply with RF exposure requirements and should be avoided. Use only the headset or an approved alternate.

This device is compliant with the essential requirements and other relevant provisions of Directive 1999/5/EC. All essential and other requirements have been complied with.

1. CAUTION: RISK OF OVERHEATING IF BATTERY IS REPLACED BY AN INCORRECT TYPE. DANGER OF USER BURNING ACCORDING TO THE INSTRUCTIONS.

2. The device complies with RF specifications when the device used at 0cm from your body. Please visit <https://www.skg.com> (or confirmed) to find our CE Declaration of Conformity (DoC) file, and printer correct info.

90*100mm 128g C2S, 印刷颜色3252C+Black 风琴4折