




Instagram: defi_smart_official


E-mail: gasian@gasian.co.kr



DEFI Smart Pillbox M

Instruction

What is included



Smart pillbox M

USB-C charging cable

Information

Product name	Smart Pillbox M
Model	GPB-2.0
Power	3.7V Lithium-ion battery
Bluetooth frequency	2.4GHz
Temperature range	-4°F (-20°C)

Safety information & Key features

Safety information ① If the product is not turned on, charge it with the enclosed charger.
② if the alarm does not go off when it's connecting the pillbox, please retry after turn off the application.

Key features ① DEFI Smart Pillbox helps people to take their medicine on time by their own alarm setting.
② Set an alarm on DEFI application and it will connect with a smartphone.
③ DEFI Smart Pillbox needs Bluetooth connection with your smartphone.
④ It extends battery life up to 30 days on a single charge.

Installing the DEFI app



① search and download "DEFI" on Google Play or APP store.



③ After connecting with your smart pillbox, set an alarm.



② Ensure that your mobile's Bluetooth is turned on for the app to connect with your DEFI Smart Pillbox and create your DEFI account.



④ Check the operation by filling up your pillbox (IR sensor) and history (open/closed).

Warranty Form

Product name	Smart Pillbox M
Model	GPB-2.0
Serial number	
Place of purchase	
Purchase date	