





Red Light Therapy Lamp

User Manual



FCC ID: 2BGFSSSTR-LT-002


Package List

Part	Part Name	Part Count
	Red Light Therapy Lamp	× 1
	DC24V-3A Adapter	× 1
	Glasses	× 1
	Wrench	× 1
	Remote control	× 1

Product Parameter

Model	STR-LT-002
Material	Metal+PC
Working Voltage	AC100-240V, 50/60HZ
Wattage	MAX 50W
LEDs QTY	78PCS
Brightness Adjustment	5-Level Brightness
Switch	Touch switch+Remote control
LED Beads	Red 660nm(78pcs)+NIR 850nm(78pcs)
Timer Function	15mins/30mins/45mins/60mins


Remote Control Instructions



Tips

- The remote control has been paired with the lamp before leaving the factory.
- If it fails to operate the lamp, please follow the steps below to re-pair it.
- Plug in and power on the lamp.
 - Immediately press and hold the switch button on the remote control to start the pairing process.
 - The light will flash once to indicate that pairing is successful.

Button Guide



Tips

- Wear eye protection when using your device.
- Ensure the target skin area is fully exposed, as clothing and accessories will impede the penetration of red light.
- Find a suitable location for your device where you can sit or stand for 20 minutes at a time.
- If you have sensitive skin and the target areas look red or feel tight, please shorten your usage time.
- Red light therapy is intended as a form of physical therapy only and is not designed for medical treatment.

Warning:

Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC Statement:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

RF exposure compliance statement:

This device has been evaluated to meet the general RF exposure requirement