


JP

- EN** Instructions for eye massager (P1) **DE** Anleitung für Augenmassagegerät (P2)
- IT** Istruzioni per massaggiatore oculare (P3) **FR** Instructions pour masseur oculaire (P4)
- ES** Instrucciones para masajeador ocular (P5) **HU** Szemmasszírozó használati útmutató (P6)
- CS** Návod k očnímu masážnímu přístroji (P7) **PL** Instrukcja masażera do oczu (P8)
- NL** Handleiding voor oogmassageapparaat (P9) **SV** Bruksanvisning för ögonmassageapparat (P10)
- EL** Οδηγίες για μασάζ ματιών (P11) **DA** Brugsanvisning til øjenmassageapparat (P12)
- JP** アイマッサーの取扱説明書 (P13)

Preparation

1. Turn on the machine to check whether the power of the massage is sufficient
2. If the power is insufficient, please use DC5V stable power supply configured by our company to charge until the power is sufficient

Charge

1. Find the charging port on the left side of the massager body
2. Insert the output plug of the USB charging line into the charging socket of the product, and then connect the USB end to the charging adaptor of the home mobile phone.
3. During the charging process, the red charging indicator is always on, and the red indicator is off after the charging is completed (the blue indicator is on after charging is completed for voice model).
4. During the charging process, then machine can not work (self-protection function, automatic shut off for the machine during charging).
5. Charge for 3 hours.

Operation Method

Press the switch for 2 seconds to turn on the machine. The default mode is automatic.

When press key shortly, the key will enter the next mode. The massage mode will be cycled in turn. Press the button for 2 seconds to switch off the machine.

Mode cycle with one button control

Automatic: Pneumatic massage + vibration + hot compress + music

Mode 1: Pneumatic massage + vibration + phone music

Mode 2: Pneumatic massage + hot compress + music

Mode 3: Pneumatic massage + music

Mode 4: Hot compress + music

The default working time of the machine is 15 minutes after it switches on. If you want to extend the using time, you need to switch on the machine again (it is recommended that the using time should not exceed 15 minutes for each time)

Music control key +

Press this key for a long time to control the volume and sound to increase. Press this key for a short time to switch to the next music

Music control key -

Press and hold this key for a long time to control the volume and sound to be smaller. Press this key for a short time to switch into the previous music.

Call answer key

Connect, press this key when call comes to answer.

How to connect your phone

Turn on the device, locate the device, and click the link. After the link is successful, you can play music from your phone and answer incoming calls

1. Operating Frequency: 2402MHz ~ 2480MHz
2. Maximum Radio-Frequency Power: 2.28 dBm

Matters needing attention

- Do not use the device after drinking alcohol
- Do not use this device while wearing contact lenses
- Do not use this product in the environment with temperature over 40 °C and excessive humidity.
- It is forbidden to wash the product with water. Please use cotton cloth soaked with alcohol for cleaning.
- Do not disassemble or attempt to repair by yourself.

Trouble and Solution

Can't switch on the power source?

- low battery, please arrange charging.
- incorrect connection for power adaptor, please connect it properly.
- broken inner lithium battery.

if can't charge, please change new battery.

Specification

Model: JB-018

Rated voltage: 5V-1

AP Product Power: 5W

Material: ABS

Reliable, safe, comfortable

Standard: GB4706.1-2005

FCC Warning

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE 1: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

NOTE 2: Any changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

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