
Homecare Device

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Manual

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# Homecare Device

---

Thank you for purchasing a product from JLU Co., Ltd.

This product is a beauty device, not a medical device.

Please note that the user experience and perceived effects may vary by individual. If any abnormal symptoms occur, discontinue use immediately.

## Table of Contents

|                                                                      |       |    |
|----------------------------------------------------------------------|-------|----|
| <b>01. Caution</b>                                                   | ----- | 13 |
| <b>02. Product Introduction</b>                                      |       |    |
| - Components & Name                                                  | ----- | 14 |
| - Product Specifications & Charging Information & Device Care Method |       | 15 |
| <b>03. Usage Method</b>                                              |       |    |
| - Usage Method                                                       | ----- | 17 |
| <b>04. Warranty &amp; A/S</b>                                        |       |    |
| - Warranty & A/S                                                     | ----- | 20 |
| - Product Warranty                                                   | ----- | 21 |


## ※ CAUTION

1. This product is a beauty device, not a medical device. Please refrain from using it for any purpose other than beauty care.
2. Remove the Protective film before using the device.
3. Before using the product, please make sure it is functioning properly.
4. Please keep out of reach of children.
5. Not intended for use by individuals with reduced physical-sensory, or mental capabilities (including children) or those lacking experience and knowledge, unless supervised or given instructions on device use by a responsible person.
6. Avoid use on pregnant-breastfeeding women.
7. Remove metal accessories such as necklaces-earrings etc., before use.
8. Avoid use on areas with skin diseases such as wounds-infections-boils.
9. It is recommended to consult with a medical professional when using with dental treatments such as implants-nerve treatments-orthodontics.
10. It is recommended to consult with a medical professional before undergoing any cosmetic procedures.
11. Avoid use on people who require health precautions, such as having pacemaker-device implants-history of severe illness-outpatient treatment.
12. This product is a skin care device intended for use on the face only. Avoid using it on any area other than the face.
13. Avoid use on pupils(eyes), around the eyes and where the thyroid gland passes.
14. During use, make sure the device is in full contact with the skin.  
Do not leave any space between the device and the skin.
15. Use carefully, as prolonged use on one area may cause skin irritation.
16. Set the appropriate power based on the skin sensitivity.
17. We are not responsible for any side effects that may occur as a result of not following our recommendations.


## ※ Product Care

1. If a detachable power supply is provided, only use the device with the power supply provided with it.
2. Do not modify or disassemble-repair the product. If the user modifies-disassembles-repairs the product, free A/S service cannot be provided within the warranty period, therefore be sure to use the designated A/S center.
3. Do not touch with wet hands, do not wash or immerse the product in water. Be careful as it may cause electric shock-malfunction.
4. Do not store the product in humid environments (bathrooms-laundry rooms-saunas etc), as this may cause malfunction. Store it in a dry place, away from moisture.
5. Avoid dropping the product or subjecting it .
6. Disinfect the head(contact area) with alcohol after use. Use a soft cloth for disinfection, as rough cloths or strong wiping can cause scratches or damage.
7. Use wet tissues or a soft cloth to clean the main body, as using alcohol may cause discoloration-contamination.
8. For detailed product care-A/S services, follow the manufacturer's warranty regulations.

Product Components


Main Body

Charging Cable  
C-TYPE

Product Name


Product Specifications

- \*This product's specifications may be subject to change without separate notice for performance improvement.
- \*The specifications of this product may vary depending on the measurement method, such as the measuring instruments and measurement locations.

|                     |                                      |                     |                   |
|---------------------|--------------------------------------|---------------------|-------------------|
| Model Name          | JTMC200B                             | Product Name        | Skin care device  |
| Product Size        | About W47 x H220 x D47 mm            | Product Weight      | About 299g        |
| Input               | DC 5V 1A                             | Battery             | 800mAh            |
| Charging Cable      | C-TYPE                               | Charging Time       | About 120 minutes |
| Usage Time          | Around 2~3,000 uses on a full charge | Manufacturer/Seller | JLU Co., Ltd.     |
| Date of Manufacture | Separate Notation                    | Customer Center     | +82 51-331-7185   |
| FCC                 | 2BHHE-JTMC200B                       | A/S Center          | +82 10-4505-7989  |

Charging Information


When connecting the charging cable to the charging port on the bottom of the device, push the connector in completely. If not, charging may not proceed smoothly. If the device does not power on, check the connection of the charging port on the bottom of the device.

### Charging Information


**Usable**

**Charger 5V 1A / 5V 2A**


**Not usable**

**Non-standard charger (5V 3A)**

1. Please be sure to use an adapter with the rated voltage (5V 1A / 5V 2A). Adapters that do not meet the specifications such as laptop chargers or ultra-fast chargers must not be used. If a non-compliant adapter is used, charging and usage time may differ from the stated information, and any malfunction or damage resulting from this will be considered consumer negligence. The manufacturer will not be held responsible.
2. Please make sure to use the included charging cable. (or a C-to-A cable)
3. Charging using a power bank is not allowed.
2. Please make sure to use the included charging cable. (or a C-to-A cable)
3. After the device is fully charged, if you turn it on and then connect the charger again, the battery indicator at the top may change from 3 bars to 2 bars.
4. This does not mean the battery is insufficiently charged. The 3 bar level is designed to correspond only to the highest voltage range, so even brief use may cause the indicator to display 2 bars.
5. The battery is still fully charged, and this does not affect the actual usage time.
6. Near full charge or right before complete discharge, the device is designed to charge slowly with a micro-current for battery safety. Because of this, charging may feel slightly longer than expected.
7. Do not operate the device while charging (e.g., pressing buttons), as this may cause malfunction.
8. If the charging port is not properly aligned, or if dust or foreign substances accumulate on the connector, poor contact may occur. Please check carefully before charging.

### Device Care Method

#### Alcohol Disinfection Allowed

(Use a soft cloth for disinfection)

\*Using a rough cloth or wiping the surface too hard may cause scratches or damage. Product damage due to user negligence is not eligible for exchange or refund.

#### Alcohol Prohibited

Use wet tissues or a soft cloth to clean the device body, as using alcohol may cause discoloration-contamination.

|                         |                                                                                                                                                                                                                                                              |
|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Usage Cycle             | Approximately every 14 days, 100 for full-face use is recommended.<br>*Depending on your skin condition, adjust the intensity and frequency of use.                                                                                                          |
| Prohibited Areas of Use | Pupils (eyes), around the eyes, and the area where the thyroid gland passes                                                                                                                                                                                  |
| Precautions             | During use, be sure to use the device while keeping it accurately and firmly attached so that there is no gap between the skin and the device.<br><br>Depending on the individual, slight bruising may occur due to suction pressure. (Temporary phenomenon) |

Temperature Sensor Function Guide

A temperature sensor safety function is applied to prevent device overheating.

※ The above operation is a normal safety control function for device protection and is not a malfunction or defect.

|                                                 |    |                                      |
|-------------------------------------------------|----|--------------------------------------|
| When the measured temperature reaches 40°C      | ▶▶ | Enters Temperature Sensor LOCK state |
| When the measured temperature is 40°C or higher | ▶▶ | Use is restricted                    |
| When the measured temperature is below 40°C     | ▶▶ | Use available                        |

Gyro Sensor that Detects Movements

To prevent malfunction and overheating, the device is equipped with a movement detector gyro sensor. If no movement is detected during operation, the power will automatically turn off for safety. To resume use, please press the power button to restart the device.

|                                    |    |           |
|------------------------------------|----|-----------|
| No movement for more than 1 minute | ▶▶ | Power OFF |
|------------------------------------|----|-----------|

### Caution


#### \*Precautions for Water·Moisture, and the Use of Skincare

1. This product is not waterproof. Please avoid storing or using it in humid environments such as bathrooms, laundry rooms, or saunas.
2. Exposure to water or moisture may cause product malfunction.
3. If skincare are applied directly onto the device, the contents may seep into gaps in the product or into the charging port, which may cause malfunction. Any damage-loss resulting from this will be considered consumer negligence, and the manufacturer will not be responsible.

#### Prevent liquid from entering the device

Head (Contact area)

Gaps in the product

Gaps in the product

Charging port

## 03. Usage Method - Quick Guide

19

### STEP 01

After cleansing, dry your face and use the device.

\*This product is not used together with cosmetics.

### STEP 02

Turn on the power and set the mode, power level, and suction level.

### STEP 03

Please focus on caring for the area below the cheekbones, such as the lower jaw and cheeks.

### STEP 04

Wait until the skin is held by suction pressure and then released again, then move to the next area and continue use.


### Usage TIP!

- \*Do not forcibly detach or move the device while suction is operating.
- \*After suction is released, please move to the next area.
- \*If used continuously on the same area, skin irritation may occur, so please be cautious.
- \*Depending on the individual, slight bruising may occur due to suction pressure. (Temporary phenomenon)

## Operation

|                              |                                                                                             |
|------------------------------|---------------------------------------------------------------------------------------------|
| <b>POWER ON</b>              | Press the power button for more than 3 seconds and the device will turn on with the screen. |
| <b>POWER OFF</b>             | When it turns on, press the power button for more than 3 seconds to turn it off.            |
| <b>Mode Setting</b>          | Press the MODE button briefly to set Mode.                                                  |
| <b>power Level Setting</b>   | Press the power button briefly to set the power level. (Levels 1~5)                         |
| <b>suction level setting</b> | Press the INTEN button briefly to set the suction level. (Levels 1~3)                       |
| <b>Memory Check</b>          | Press and hold the MODE button to check the accumulated usage.                              |
| <b>Pause</b>                 | It is possible to pause with the start/stop button.                                         |

## Display

- ① Currently selected mode
- ② Suction level
- ③ Power level
- ④ Battery level
- ⑤ Temperature status
- ⑥ Temperature LOCK icon

## 04. Product Warranty & A/S

20

This warranty policy applies only within the warranty period and conditions specified by our company.  
For A/S service, please contact the place of purchase or our headquarters A/S center.

### ※ Warranty Period

Within 1 year from the date of production

if the product is registered as genuine and a manufacturing defect of the main device is confirmed during the free maintenance warranty period, free A/S service will be provided.

### ※ Paid A/S

#### 1) Products outside the warranty scope

- Products past the warranty period
- Battery performance degradation and wear of consumables due to normal use

#### 2) Consumables and accessories

- Consumables (head, cartridge, cable, adapter, handpiece, cart, etc.)
- Separately provided accessories (stand, cradle, etc.)

#### 3) Customer negligence or improper management

- Failures caused by dropping-impact-disassembly-modification-repair
- Issues due to unofficial software-firmware changes
- Damage caused by natural disasters (fire, earthquake, flood, etc.)

### ※ Exceptions (A/S is not available)

Lost or stolen

### ※ A/S Shipping Fee

**For services conducted via courier, shipping costs are shared between the seller and the customer.**

- Before repair : Customer bears the cost
- After repair : Seller bears the cost

### ※ A/S Center

- Contact Number: +82 10-4505-7989
- Operating Hours: 9:00 AM – 6:00 PM Closed on Sat/Sun/public holidays (Lunch break 12:00 PM – 1:00 PM)
- Address: JLU Co., Ltd. 411 office, Busan Knowledge Industry Center, 111, Hyoyeol-ro, Buk-gu, Busan, Republic of Korea

# Product Warranty

Thank you for purchasing JLU products.

This product is a genuine product manufactured by the headquarters and shipped after passing strict quality control and inspection processes.

Free service within the warranty period is available if you complete the product registration through the homepage.

A/S Center : +82 10-4505-7989

### [How to Register the Product]

Visit the official website ([www.jlu.kr](http://www.jlu.kr)) → Go to the Product Registration page  
→ Log in → Enter your product serial number

\*S/N : Separate notation on the product

## FCC Compliance Statement

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

## FCC Interference Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates,uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## FCC Caution

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.