



Model: MGE-2302

CE RoHS
Made in China

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

This device and its antenna(s) must not be co-located or operating in conjunction with any other antenna or transmitter.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement.
The device can be used in portable exposure condition without restriction

Important	01
Important Safety Precautions	01
Product Diagram	02
Before use	02
Indicator light description	02
At use	03
Maintenance	04
How to clean	04
Storage	04
Error Diagnosis	05
Product Specifications	06
Warranty & Customer Support	06

Please keep this manual properly for future reference. Before using this appliance, please read all the instructions carefully.

- Do not wear accessories around your eyes when using this product.
- For those with symptoms of eye disease, patients with major eye diseases such as eye surgery, glaucoma, cataract, and retinal detachment, please DO NOT use this product.
- For those who have medical electronic devices implants, such as cardiac pacemakers, or patients with malignant tumors, heart disease, or severe hypertension, please consult your physician before using this product.
- Children and mentally/physically disabled can only use the product with special supervision of a responsible adult.
- If you experience bodily pain and discomfort during use, please stop using it immediately and consult your doctor.
- This product is not a toy and can not be given to children for fun.
- This product is for single-person use. Do not operate simultaneously on more than 1 person each time.
- This product is only for home use. It is not a piece of medical equipment and is not intended for self-diagnosis, treatment, or other commercial purposes.
- Please charge the product using only accompanied type-c cable and please be careful while you charge.
- Please cut off the power when the product exhibit malfunctioning. Do not disassemble or fix the product on your own.
- Do not use the product if you have a fever or body temperature higher than 38°C.
- Please do not disassemble the product on your own to prevent product malfunctioning.
- Do not use a power source other than the rated power to prevent fire, physical injuries, or product damages.
- Do not damage, destroy, over-bend, pull, twist, or tie the type-c cable.
- Do not expose the product in the environment of high-temperature or hot water.
- The product should not be discarded with regular home disposal. Please follow the local guideline for disposing of the product.
- The direction of connecting the product or the type-c cable depends on the product model. Please do not force to connect to the product in case of damaged or shorted connected that may lead to fire or product malfunctions.

Diagram illustrating the front view of the smart glasses with the following components labeled:

- Power button
- Music button
- Type-c Charging port

- (1) Please check the integrity of the machine and accessories
- (2) Please remove the accessories around the eyes during use to avoid affecting the massage area.
- (3) Adjust the elastic headband to the size that fits your head shape.

Battery indicator light			
	Color	Status	Description
Normal working status	Orange	on light	Heating on
	Green	on light	Heating off
Charging status	Red	breathing light	0%-20%
	Green	breathing light	20%-99%
	Green	on light	Fully charged
Low battery reminder	Red	flash light	Low battery

At use

*Without connecting Bluetooth status, massage mode runs for 15 minutes and ends, connecting Bluetooth status, massage mode runs for 15 minutes and ends, music function runs for 30 minutes and ends

How to clean

- (1) The product shell can be cleaned with a clean cloth;
- (2) Stains can be removed by rubbing alcohol on a clean rag;
- (3) For deep dirt, wet the cloth with neutral detergent and gently wipe it off.
- (4) Do not pour water to the product
- (5) Leave it to dry naturally.

Please place the product inside of the storage bag with its original package and store it in a safe place.

Do not put heavy weights on top of the product.

Please take moisture, rodent and pollution prevention measures before storing the product.

For product in more than 3 months of storage, please recharge before use.

Do not store the product in places with direct sunlight, or corrosive gas.

Error	Reason	Solution
Can not turn on.	Insufficient battery.	Please recharge.
Air-pressure strength is weak.	The strap is too loose.	Adjust the length of the strap , make the massager comfortable to wear
Stopped working during operation.	Massage timer is off.	Press the Power button to reactivate the product.
	Battery is ran out.	Recharge the battery.
Cannot connect Bluetooth	Mistakenly touched the music button, Bluetooth is turned off	Long press the power button for 2 seconds to turn off, and the Bluetooth function will be restored after restarting.
Bluetooth playback music sound is low	Mobile phone playback volume is low	Adjusting the phone volume on your phone

Product name	Eye Massager
Product model	MGE-2302
Product power	4W
charging voltage	5V
battery capacity	1200mAh
Product size	215×117×80mm
Net weight	340g

*Note: Products purchased directly from NAIPO benefit from our after-sales service. If the purchases are made from other selling teams, we recommend you contact them for warranty issues.