


Relax Your Eyes Unwind Your Mind.



Charging your Restnature

Charge the device using the provided type-c charging cable and a DC 5V/1A USB charger (not included). Note that QC and PD fast charger plugs cannot charge the product.




Charging Indicator
Flash Red : Low Battery
Light Red : Charging
Light Green : Fully Charged

Low Battery Alarm : When the battery is low, the red light will flash and a warning sound will play. The device will shut off when the battery runs out.

1

Wearing Your Restnature

Remove any eye accessories (such as eyeglasses, contact lenses, or false eyelashes). Put on the device and adjust the buckles for a comfortable fit.




Tips : For a better eye massage experience, you can slightly loosen the head band when you feel the pressure is too strong, or slightly tighten the head band if the massage is weak.

2

Using Your Restnature

Power Button

- Light Orange : **Heating On**
- Light Blue : **Heating Off**
- Light Off : **Power Off**



Music Button

- Light White : **Music On**
- Flash Blue : **Bluetooth Pairing**
- Light Blue : **Bluetooth Connected**
- Light Off : **Music Off**

- Long Press:** Power On / Off
- Short Press:** Switch Mode (Strong Mode → Soft Mode → Heat Only → Sleep Mode)
- Double Press:** Heating On / Off


- Long Press:** Bluetooth On / Off
- Short Press:** Next Track
- Double Press:** Music On / Off

Tips : To maximize your relaxation, keep your eyes closed for 1-3 minutes after using the eye massager before opening your eyes.

3

Maintaining Your Restnature

To keep the product clean, use a soft cloth with mild detergent to wipe it. Avoid using alcohol and ensure that no water or other liquids get inside. Keep the product in the storage bag to protect it from dust.



Notes : Battery life varies by use. If not used for a long time, turn it off and charge it once a month to prolong its lifespan.

4

Need Help ?

Contact us at support@restnature.co
We are more than happy to serve you in 24 hours.

Follow Us !



Orestnature.co



O@restnature_official



O@Restnature

5

FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This device and its antenna(s) must not be co-located or operating in conjunction with any other antenna or transmitter. NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.	However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures: -- Reorient or relocate the receiving antenna. -- Increase the separation between the equipment and receiver. -- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. -- Consult the dealer or an experienced radio/TV technician for help. The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction	EN This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions: 1. This device may not cause interference. 2. This device must accept any interference, including interference that may cause undesired operation of the device. FR L'émetteur/récepteur exempt de licence contenu dans le présent appareil est conforme aux CNR d'Innovation, Sciences et Développement économique Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes: 1. L'appareil ne doit pas produire de brouillage; 2. L'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement. The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction
--	---	---