


BLUETOOTH EYE MASSAGER

Operating instruction



Please read the instruction carefully before use

PLEASE READ THIS INSTRUCTION BEFORE USE

Thank you for buying this eye massager. Please read the operation manual carefully before use, and pay special attention to safety matters, so as to operate and use this product correctly

This manual provides users with a safe and effective way of operation. Please keep it properly after reading it for future reference

ORDERS TO RECORD

1. Safety precautions

2. Product profile

3. Please confirm the commodity content

4. Name and function of each part

5. product usage

6. Cleaning, maintenance and storage

7. Faults and solutions

8. Product specifications

2

2-3

3

4

5

6

7

7

1

SAFETY PRECAUTIONS

If you wear contact lenses, please remove contact lenses before using this product. If there is any dirt or foreign matter in the mask, please wipe it with a wet cloth before using it to avoid accidental eye.

During the charging process, the product is not supported to start up (charging is conducted under the working state, the product starts the self-protection program and stops working automatically).

Do not use the eye massager during outdoor walking or driving.

Do not use eye massager in hot and humid places.

PRODUCT INTRODUCTION

The eye massager uses the microcomputer chip control technology for the entire human nature simulation control, the application of air pressure that a whole new concept to the eye care, under the control of air pressure in the computer to a full range of muscles around the eye pressure, vibration and coupled with the internal function and the function of the far infrared heating, the eye acid bilges, the eye pain, relieve eye fatigue.

1.Using the physical principles of Pressure Therapy and Thermotherapy, it perfectly fits the eye area, deeply massages the eye acupuncture.

2.Shiatsu, Vibration, Air-Press, Hot Compress and Auto—Five Combination for choice, fits every needs and maximize relieve your eye fatigue.

3.15 minutes time setting function and intelligent thermal induction anti-overheated. Once overheated the device will auto stop.

4.Include with Bluetooth function. Enjoy your relaxing time with music. Adjustable volume.

5.Full charge in 2 hours. Long time battery life good for travel and trip.

6.180°folding design together with storage bag. Easy to carry.

7.It's highly recommend to people who suffer from eye fatigue.

2

2. SCOPE OF APPLICATION:

For the eyes

3. CONTRAINDICATIONS:


Disables eye injuries or inflammation.

If there are symptoms of eye diseases can not be used, general implementation of eye surgery, glaucoma, cataract, retinal detachment and other major eye diseases patients should not be used.

For patients with deep stroke, hypertension and hypotension, please consult a professional physician first. If you feel any discomfort during use, please stop using immediately. Please avoid using the eye massager during drinking or physical discomfort

PLEASE CONFIRM THE COMMODITY CONTENT

Eye massage device, USB charging cable, Instruction manual





3

ONE KEY CONTROL LOOP

Integrated mode: Air pressure + Vibration + Hot compress

Clear mode: Air pressure + Hot compress

Sleeping mode: Air pressure

Vitality mode: Air pressure + Vibration

Smart mode: Hot compress +

NAME AND FUNCTION OF EACH PART

1. One-click operation

(1) Long press turn on/off key for 3 seconds to start/close massager

(2) Short press up/down key to change songs

(3) Long press up/down key to change volume

(4) Double click turn on/off key to open/close music

(5) Single click turn on/off key to change massager mode


(6) The host starts bluetooth automic open mobile search

corresponding type connection

Turnon/off key

Up Key

Down Key



4

PRODUCT USAGE


Method of use


1. Preparation before use

(1) take out the eye massager gently

(2) put the eye massage device on the eye position

2. Usage (as shown)





5


PRODUCT CLEANING, MAINTENANCE AND STORAGE


Product cleaning, storage and maintenance

1. Product deaning and maintenance

If the product is dirty, please wipe it with a dry soft cloth.

If it is particularly dirty, you can wipe it with medical alcohol.





6

2. Keep

Keep it dry, ventilated and clean

Please keep out of reach of children

Do not disassemble the main engine without permission

Bluetooth Connection

1.Open the Bluetooth on your mobile phone and search the Bluetooth Name: YBY. The device will auto connect.

2.This product could only connect one Bluetooth device. If want to change another device, please turn off the Bluetooth on the phone you are currently connected to.

3.Please connect the Bluetooth within 10s when turn on the device. If you want to connect during use, please restart the device.

4.After connecting, you could play music on your phone and enjoy that.

Tips: The bluetooth cannot connect successfully in low battery. First connection if fail to search Bluetooth "YBY", please long press the "Turn on/off" button to turn off the device then restart and also try to turn off the bluetooth on the phone and restart for several times.

Once connect successfully, the device will save the data and next time it will be connect at once and you don't need to repeat above process.

PRODUCT SPECIFICATION

The product name	Eye massager
Automatic timing	15 minutes
The rated voltage	DC 5V
Rated power	3.4W
Power adapter	AC 100-240V DC5V

7

COMMON FAILURES AND SOLUTIONS

Q: Is there a time setting function?

A: It has 15 minutes setting function

Q: What should I do if I want to turn off the sound after I turn it on?

A: The volume is adjustable. Long press 3s "Up/Down Key" turn up/down the volume. Double presses the "Turn on/off Key" turn off the sound.

Q: Why I can't charging?

A: Please use an new USB cable or AC adapter.

Q: Why the air press massage is not strong enough?

A: Low in power. Please charge.

Q: Why I can't connect the bluetooth?

A: Please restart the product and the bluetooth on the mobile phone to reconnect. This product could only connect one mobile phone. If want to change another device, please turn off the Bluetooth on the phone you are currently connected to.

Q: Is there a mode without background music?

A: YES. Double presses "Turn on/off Key" button turn off/ on the music. If during the condition of connecting the Bluetooth, please turn off the music on the phone.

After-Sales Service

If you meet any problems during use, please feel free to contact us at anytime.

Our Customer Support Team is ready to fix your issues.

8

MESSAGE BRINGS WISDOM TO LIFE

FCC Warning:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 0cm between the radiator and your body.