



DrKumo

Quick Start User Guide

Model: KumoPal™ XL100

Part Number: HUB-101A



Welcome! Let's Get Started

This Quick Start Guide introduces the main components of your KumoPal Hub, provides instructions to power it on, and includes important care and safety information. Follow the steps on the next pages to get your device ready for use.

What's in the Box

1. KumoPal Hub
2. Quick Start Guide
3. USB-C Charging Cable
4. USB-C Charging Adapter

Getting to Know Your Hub



1. Power Button
2. Volume + Button
3. Volume - Button
4. USB 2.0 Type-C Port
5. Speakers
6. Microphone
7. Headphone Jack
8. Front-facing Camera
9. Rear-facing Cameras
10. Flashlight

What each component does

1. Power Button

- Power On: Press and hold the Power button to turn on the Hub.
- Power-off: Press and hold the Power Button for a few seconds. Then, on your screen, tap Power off.
- Forced restart: Press and hold the Power button for about 8-9 seconds, or until your phone restarts.

2. Volume + Button

- Increase the volume of your Hub

3. Volume - Button

- Decrease the volume of your Hub

4. USB 2.0 Type-C Port

The USB-C charging port is where you plug in a charging cable connected to a power adapter to charge your device's battery or transfer data.

5. Speakers

Provides clear and immersive stereo sound directly from your tablet.

6. Microphone

Captures audio clearly for video

calls, voice commands, or audio recording.

7. Headphone Jack

Connect wired headphones or external audio devices for private listening or enhanced audio quality.

8. Front-facing Camera

Ideal for selfies, video conferencing, and facial recognition features.

9. Rear-facing Cameras

Captures high-quality photos and videos from the rear side of your tablet.

10. Flash Light

Provides illumination for capturing clearer photos/videos in low-light conditions or serves as a flashlight.

Starting your Hub

Step 1: Charge Your Hub

1. It is best to fully charge your Hub before you use it for the first time. This can take a few hours.
2. Connect the small end of the Charging Cable to the Charging Port on your Hub.
3. Connect the large end of the Charging Cable to the Power Adapter.
4. Plug the Power Adapter into a wall outlet.
5. A light or symbol on the screen will show you that it is charging.

Step 2: Turn On Your Hub

1. Once your Hub is charged, you are ready to turn it on.
2. Press and hold the Power Button for a few seconds. The screen will light up. It may take a minute to start.

Welcome to the DrKumo App

It will start and open the DrKumo program for you automatically.

You will not see other screens or apps like on other Hubs. This is normal.

Your care team has already set everything up for you.

Sleep and Wake Your Hub

- To save battery, you can put the screen to sleep.
- Quickly press and release the Power Button once.
- The screen will go dark.

Waking the Screen Up

- Quickly press the Power Button again.
- The screen will turn back on.

If You Need Help

What if my Hub won't turn on?

- Make sure it has been charging for at least 30 minutes. Use the cable and power adapter that came in the box.

What if the screen is dark?

- The screen may be asleep. Quickly press the Power Button once to wake it up.

What if I am still having trouble?

- Tap the  icon on your app for additional support.

Important Safety & Care

1. Please only use the charger that came in the box.
2. Keep your Hub away from water and liquids.
3. Do not try to open the Hub or replace the battery yourself.
4. If you're not using your Hub, lock the screen.
5. Lightly touch the screen. Do not press the screen hard.
6. Hearing Damage Risk: Do not listen at high volume levels for a long period.

7. Properly place / store it to avoid other objects falling on the Hub.
8. Do not place the Hub directly on the bed before going to bed because heavy pressure is likely to damage the screen during deep sleep

Certification Statements

FCC-USA

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions.

This device may not cause harmful interference.

This device must accept any interference received, including interference that may cause undesired operation.

FCC ID: 2BP5Q-XL100.

CAUTION: Modification or dismantling of this equipment not expressly approved by Manufacturer may void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with

the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

SAR Information Statement

Your wireless phone is a radio transmitter and receiver. It is designed and manufactured not to exceed the emission limits for exposure to radiofrequency (RF) energy set by the Federal

Communications Commission of the U.S. Government. These limits are part of comprehensive guidelines and establish permitted levels of RF energy for the general population. The guidelines are based on standards that were developed by independent scientific organizations through periodic and thorough evaluation of scientific studies. The standards include a substantial safety margin designed to assure the safety of all persons, regardless of age and health. The exposure standard for wireless mobile phones employs a unit of measurement known as the Specific Absorption Rate, or SAR. The SAR limit set by the FCC is 1.6 W/kg. * Tests for SAR are conducted with the phone transmitting at its highest certified power level in all tested frequency bands. Although the SAR is determined at the highest certified power level, the actual SAR level of the phone while operating can be well below the maximum value. This is because the phone is designed to operate at multiple power levels so as to use only the power required to reach the network. In general, the closer you are to a wireless base station antenna, the

lower the power output. Before a phone model is available for sale to the public, it must be tested and certified to the FCC that it does not exceed the limit established by the government adopted requirement for safe exposure. The tests are performed in positions and locations (e.g., at the ear and worn on the body) as required by the FCC for each model. The highest SAR value for this model phone when tested for use on the body, as described in this user guide, is 1.160W/Kg(Body-worn measurements differ among phone models, depending upon available accessories and FCC requirements). While there may be differences between the SAR levels of various phones and at various positions, they all meet the government requirement for safe exposure. The FCC has granted an Equipment Authorization for this model phone with all reported SAR levels evaluated as in compliance with the FCC RFexposure guidelines. SAR information on this model phone is on file with the FCC and can be found under the Display Grant section of <http://www.fcc.gov/oet/fccid> after

searching on

FCC ID: 2BP5Q-XL100 Additional information on Specific Absorption Rates (SAR) can be found on the Cellular Telecommunications Industry Association (CTIA) web-site at <http://www.wow-com.com>. * In the United States and Canada, the SAR limit for mobile phones used by the public is 1.6 watts/kg (W/kg) averaged over one gram of tissue. The standard incorporates a substantial margin of safety to give additional protection for the public and to account for any variations in measurements.

Body-worn Operation

This device was tested for typical body-worn operations. To comply with RF exposure requirements, a minimum separation distance

of 10mm must be maintained between the user's body and the handset, including the antenna.

Third-party belt-clips, holsters, and similar accessories used by this device should not contain any metallic components. Body-worn accessories that do not meet these

requirements may not comply with RF exposure requirements and should be avoided. Use only the supplied or an approved antenna.

FCC Radiation Exposure Statement:

The device has been tested and comply with FCC SAR limits

This device has been tested and comply with IC SAR limits.

Operation of this device in the band 5150-5250 MHz is restricted to indoor use only. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.



ENERGY STAR® Compliance



The DrKumo Hub meets the ENERGY STAR guidelines for energy efficiency. Reducing energy consumption saves money and helps conserve valuable resources. For more information about ENERGY STAR, go to energystar.gov.

For more detail, visit:

<https://hub.drkumo.com>

