

Oral Device Instruction for Use

The Oral Device is intended for use by the orthodontic patient during treatment with aligners to facilitate minor anterior tooth movement. The vibrations generated by the device gently remodel the alveolar bone to accelerate up anterior tooth movement when the teeth are under force and to secure the teeth in retention.

Before first-time use: DEVICE MUST BE FULLY CHARGED. IT WILL NOT ACTIVATE UNTIL IT IS CHARGED.

1. Oral Device is shipped in “shipping mode” to save battery life during transport.
2. To charge the device, plug the USB-C cable into the back of the device.
3. When fully charged, unplug the Oral Device device and follow the General Instructions.
4. Note: The recommended charging time is once per week.

How to use Oral Device?

Remove the Oral Device device from the case.

1. Rinse the mouthpiece.
2. Gently bite down on the mouthpiece.
3. Press the “On” button to begin the session. (There is an automatic shut-off feature after 5 minutes of use.)
4. Rinse and dry the mouthpiece after use.

How often to use Oral Device?

Use the Oral Device device as indicated by your physician.

Charging and Cleaning

Charging

1. Charge by connecting the charging cable directly to the back of the device and to a USB power source.
2. Oral Device's LED lights will flash once as charging begins. The number of LEDs illuminated will correspond with the percentage the device is charged. When charging is complete, the LEDs will automatically turn off. The recommended charging time is once per week.
3. A rechargeable lithium-ion polymer battery pack internally powers the device. The device is inoperable while charging is occurring.

Storage

1. Proper care should be taken to ensure that the Oral Device device is not damaged. The device should be stored at room temperature.

2. Storage/Operating Temperature: -20°C to 45°C Storage/Operating Humidity: 15 to 90% Storage/Operating Atmospheric Conditions: 700 kPa to 1060 kPa

3. The device is water-resistant, NOT waterproof:

- DO NOT submerge in water
- DO NOT use in the shower
- Keep USB port clean and dry
- Wipe dry before charging after cleaning


Mouthpiece Care

Cleaning

The mouthpiece should be rinsed with water after each use. For a deeper cleaning, brush the mouthpiece with a toothbrush during rinsing.

Mouthpiece Replacement

1. The mouthpiece can be detached from the actuator for cleaning and replacement.
 - (A) Turn the mouthpiece one-quarter counterclockwise and gently pull to detach.
 - (B) To re-attach, place the mouthpiece on the actuator and turn one-quarter clockwise.



Contraindications

- Do not use Oral Device if you have conditions such as vertigo to avoid symptoms such as headache, dizziness, and aching facial pain.
- Do not use Oral Device if you are pregnant.

Factory Reset

In the event your Oral Device device becomes non-responsive, you may perform a factory reset by pressing and holding the “On” button for 10 seconds. The device will automatically reset.

Device Notifications

LEDs slow breathing	Device is Charging
Four solid LEDs	Battery charged
Three solid LEDs	Battery 60%-80% charged
Two solid LEDs	Battery 40%-60% charged
One solid LEDs	Battery 20%-40% charged
One blinking LED	Low Battery



Precautions

- It is important to read the Oral Device “Instructions for Use” before operating the device to ensure proper use, maintenance, and care.
- Oral Device is intended to be re-used.
- Before use, inspect the Oral Device device and packaging to verify that no damage has occurred due to shipping and handling. Do not use damaged products.
- Always recharge Oral Device according to the “Instructions for Use” before use. Refer to the CHARGING section for more information.
- Never use any other alternative power sources or instruments in conjunction with Oral Device.
- Use Oral Device only as instructed.
- Oral Device is water resistant, not waterproof.
- DO NOT submerge in water.
- DO NOT use in the shower.

- Keep USB port clean and dry.
- Wipe dry before charging & after cleaning.



Warnings

- Discontinue usage if you experience headache, dizziness, discomfort, or pain.
- Do not make any modifications to the upLift 5 device. Any modification may impair the functionality of the device and void the warranty.



Disposal

Used electrical and electronic equipment (WEEE) should not be mixed with general household waste. For proper treatment, recovery, and recycling, please take this product(s) to designated collection points where it will be accepted free of charge. Disposing this product correctly will help save valuable resources and prevent any potential negative effects on human

health and the environment, which could otherwise arise from inappropriate waste handling. Please contact your local authority for further details of your nearest designated collection point/household hazardous waste collection point. Penalties may be applicable for incorrect disposal of this waste, in accordance with your national legislation.



Exposure & Interference

- This equipment has been tested and found to comply with the limits for a Class B digital device pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in an installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures: Reorient or relocate

the receiving antenna. Increase the separation between the equipment and the receiver. Connect the equipment to an outlet on a circuit different from that to which the receiver is connected. Consult the dealer or an experienced radio/TV technician for help.

- Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
- This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) the device must accept any interference received, including interference that may cause undesired operation.



Replacement

If your Oral Device device and/or any of its components no longer work properly or require repair during the warranty period, contact your doctor. Do not attempt to service or perform maintenance on the device. No modification of the device is allowed.



Service Life

The recommended service life of Actuator is two years from the date of manufacture.



Limited Warranty

This policy is applicable to domestic US purchases only: Oral Device will replace authentic devices and attachments eligible for warranty service with merchandise of equal or better quality. Every new Oral Device device sold includes a two (2) year warranty from date of purchase covering any manufacturer-

related defects resulting from everyday use of the product. Warranties are only valid for items purchased through authorized sellers, returned within the warranty dates, and accompanied by a physical product purchase receipt from the location purchased. This warranty is void if the product has been subject to an accident, misuse, abuse, improper maintenance or repair, unauthorized modification, or any other use not found in the printed instructions. For devices purchased outside of the US, International Returns and International Warranties: All product returns outside the United States are per your place of purchase: spa, retailer, or distributor. Please consult your receipt or place of purchase for their return and warranty policy.

Technical Specification & Symbols

Warning: No modification of this equipment is allowed.

Power supply	3.7 Volt Battery, Internally Powered
Equipment Type	 Bluetooth Pairing Active
Ingress Protection Rating: IPx 3	Water falling as a spray at any angle up to 60° from the vertical shall have no harmful effect.
	Bluetooth Enabled
	FCC Declaration of Conformity
	Waste Electrical and Electronic Equipment Directive (WEEE)
	Restriction of Hazardous Substances Directive (RoHS)
	Double Insulated Electrical Appliance
	Caution
	The official mark of Europe Certificate
	Medical Device
	Single patient multiple use
	Use by prescription
	Manufactured Date
	Batch code
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation
	Type BF applied part
	Non-sterile



For: Medina Management, LLC
18975 Collins Ave
Unit 3905
Sunny Isles Beach, FL, 33160
United States
DUNS: 144958756