

Air Compression Recovery System



Explore
More with Us



Scan for
Customer Support

Stay Connected and Supported.

S9083

Thank you for choosing the Air Compression Recovery System.

Established in 2001, CINCOM is a leading global innovator in healthcare and wellness solutions, specializing in advanced dynamic air compression massage devices. Our touchscreen leg recovery system features intelligent air compression, single leg massage, and precision independent zone massage to provide a relaxing spa experience. The technology is designed to simulate reflexology, providing targeted recovery for your feet, calves, thighs, and entire legs.

With a proven track record of enhancing wellness for millions of customers worldwide, CINCOM remains committed to delivering premium products and exceptional service.

Air Compression Recovery System



Explore
More with Us



Scan for
Customer Support

Table of Contents

What's in The Box?	1~1
Terms of Use	2~2
Precautions	3~4
Warning	4~4
Charging Cautions	5~5
Product Description	6~7
Operations	8~12
Bluetooth Pairing & Disconnection	13~13
Specification	14~14
Troubleshooting	15~16
Maintenance & Cleaning Instructions	16~16
FCC Approved Equipment Authorization	17~17
Warranty	18~19
Labels	20~20



There is a single instruction manual
for this product in all sizes.

What's in The Box?

2 x Massage Sleeves(with Controllers) →



1 x Power Adapter / DC 5V3A →



1 x Type-C Charging Cable →



1 x Storage Bag →



1 x User Manual →



1 x Quick Start Guide →



1 x Shoulder Strap →



Terms of Use

This product is provided by CINCOM to the purchaser on the condition that the purchaser and user acknowledge and accept full responsibility for its proper and safe use. Improper, careless, or unsafe use of the CINCOM Air Compression Recovery System or any of its accessories may result in serious injury or death.

By using this product, the user assumes all risks and liabilities associated with its use, including but not limited to misuse, improper handling, failure to follow instructions, or use in violation of applicable laws or regulations. CINCOM shall not be held liable for any personal injury, property damage, or loss of life arising from the use or misuse of this product under any circumstances, whether such use is intentional, reckless, negligent, or accidental.

To the fullest extent permitted by law, CINCOM disclaims any liability for special, indirect, incidental, or consequential damages arising out of the assembly, installation, operation, or use of this product, including any accessories or materials required for its use. By operating or using this product, the user agrees to release and hold harmless CINCOM from any and all claims, liabilities, or damages related to its use.

If you do not agree to accept full responsibility for the use of this product, please do not operate it. Do not open or use any of the enclosed materials. Return the product to CINCOM in new and unused condition. Please note that products that have been used or are no longer in "as-new" condition cannot be returned or exchanged.

All information contained in this User Manual is subject to change without prior notice.

CINCOM reserves the right to make changes or improvements to its products at any time without obligation to update or retrofit previously sold units. If you have any questions regarding your product or its operation, please scan the QR code provided to contact CINCOM Customer Support.

Precautions

The CINCOM Air Compression Recovery System is designed as a personal care product. It is not intended for diagnosis/treatment of any health condition. Operate the CINCOM Air Compression Recovery System sensibly and with care, and it will be exciting, safe, and relaxing for you. Failure to operate your model in a safe and responsible manner may result in property damage and serious injury. We strongly recommend that people who have very sensitive legs use the lowest intensity. Please contact your physician before using the CINCOM Air Compression Recovery System if you:

1. Never charge the device that exceeds AC100-240V as it could cause fire.
2. Make sure the power adapter fits entirely into the port.
3. Avoid placing any full glasses or liquid bottles next to the device as they could spill and damage the device.
4. If the adapter cable is damaged, do not charge the unit.
5. Do not leave the product unattended when plugged in. Unplug from outlet when not in use.
6. Do not operate after taking pain killers, sedatives or alcohol.
7. Individuals with pacemakers should not use this product. Always consult a doctor before using it.
8. Keep cord away from heated surfaces.
9. Do not use with wet hands or in humid places. Wetness and humidity can cause electric shock or mechanical malfunctioning.
10. Do not attempt to disassemble or repair the unit, charging cord or screen. (doing so voids the warranty).
11. If you have a medical condition like diabetes, inflammation, suffer from blood clots, recently had surgery and/or are pregnant, suffering from heart disease or osteoporosis, we recommend you consult your doctor before using the massager.
12. Never use the device more than 1 hour in a row.
13. Never cover the massager while using it, this may result in overheating of the unit.
14. Never use the massager in a wet/humid room, in high temperature or next to corrosive gas.
15. Don't forget to turn off the device once you finish the massage, in order to protect the device.

16. Use this appliance only for its intended use.

17. Most importantly, use good common sense at all times.

©2026 CINCOM. All Rights Reserved. (Manufactured by Shenzhen Dongjilian Medical Tech Co., Ltd.)

WARNING

Do not use the device to supplement or replace any medical treatment that you are currently receiving. Please consult your physician before use if you are experiencing any unexplained pain, already receiving medical treatment for circulation or vascular issues, or using at home medical equipment.

Please contact your physician before using the CINCOM Air Compression Recovery System if you have physical health concerns or have any of the following conditions:

- You are pregnant.
- You have a pacemaker or other medical implant.
- You are suffering from heart or vascular disease.
- You have irregular blood pressure or are receiving treatment to normalize your blood pressure.
- You are suffering from osteoporosis.
- You are suffering from cancer.
- You have recently undergone surgery.
- You have concerns regarding areas of your health.

NEVER leave the appliance unattended, especially if children are present. We recommend you limit the use of the massager to 30-minute sessions and for no more than an hour (60 minutes) consecutively.

NEVER use this product directly on swollen or inflamed areas of skin or open wounds.

NEVER be used by an individual suffering from any physical ailment that limits the user's mobility.

NEVER be used on persons who have sensory deficiencies affecting their hands and appendages.

Charging Cautions

This massager is capable of extreme performance. Responsibility, maturity, and common sense are required. This model must be treated with respect and caution to ensure your safety.

Operating the massager in a careless, unsafe manner, without adequate care and preparation can result in serious injury. The massager includes a lithium battery. Batteries have a severe risk of fire and injury if they are improperly handled, abused, or misused. The user must read and understand all included instructions, warnings, and precautions regarding the handling, care, and use of batteries.

- BEFORE you charge, ALWAYS confirm that the charger/plug settings exactly match the type (chemistry), specification, and configuration of the battery to be charged. DO NOT exceed the maximum manufacturer recommended charge rate which is 5V/3A.
- ALWAYS charge batteries in a well-ventilated area.
- Charge the battery in a safe area away from flammable materials. Monitor the charging process and never leave batteries unattended while charging.
- DO NOT allow children to charge the massager.
- DO NOT disassemble, crush, short circuit, or expose the battery to flame or other source of ignition.
- DO NOT leave the massager unattended while charging, discharging, or anytime the charger is ON with a battery connected. If there are any signs of a malfunction unplug the power source and/or stop the charging process immediately.
- This device cannot be used while charging. Please unplug the charger once the battery is fully charged.
- If the massager gets hot to the touch during the charging process (temperature greater than 110°F / 43°C), immediately disconnect the battery from the charger and discontinue charging.
- Keep a Class D fire extinguisher nearby in case of fire.
- DO NOT expose the unit to water or moisture.



WARNING: RISK OF INJURY!

We strongly recommend that people who have very sensitive legs use the lowest intensity.



CHARGING OR DISCHARGING BATTERIES POSES A NUMBER OF SIGNIFICANT HAZARDS WHICH CAN INCREASE THE RISK OF FIRE, EXPLOSION, SERIOUS INJURY, AND PROPERTY DAMAGE IF NOT PERFORMED AS PER THE INSTRUCTIONS. BEFORE USE, READ AND FOLLOW ALL MANUFACTURER'S INSTRUCTIONS, WARNINGS AND PRECAUTIONS.

NEVER ALLOW CHILDREN UNDER THE AGE OF 14 TO USE THE MASSAGER WITHOUT THE SUPERVISION OF A RESPONSIBLE, KNOWLEDGEABLE ADULT.

THIS IS NOT A TOY.


Product Description

CONTROL UNIT




Type-C Charging Port

Operations



LEG SLEEVES



INDICATIONS FOR USE

The Air Compression Recovery System is indicated for the temporary relief of minor muscle aches and pains and for temporary increase in circulation to the treated areas in people who are in good health. Air Compression Recovery System kneading and stroking of tissues by using an inflatable garment.

- After unpacking the massager, please check the integrity of the device accessories and look for any damage.
- Before first use, charge the massager for 3-4 hours.
- The CINCOM Air Compression Recovery System is an air pressure massager intended to temporarily relieve minor muscle aches and/or pains, and to temporarily increase circulation to the treated area. Based on the principles of human engineering, a unique design and the distribution of acupuncture points, the CINCOM Air Compression Recovery System meets the demands of your leg according to its structure. It will be like having a professional physiotherapist at home.
- Please be sure to remove all sharp jewelry. Do not cross your legs.
- If the CINCOM massage feels uncomfortable, you can reduce the intensity or stop the session.

BEFORE USE: FULLY CHARGED

Use the included 5V/3A power adapter to charge the device before first use. White light blinks - Charging in progress; White light illuminates steadily - Fully charged.

START BUTTON

***After the device is powered on, it will enter standby mode. Press the Start Button  to begin operation.**

START:

Press the Start Button . You will hear a single "beep", and the massager will start working.

PAUSE/POWER OFF

Press the Start Button  again to pause the massage. You will hear beeping sounds, and the air chambers will gradually deflate. If there is no operation within two minutes, the device will automatically shut off.



TIME BUTTON



The start on default time is 30 minutes. The time setting range is 15-60 minutes. Touch the " < " or " > " button to set the time as needed, adjusting ± 5 minutes at a time. After the time has expired, the device will shut off.

STRENGTH BUTTON



The default starting pressure is 50 mmHg. The pressure can be adjusted within a range of 50-150 mmHg. Press the " < " or " > " buttons to increase or decrease the pressure in 10 mmHg increments, allowing you to customize the intensity to your comfort level. For reference, the pressure levels are generally categorized as follows:

- **Low: 50–80 mmHg**
- **Medium: 90–120 mmHg**
- **High: 130–150 mmHg**

If you have sensitive legs or are using the device for the first time, we strongly recommend starting at the lowest pressure setting (50 mmHg) and gradually increasing the level as your body adapts.


MODE BUTTON




The default program mode is Combination, in which all 5 air chambers (zones) operate simultaneously. Using built-in pressure sensing technology, the system applies and releases pressure across different areas of the leg based on its natural shape and contours. As a result, each user may experience a slightly different, personalized massage. Press the Mode Button (MODE) to cycle through and select your preferred program mode.

LEG MASSAGER MODES


Mode 1: Combination Mode




Mode 2: Circulation Mode



Mode3: Sequence Mode



CHARGING & POWER MANAGEMENT



Battery Performance: Full charge provides 1.5-2.5 hours of typical usage.
For best results: Complete initial charge before first use.

⚠ CAUTION

- Make sure both users complete the unpairing steps before starting a session.
- If Bluetooth remains connected, settings adjusted on one unit will continue to affect the other unit.
- To restore synchronized dual-leg operation, power on both units and follow the Bluetooth pairing instructions.

⚠ SAFETY WARNING - Emergency Stop Switch

Before use, make sure the switch is in the locked position .
To activate the emergency stop, slide the switch in the direction of the arrow.



Emergency Stop Switch




This function is intended for emergency use only. During normal operation, the switch should remain in the locked position. If the system malfunctions or the air pressure becomes excessively high, move the switch to the arrow icon position to immediately shut down the device.

When the emergency stop is activated, the device will emit 6 beeps, the display will show "— —", and the system will automatically power off.


LED POWER INDICATORS

STATUS	LIGHT INDICATOR
Charging in progress	White light blinking
Fully charged	White light steady
Battery >50%	White light on
Battery 20%-50%	Orange light on
Battery 10%-20%	Red light on
Battery <10% (Critical)	Red light blinking
Low Battery Alert	When red light blinks, recharge soon to maintain optimal performance.


HOW TO FOLD THE SLEEVES?




① Fold the sleeve in half



② Fold the sleeve upwards from the heel area



③ Fold up 2-3 times




④ Place the leg massager in the storage bag

Bluetooth Pairing & Disconnection


PAIR DEVICES

Touch and hold the Mode Button (MODE) on both controllers for 5 seconds until you hear a beep. The Bluetooth indicator lights will flash for approximately 30 seconds. Once the lights stop flashing, the devices are successfully paired.



UNPAIR DEVICES

Touch and hold the Mode Button (MODE) on one of the controllers for 5 seconds until you hear a beep. The Bluetooth indicator light will flash for a while before turning off, indicating that the devices have been successfully unpaired.



*If two users wish to operate the leg massagers independently, the Bluetooth connection must be disconnected first.

Specification

Product Name	Air Compression Recovery System
Model Number	S9083
Dimensions	418×195×235mm /16.5×7.7×9.3 in
Battery Capacity	7.4V 2600mAh
Charging Time	3-4 Hours
Adapter	AC Input: 100~240 Volts, 50/60HZ, DC 5V/3A
Timer	15Mins-60Mins, Default Value is 30Mins
Strength	50-150mmHg, Default Value is 50mmHg
Compression Mode	3 Modes: Combination, Circulation, Sequence
Sleeve Type	5 Chambers Sleeve
Material	• Sleeves: Durable Nylon • Unit Control: Fireproof ABS
Operating Conditions	• Temperature: +5°C to 35°C (41°C to 95°C) • Humidity : 5% to 80% non-condensing • Atmospheric Pressure : 70 kPa to 106 kPa
Storage Conditions	• Temperature: -20°C to 40°C • Humidity: 30% to 85% non-condensing • Please keep dry and avoid direct sunlight exposure

Troubleshooting

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
Noise during operation	1. "Clicking" sound from airbag inflation. 2. "Hum" sound from the air pump.	These are normal operational sounds produced by the airbag and pump during inflation and deflation.
Cannot charge or indicator light not on	1. The Type-C connector is not properly connected to the control unit. 2. The adapter plug is not securely inserted.	Ensure the charging cable and adapter plug are firmly connected.
Device stops working during operation	Timer has reached the preset duration.	To continue using, restart the device.
Charging indicator flashes red, beep sounds ten times, and device shuts down	Low battery – please recharge the control unit before use.	Please fully charge the control unit before using the device. (Note: The device cannot operate while charging.)
Beeps six times and shuts down/continuous beeping	The Emergency Stop Switch is in the "→" position.	Set the Emergency Stop Switch to the "🔒" position.
Device does not operate after powering on and shuts down after 2 minutes.	The Start/Pause button (▶) was not pressed after powering on.	After powering on, tap the Start/Pause button (▶) to begin operation.
Device shuts down automatically	The device was charging while turned on.	Ensure the device is turned off during charging.

Bluetooth indicator not on / left and right leg sleeves do not operate synchronously.

The units are not successfully paired.

Re-pair the units following the instructions provided on P13 of this manual.

If you have any questions, concerns or if you would like to order any spareparts or replacements. Please scan the QR code to contact us and we will get back to you as soon as possible.

Maintenance & Cleaning Instructions

MAINTAINING THE SYSTEM

The CINCOM Air Compression Recovery System requires no routine maintenance or service except for the care in this section.

STORING THE SYSTEM

Store your Air Compression Recovery System in a clean, dry location.

CLEANING THE DEVICE

To clean the Control Unit:

- Wipe down the device with a damp, clean cloth.
- Dry thoroughly with a clean cloth.

CLEANING THE LEG ATTACHMENTS:

- Wipe down the leg attachments inside and out with a damp, clean cloth.
- Dry thoroughly with a clean cloth.
- Do not machine wash or dry.
- Do not dry clean.

FCC Approved Equipment Authorization

FCC ID: 2BQV7-S9083

FCC Warning

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

NOTE 2: Any changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Warranty

CINCOM has made every effort to ensure your product is as strong and durable as possible. However, it is still possible that some mishaps could happen and we take full responsibility.

To guarantee you have the best experience we offer a 1-year warranty covering defects in materials, workmanship and assembly on all our products. To ensure that you don't invalidate the terms of the warranty please read and understand this booklet.

COVERED BY WARRANTY

- If a defect in materials, workmanship or assembly is found BEFORE the product has been used, we will exchange the entire product, if the product is returned complete and in perfect condition in its original packaging.
- 1-year warranty for the Air Compression (excludes damages caused by the user. Including but not limited to incorrect maintenance.)
- 6-month warranty for the Battery (excludes damages caused by the user, including but not limited to incorrect maintenance.)

**When the air compression, the battery (or the programs) start to malfunction, please contact us at service@cincomhealth.com and we will provide a replacement.

THIS IS A COMPLEX HIGH PERFORMANCE AIR COMPRESSION RECOVERY SYSTEM AND MUST BE TREATED WITH CARE AND RESPECT.

NOT COVERED BY WARRANTY

- If this warranty certificate and purchase records cannot be presented.
- Improper use caused by artificial damage, such as access to improper power supply, usage of unsuitable accessories, improper installation, and not following the instructions.
- Damage caused by transportation and other accidents. Repair and modification not approved by the Company. Further damage caused by force majeure (such as natural disasters, voltage abnormalities, etc.)
- Regular use of the product is caused by aging but does not affect the normal use of the product.
- This free service does not include the needs to replace the consumable parts (batteries, chargers, power cord, massage head etc.) and accessories (cloth bags, storage bags, etc.), auxiliary devices.
- CINCOM DOES NOT operate a new for old warranty once the product has been used.

**Warranty is only valid if the items were purchased directly from CINCOM (purchases from third party sellers will not be covered under the warranty)

**IT IS THE PARENTS' OR GUARDIANS' RESPONSIBILITY TO ENSURE MINORS RECEIVE APPROPRIATE GUIDANCE AND/OR SUPERVISION WHEN USING THIS PRODUCT.

If you have any questions, concerns or if you would like to order any spareparts or replacements. Please feel free to contact us at service@cincomhealth.com and we will get back as soon as possible.

Labels

The following labels and symbols appear on the product and/or packaging.

LABEL	DESCRIPTION	LOCATION
	Read instructions before use	On packaging label
	Warnings, read instructions before use	On packaging label
	Double Insulation. Equipment consists an extra layer of supplementary insulation over the basic insulation, or a single layer of special reinforced insulation directly over the live parts.	On packaging label
	Shenzhen Dongjilian Medical Tech Co., Ltd. Building B1, No. 70 Tongxin Road. Baolong Street, Longgang District, Shenzhen, China	On manual label
	Separate collection for waste electrical and electronic equipment.	On product label
	In accordance with Directive 2014/35/EU electrical equipment designed for use within certain voltage limits, and Directive 2014/30/EU electromagnetic compatibility	On product label
	The electromagnetic radiation from the device is below the limits specified by the Federal Communications Commission.	On product label
RoHS	A directive in the European Union that regulates the use of certain hazardous substances in electrical and electronic equipment (EEE).	On product label