



LEAD A HEALTHY NEW LIFE



+ FOLLOW US



✉ service@quinear-ec.com

🌐 www.quinear-ec.com



Air Compression
Recovery System

S9086

Contents



Safety Precautions	1
FAQS	2
Battery Charging Instructions	4
Specifications	5
Troubleshooting	5
Product Description	6
Operating Instructions	8
How to use the Bluetooth	10
Storage & Maintenance	11
Cleaning	11
Disposal	11
Package includes	11
Contact Us	12
FCC Approved Equipment Authorization	12
Labels	13



Please take a moment to thoroughly read the user manual before enjoying your QUINEAR

Safety Precautions



Please read the following safety precautions carefully to ensure safe and correct use of the QUINEAR Air Compression Leg Recovery System.

■ Warnings

No modification of this equipment is allowed.

1. If you experience severe pain, any unusual symptoms, or want to remove the sleeves in an emergency during use .

- Press the power button to stop the device or turn on the emergency brake switch.
- Remove the sleeves and consult your doctor.

2. Consult your doctor before using this device if you have any of the following conditions.

- Experiencing inflammation, an infection, pain of unknown origin, at or near the site of application.
- Have bone fractures or dislocations at or near the site of application.
- Receiving medical treatment, or have a contra indication requiring the use of any medical control unit.
- Using pacemaker or other medical devices which are susceptible to electrical interference.
- Suffering from malignant tumors.
- Suffering from cardiac diseases.
- Suffering from osteoporosis.
- Having serious peripheral neuropathy dysfunction or sensory disturbance caused by diabetes.
- Being unsuitable to do the massage because of traumas on the body.

3. Keep it out of the reach of infants, children and people without the ability to use it independently.

4. It is recommended to use the original power adapter.

5. Don't scratch, damage, process, excessively bend, pull or twist the power cord of power adapter.

Otherwise, it may cause fire or electric shock.

6. Don't plug or unplug the power adapter with wet hands.

7. Don't put the air pump controller in the quilt or use the machine in high temperature condition.

8. It is forbidden to remodel, disassemble or repair the machine without permission.

■ Cautions

1. Only use the power supply provided with the system. Using a different power supply may cause the system to operate incorrectly.
2. Don't use it in the bathroom or other humid places.
3. Keep it out of the reach of children.
4. Don't walk around when you are using this device or wearing the sleeves.
5. Do not leave the system, power supply, or any accessories where they could be damaged by children, pets, pests, or liquids.
6. The attachments are designed to be used by only one person at a time.

FAQS


Q1: How dose this product work?

A1: This is an air compression device, 5 air chambers(foot, lower calf, upper calf, lower thigh, upper thigh) inflated and deflated with a sophisticated sequence to stimulate hand kneading and stroking, maximizes blood flow throughout the body to help with reduced muscle soreness and/or pains to the treated areas, and improve circulation.


Q2: How many massage modes, and what's the difference?

A2: There are 3 massage modes, please check the difference as below:


Mode 1: Combination Mode



Mode2: Circulation Mode



Mode3: Sequence Mode



Q3: What should I do if I feel the massage strength is too weak or too strong?

A3: There are 11 levels of massage intensity available, 50-150 mmhg Adjustable. The default setting is 50mmhg at startup, and 10mmhg at each adjustment, so you can flexibly choose the strength you want. If you still feel it's too strong at the lowest setting, please stop using the device, we suggest to consult your doctor before next use.

Q4: How long should I use it?

A4: It is recommended to use it once or twice a day, no more than 30 minutes at a time. But you can use it longer if you feel too tired and want it to massage longer to work on the soreness.

Q5: Can I only apply on one leg?

A5: Turn on one leg and do not turn on the other leg. The Bluetooth indicator light will not light up. (Please check P9-P10 to see more).

Q6: Why does the air pump controller make sound?

A6: The sound comes from the working air pump in the control unit, providing air continuously to the airbags in the sleeves, it's a low sound and that's a normal phenomenon.

Q7: When should I use this product?

A7: This product pushes the blood to flow, effectively remove muscle soreness and fatigue, you can use it anytime you want. You can use it before workout to warm up your muscles, use it after exercise for better and faster recovery, or use it when you back home after a long day work.

Q8: Can I synchronize the massage settings for both leg covers

A8: Press the MODE button on the host for 5 seconds at the same time. After hearing a beep, release it at the same time. At this time, both Bluetooth indicator lights will flash. After about 30 seconds, both Bluetooth indicator lights will light up at the same time, indicating a successful connection.(Please check P10 to see more).

Battery Charging Instructions

Battery Capacity	2600mA/h
Charging Time	2H
Using Time with One Charging	2-4H (Depending on intensity and mode)
When First Charging	Please make sure the machine is fully recharged for first charging

Status	Battery Indicator Light	Description
Charging Reminder	1. Red flashing; 2. Battery indicator light flashes red, emits a beep sound, and automatically shuts off.	Low battery (below 10%) or critically low: Unable to power on. Please charge immediately.
Low Battery Reminder	Red	Battery is below 20%.
Medium Battery Reminder	Yellow	Battery is below 50%.
High Battery Reminder	White	Battery is above 50%.
Charging	White flashing	Charging in progress.
Fully Charged	White solid	Charging complete.

Specifications

Model	S9086
Name	Air Compression Recovery System
Adapter	AC Input: 100-240 Volts, AC 50/60Hz DC Output: 5V 3A
Rated Power	15W
Compression Mode	3 Modes: Combination, Circulation, Sequence
Intensity Levels	11 Levels: 50-150 mmhg Adjustable
Operating Time	15Mins--60Mins, Default value is 30Mins
Product Weight	4.25kg / 9.37 lbs
Sleeve Type	5-Chambers Sleeve - Foot / Low Calf / Upper Calf / Low Thigh / Upper Thigh
Leg Boots Size Small	Inseam: 28.4 in; Thigh Circumference: 28.2 in; Foot Length: 10.6 in
Leg Boots Size Medium	Inseam: 29.5 in; Max Thigh Circumference: 32.5 in; Max Foot Length: 12.8 in
Leg Boots Size Large	Inseam: 32.3 in; Thigh Circumference: 34.1 in; Foot Length: 13.4 in


Troubleshooting

Problems	Possible Causes	Solutions
Doesn't work and screen is off	1. Power is not turned on. 2. The controller's out of power.	1. Press the power button to turn the controller on. 2. Charge the controller fully before the first time using.
After booting up, the host beeps 6 times and shuts down	The emergency brake switch is in the "arrow" position	Push emergency brake switch to the "locked" position
After turning on the product, it does not work. After 2 minutes, it will beep once and turn off	Not pressing the start button after booting up	After turning on, you need to touch the start button to work
Unsynchronized leg massage	The Bluetooth is not connected	Re pair devices (Please check P10 to see more)
The Bluetooth lights on both legs are not on	The Bluetooth is not connected	Re pair devices (Please check P10 to see more)
The host cannot charge or the power light is not on	1. The charger plug is not securely plugged into the power supply 2. The power cord is not securely plugged into the socket	Please connect the charger and power cord correctly
Stop working during product operation	Time is up	Please restart the device to continue using
The power light flashes red and shuts down after 10 beeps	Low battery, remind to charge	Please charge the product


* You're welcome to reach us at service@quinear-ec.com if further assistance is needed.

Product Description

● Air Pump Controller



● Massage Sleeve



Operating Instructions

Please read safety precautions carefully before using the product. In order to get better massage, please use it according to the instruction.

Tips

* Please make sure the zipper is zip up completely before wearing. Don't try to use the system with the zipper partially or totally unzipped.


* Please note that the product cannot be used while charging. It can only be charged when it is turned off, and cannot be charged when it is turned on.

* Please ensure that the emergency brake switch of the device is locked before use.

1. Put the sleeves on your leg properly.




2. Confirm that the emergency brake switch is in the locked state



Emergency Brake Switch

3. Press 2s on the "⏻" Power button to turn the power on,




4. Touch the "▶" to begin massage.




5. Find a comfortable position lying down, you can begin the massage. It will begin with default setting of: TIME-30 Mins, INTENSITY-50 mmHg, MODE-COMBINATION.

6. Touch the "M" mode button to switch the massage mode. There are 3 modes that can be switched at will: Combination mode, Circulation mode, Sequence mode.



Pattern selection


7. Touch the "<" or ">" intensity button to increase or reduce the compression intensity, 11 Levels: 50-150 mmHg adjustable, each touch +/- 10 mmHg.



Intensity reduction

Intensity increase

8. Touch the "<" or ">" to increase or reduce the massage time. 15Mins-60 Mins adjustable, default value is 30 Mins, it will adjust 5 Mins every time you adjust. The time can be adjusted anytime while the session is running.



Time reduction

Time increases

The minimum limit is 15 minutes, the maximum limit is 60 minutes, and the time range for each adjustment is 5 minutes

(Minimum limit 15 minutes)
(Maximum limit 60 minutes)

9. The device will automatically shut off after the treatment time is up, you can press the power button and touch the start button if you wanna enjoy the massage again at any time.



Start/Stop

Single-Leg Massage Operation Introduction

Just Turn on one leg sleeve, the bluetooth will not connect.



No Light


How to use the Bluetooth

Pair devices

1. Simultaneously touch the MODE button for 5 seconds, beep once, and then release it at the same time. At this time, the link indicator light will flash.




2. The first connection takes about 30 seconds, and the Bluetooth indicator light will light up, indicating a successful connection.



At this point, only one host needs to be operated, and the two hosts will synchronize their changes.

Unpair devices

1. Touch any of the MODE buttons for 5 seconds. After hearing a beep, the pairing will be canceled. (At this time, the Bluetooth indicator light will flash for a while and then turn off)



Storage & Maintenance

1. Don't disassemble the device by yourself. It has no user-serviceable parts. When service or repairs are required, please contact customer service at service@quinear-ec.com.
2. Keep it out of the reach of infants, children and people without the ability to use it independently.
3. Don't place it in high temperature and humidity condition.
4. Avoid direct sunlight.
5. Avoid needles puncture the sleeves.
6. Don't place heavy stuff on it.
7. Keep the open ports of the control unit, and power in let free of debris.
8. Do not allow lint or dust to accumulate on the control unit. If lint or dust accumulates, wipe down the system with a dry cloth before use.

Cleaning

Make sure to cut off the power before you clean the device.

1. If dirty, please wipe the air pump controller, sleeves with a damp clean cloth.
2. Dry thoroughly with a clean cloth.
3. Don't use gasoline, alcohol, diluent and other irritating liquid to wipe the machine in case it causes malfunction or the components are damaged or discolored.
4. Don't dry clean.

Disposal

Please adhere to the local regulations when you dispose of the wasted

Package includes

- 2 x Massage Sleeves
- 1x User Manual
- 1 x Power Adapter
- 1 x Portable Storage Bag

Need Help? Contact Us!

For any issues during use, reach out to our dedicated support team anytime.

✉ Email: service@quinear-ec.com

Important:

- Include your order number and a brief description of the issue.
- Attach photos or videos (if applicable) for faster assistance.

We're here to help!

FCC Approved Equipment Authorization

FCC ID: 2BQV7-S9086

FCC Warning

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

NOTE 2: Any changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Labels

The following labels and symbols appear on the product and/or packaging.

Label	Description	Location
	Read instructions before use	On manual label
	Warnings, read instructions before use	On manual label
	Double Insulation. Equipment consists an extra layer of supplementary insulation over the basic insulation, or a single layer of special reinforced insulation directly over the live parts.	On manual label
	Shenzhen Dongjilian Medical Tech Co., Ltd. Building B1, No. 70 Tongxin Road, Baolong Street, Longgang District, Shenzhen, China	On manual label
	Separate collection for waste electrical and electronic equipment.	On package and product label
	In accordance with Directive 2014/35/EU electrical equipment designed for use within certain voltage limits, and Directive 2014/30/EU electromagnetic compatibility	On package and product label
	The electromagnetic radiation from the device is below the limits specified by the Federal Communications Commission.	On package and product label
	A directive in the European Union that regulates the use of certain hazardous substances in electrical and electronic equipment (EEE).	On package and product label
	For environmental protection, please recycle it instead of throwing it away as trash.	On package and product label