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AZLG®

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⑥

英文说明书 OraMemo Smart Ring User Manual Please read carefully before use and keep for future reference. ①	Package contents Smart Ring Charging Dock User Manual Charging Cable OraMemo Ring parameters Input: 5V DC 30mA Capacity: 3.92V DC 16.5mAh ②	Download the MemoHealth App iOS: Open the App Store on your device, search for "MemoHealth APP" and install it. Android: Open Google Play Store, search for "MemoHealth APP" and install it. First-time activation 1.Put the ring into the charging dock; it will begin charging automatically. 2.On its very first charge the ring flashes yellow once to show it has been activated. 3.While charging, the inner red LED on the ring blinks, and the indicator on the dock glows green. 4.We recommend fully charging the dock before first use to ensure it has enough power to activate and supply the ring. ③	Pair the ring 1.Turn on Bluetooth. 2.Open the MemoHealth APP → Me tab → Add Device. 3.Choose the ring; pairing completes automatically. ④	Wearing the ring: Slip the ring onto your finger, and for more accurate measurements, rotate it so the sensor faces the inside of the finger. The first time you use it, check the ring's battery level in the app. Charge it promptly if the level is low so you can wear it for extended periods. Charge the ring: 1.Place the ring into the dock and rotate it until it magnetically snaps onto the charging contacts. 2.Charging starts automatically; the ring's battery level can be viewed in Me → My Device in the app. ⑤	Here's what the indicator light on the charging case means while your ring is charging: Status ①: If the charging case shows a solid green light, the ring is currently charging. Status ②: If the green light is flashing, the charging case itself has low battery and needs to be charged. Status ③: Once the ring is fully charged, the green light on the charging case will turn off. ⑥	Charge the dock: While the dock itself is charging, the indicator blinks red. When the dock is fully charged, the red light stays on. Warm tip: If you won't be wearing the ring for a long time, charge it at least once every two months to maintain battery performance. ⑦	APP Download Customer support: If you have any questions, concerns, or need assistance regarding your product, please don't hesitate to reach out. You can contact Amazon Customer Service directly through your order page or email us at service@OraMemoai.com. OraMemo is committed to providing timely and helpful support for all your inquiries—we're here to help! ⑧	FCC STATEMENT: This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: ① This device may not cause harmful interference, and ② This device must accept any interference received, including interference that may cause undesired operation. Warning: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, ⑨	may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures: ● Reorient or relocate the receiving antenna. ● Increase the separation between the equipment and receiver. ● Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. ● Consult the dealer or an experienced radio/TV technician for help. RF warning statement: The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction. FCC ID: 2BQXQ-ORAWELL10 ⑩
--	--	---	---	---	--	---	--	--	--