2-BAND-100

User Manual

EN

LI

ES

FR

DE

IT

NL

SV

DA

PT

NO

Copyright © 2025 Samphire Neuroscience Ltd

User Manual English

Your Nettle

Your Nettle is a medical device. Please treat it like one. Its intended purpose is to reliably deliver transcranial direct current stimulation (tDCS) to the motor cortex and dorsolateral pre-frontal cortex to assist in the management of pain relating to menstruation and mood symptoms relating to premenstrual syndrome.

Appropriate use of Nettle will lead to the following clinical benefits:

- An improvement in pain related to menstruation, as measured by the pain VAS score on the Visual Analogue Scale (VAS); and
- An improvement in mood symptoms related to premenstrual syndrome, as measured by the Positive and Negative Affect Scale.

Like with all medical devices, the extent to which these benefits are experienced varies across individuals.

Nettle works by creating targeted and tightly controlled waves of energy that assist the natural communication between neurons in your brain. It is drug free and does not interfere with your hormones or natural cycle.

- 1 Button
- 2 Light
- 3 Charging Port
- 4 Electrodes

- | | |
|-----------------------------|--|
| 5 Medical Charger (LXCP12) | 9 Single Use Salt Sachets (SNSAL01) |
| 6 Solution Bottle (SNBOT01) | 10 Single Use Electrode Covers (SNCOV01) |
| 7 Towel (SNTOW01) | 11 Accessories Case (SNCAS01) |
| 8 Storage Pouch (SNPOU01) | 12 Replacements Pouch (SNREP01) |

You can order more single use consumables through the Samphire mobile application or website.

Before proceeding, check to make sure you have all the required accessories, and your Nettle has arrived to you intact and without damage. If something is missing or damaged, contact us.

The Samphire Mobile Application

Your Nettle is controlled via the Samphire Neuroscience mobile application. If you do not already have this application installed on your mobile device, follow the steps below.

EN

LI

ES

FR

DE

IT

NL

SV

DA

PT

NO

1. Install the Samphire Neuroscience mobile application on your mobile.

Search for "Samphire Neuroscience" in the relevant mobile application store. The application is compatible with iOS 16+ or Android 11+ devices.

2. Create a Samphire account through the mobile application.

You'll need a functioning email address and an internet connection.

When to use Nettle

In most cases, the symptoms your Nettle is intended to address are cyclical. Typically, Nettle works best when it is used prior to and during the time when you feel the most symptoms.

By default, use Nettle once per day, for five days, starting five days before the start of your symptoms or cycle. If your symptoms are severe, you can use Nettle more frequently, up to two times per day.

If your cycle or symptoms start unexpectedly, start and continue daily use until you feel relief, or seven days elapse.

You can use your Nettle ad hoc throughout your cycle if intermittent symptoms occur. There are daily and monthly limits on how often you can use your Nettle. If you exceed them, you'll be informed through the mobile application, and prevented from running a session.

The Samphire mobile application includes features such as menstrual cycle tracking and reminders, that may make using your Nettle more convenient, but they are optional.

Running your first Nettle session

Before running your first Nettle session, you'll need to complete a quick orientation through the Samphire Neuroscience mobile application. In the orientation, you'll learn about how your Nettle works, what it feels like, and other Nettle users' experiences.

When it's time to run your first Nettle session, tap to start the session in the Samphire Neuroscience mobile application.

You'll be guided through the entire process, with helpful videos and animations shown along the way. Here is a summary of what you will do.

Setting up for your Session

1. Gather required items

These are your Nettle (fully disconnected from its charger), four Electrode Covers, the Towel, and the Bottle. If the Bottle is empty, prepare more Salt Solution by adding one Salt Sachet to tap water, and shaking vigorously.

2. Prepare your Nettle

Insert one Electrode Cover into each of the four electrodes on your Nettle, and fully hydrate them using the Salt Solution.

3. Connect your Nettle

Press your Nettle's button, and connect to it from within the mobile application. You may need to allow permissions to access your mobile device's bluetooth functionality.

4. Put on your Nettle

It sits like a pair of headphones, angled forward slightly.

EN

LI

ES

FR

DE

IT

NL

SV

DA

PT

NO

Make sure your Nettle's Button is at the back of your head, just above and behind your right ear.

5. Verify your Nettle's placement and setup

The mobile application will run some quick checks to make sure your Nettle has good contact with your head. Depending on your hairstyle, some adjustments may be required. The mobile application will guide you.

During your session

1. Session starts

Your Nettle slowly increases to its target strength. In your first session, this happens slower than usual.

2. Session continues

A session lasts for twenty minutes. You can pause your session (for up to five minutes) by tapping the pause button in the mobile application.

3. Session ends

When time is almost up, your Nettle will slowly decrease its strength to zero, before ending the session.

After your Session

1. Nettle disconnects

Your Nettle will automatically disconnect and power off at the end of the session.

2. Dispose of the Electrode Covers and dry your Nettle

Electrode Covers are single use only. If there is moisture on your Nettle, dry it off before storing it away.

DO NOT reuse Electrode Covers. Doing so may increase the risk of skin irritation.

Cleaning your Nettle

Neither your Nettle nor its accessories are sterilised or need to be sterilised. After each session, you should wipe and dry your Nettle. When doing so, pay attention to the Electrode areas, and make sure they are carefully dried.

DO NOT clean your Nettle while it's charging.

Storing and transporting your Nettle

When you are not using your Nettle, you should store it in the provided carry bag, along with the provided accessories. Before storing it away, make sure to:

- dispose of the Electrode Covers (they are single use);
- clean and dry your Nettle; and
- keep your Nettle in a safe place, away from children and pets.

DO NOT store your Nettle in a way that compresses or stretches it away from its natural shape. Doing so may cause damage.

Future sessions

You will run future sessions the same way as you ran your first

EN

LI

ES

FR

DE

IT

NL

SV

DA

PT

NO

session. As you get more familiar with how your Nettle works, starting a session gets quicker and easier.

You are able to review the setup instructions at any time in the mobile application.

Charging your Nettle

You should not have to charge your Nettle more than once per month. Charge your Nettle if it is unresponsive, or if the indicator light flashes red.

To charge, connect your Nettle to the provided Medical Charger. Re-charging a fully flat battery takes about five hours. When the indicator light is white, your Nettle is fully charged.

If stored and used correctly, your Nettle's battery will last up to 500 charging cycles.

DO NOT connect the charger, or charge your Nettle, while it is on your head. Doing so could lead to electric shock.

DO NOT use a charger other than the provided charger. Doing so could damage your Nettle. The provided charger is a medical grade charger, and has special features to protect both the Device and you. If you lose your charger, contact us to obtain a replacement charger, using the provided contact details.

Precautions

Users should take the following precautions to ensure safe use of the device:

- Unless directed by a physician, do not use the device on children
- The safety of electrical stimulation during pregnancy has not been established.

- You may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (saline).
- If you have suspected or diagnosed heart disease, you should follow precautions recommended by your physician.
- Don't use your device on any body parts other than on the head, positioned as directed.
- Consult with your physician prior to using the device after a recent surgical procedure.
- Consult with your physician prior to using the device if you have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device.
- Use caution if stimulation is applied over areas of skin that lack normal sensations.
- Care for, clean, transport and store your device carefully and out of direct sunlight, and as directed.
- Keep your device out of the reach of children and pets.
- Use this device only with the accessories recommended by the manufacturer.

Warnings

- If you are under the active care of a physician, consult with your physician before using this device.
- Apply stimulation on the head, only as directed.
- Do not apply stimulation over open wounds or rashes, or swollen, red, irritated, infected, or inflamed areas or skin eruptions.
- Do not use in the bath, shower or other water bodies.
- Do not apply stimulation while driving, operating machinery, or during any activity in which you are otherwise at significant risk of injury or responsible for other people's safety.
- Do not use this device while intoxicated or incapacitated.

EN

LI

ES

FR

DE

IT

NL

SV

DA

PT

NO

- Do not use your device in areas with strong electric or magnetic fields.
- Do not use your device without inserting and hydrating the Electrode Covers.
- Do not use your device with dried or rehydrated Electrode Covers. They are single-use only.
- Do not charge your device with a third-party charger.
- Do not apply physical alterations to your device or attempt to repair any physical damage. In the case of any damage, immediately discontinue use and contact the manufacturer.
- Do not use the device while sleeping.
- Do not perform rigorous activities such as running or dancing while wearing the device. Gentle movements such as walking around or doing other tasks are OK, so long as the device does not dislodge from your head during operation.
- Do not allow children to play with or touch the charging cable of the Device. The charging cable has an associated risk of asphyxiation or strangulation.

Contraindications

There are important contraindications that apply to using Nettle. If any of the below contraindications apply to you, you should not use Nettle.

- Persons who are under 18 years old.
- Persons with a history of seizures or epilepsy.
- Persons who are pregnant or may be pregnant.
- Persons experiencing active suicidal thoughts.
- Persons with a pre-existing neurological or neuropsychiatric condition.
- Persons with a lesion, tumour or other defect in your skull, cranium or brain.
- Persons with an implant inside their skull, cochlear implant or implanted hearing aid.
- Persons with implanted medical devices, such as a cardiac pacemaker

or neurostimulation devices, such as spinal cord stimulators, vagal nerve stimulators, auricular stimulators, or deep-brain stimulating electrodes.

If you are unsure about whether any contraindications apply to you, you must consult with your physician before continuing, especially if you are under the care of a physician for any conditions related to the central nervous system (e.g. neurology, psychiatry).

Side effects and adverse reactions

No serious adverse events are expected with the use of Nettle. However, there are some minor side effects that may occur. These include skin irritation or redness, tingling, itching or burning sensations, tinnitus (noise or ringing in the ears) and transient headaches.

These typically resolve within an hour of use and are often the result of increased blood flow or brain activity. Should they become severe, uncomfortable, or painful, stop using your Nettle.

Side effects are less likely if the Electrode Covers are properly hydrated, and Nettle is firmly positioned on the head. The mobile application guides you through these things.

If you are unsure about whether any contraindications apply to you, you must consult with your physician before continuing, especially if you are under the care of a physician for any conditions related to the central nervous system (e.g. neurology, psychiatry).

Disposing of your Nettle

Dispose of your Nettle responsibly. Nettle contains a battery and other electrical components, so it should only be disposed of in accordance with the electrical waste disposal regulations that apply in your area.

EN

LI

ES

FR

DE

IT

NL

SV

DA

PT

NO

Technical specifications

Model Number	BAND-100
EU MDR Class	Class IIa (non-invasive)
Ingress Protection	IP42
Maximum Current Strength	2 mA $\pm$ 20%
Battery	Rechargeable LiPo 3.7 V, 400 mAh
Part Classification	Type BF Applied Part
Device Size	15 cm x 15cm x 10cm
Device Weight	185 grams
Individual Electrode Size	10.5 $\pm$ 0.5 cm ²
Operating Environment	Home healthcare environment
Operating Ranges	Temperature: +10°C to +40°C Humidity: 15% to 90%, non-condensing Altitude: 0-2000 metres Pressure: 700 hPa to 1060 hPa
Storage and Transport Ranges	Temperature: -20°C to +60°C Humidity: 10% to 90%, non-condensing
Indicator Light	Green (Solid) = Stimulating Green (Flashing) = Paused Blue (Solid) = Connected, Not Stimulating Blue (Flashing) = Waiting to Connect Red = Charging White = Fully Charged Purple = Charger inserted but not charging Red (Flashing before session) = Low Battery White (Flashing during session) = Contact error Red (Flashing during session) = Other error
Charging Parameters	5 Volts and 100 mA
Provided Medical Charger	Model LXCP12 (Rated Voltage 5 V; Rated Current 2.0 A; Maximum Power: 10.8 W, Minimum Power: 9.4 W, Protection Class: II, Duty: Continuous)

Symbol explanations

	Separate collection for electrical and electronic equipment (EC Directive 2012/19/EU, Waste Electrical and Electronic Equipment (WEEE)).
	Type BF Applied Part
	European Union mark of conformity
	Bluetooth
	Consult Instructions for Use (IFUs)
	Consult Accompanying Documents

How to Get in Touch

For any questions regarding your Nettle, or about Samphire Neuroscience more generally, you can contact **support@samphireneuro.com**

Samphire Neuroscience Ltd, 85 Great Portland Street, First Floor, W1W 7LT, London, U.K.

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FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

This device and its antenna(s) must not be co-located or operating in conjunction with any other antenna or transmitter.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction
FCC ID: 2BR5Q2BAND100

**Samphire
Neuroscience**

CE 1282

Samphire Neuroscience Ltd, 85 Great Portland Street, First
Floor, W1W 7LT, London, U.K.