

Thank you for purchasing
the Lucenia Dual Shine.

Thank you for purchasing
the Lucenia Dual Shine.

01	Dual Shine Overview	
	Product Components / Names	22-23
	LCD Screen Guide	24
	How to Use	24
	Function Guide	25
	Mode Sequence	25
02	Recommended Modes	
	Morning Mode	27
	Night Mode	28
	Scalp Mode	29
03	Basic Modes	
	RF Shot Mode	30
	Toning Mode	31
	MT Shot Mode	32
	V Shot Mode	33
04	Cautions / Warranty	
	Precautions During Use	34-35
	Precautions for Storage	36
	Product Warranty	37

Scan the QR code on the back
to view detailed instructions on how to use the product.

*Please read the manual and safety precautions carefully before use.
Images in this manual are simulated for illustrative purposes.

01

I.st LUCENIA DUAL SHINE

Dual Shine Overview

Product Components / Names

Main body, USB charging cable, Dedicated cradle, User manual



Head for Scalp

Brush-type electrodes deliver energy to the scalp gently without irritation.



Head for Face

Spiral electrodes adhere closely to the skin without lifting.



LCD Screen

LEVEL

Level Control Button

ON/MODE

Power / Mode Button

I.ST LUCENIA DUAL SHINE

Cradle



Charging Port **Rapid Charging**

Use the charging port located at the bottom of the device.

Cradle Charging Port **Standard Charging**

A dedicated cradle is included for both storage and charging.


You can safely store and charge the device by connecting the power cable directly to the cradle.

Charging Indicator




LCD Screen Guide

The intuitive LCD screen allows for easy, at-a-glance operation.



How to Use



Align the dot on the device with the dot on the head.

Twist the head to the left to lock it in place.


- STEP 1.** Press and hold the ON/MODE button to turn on the device.
- STEP 2.** Press the ON/MODE button briefly to switch modes.
- STEP 3.** Press the LEVEL button briefly to select your preferred level (1 to 5).
- STEP 4.** The device will begin automatically once a level is selected and the head is placed against your skin. If you lift the device from your skin, it will pause temporarily. The device will automatically turn off after 5 minutes.
- STEP 5.** To stop the current mode and reset the timer, press the ON/MODE button briefly.
- STEP 6.** To turn off the device from any state (active or standby), press and hold the ON/MODE button.

*The device will also turn off automatically once the preset care time is complete.

Function Guide

Recommended Modes	Morning	EP + Laser + Pink LED + IR / Duration: 3 min.
	Night	RF + Laser + Red LED + IR / Duration: 7 min.
	Scalp	MC + Red LED + IR + Vibration / Duration: 5 min.
Basic Modes	RF Shot	RF + Red LED + IR / Duration: 5 min.
	Toning	Laser + EP + Yellow LED + Vibration / Duration: 5 min.
	MT Shot	EP + Blue LED + Vibration / Duration: 5 min.
	V Shot	MC + Pink LED + Vibration / Duration: 5 min.

Mode Sequence




When the device is turned on, it starts in Morning Mode.
Modes run in the following order
: Morning > Night > RF Shot > Toning > MT Shot > V Shot



Scalp Mode starts automatically when the scalp head is attached.
*This mode does not activate when the skin head is attached.

LED Color Therapy for Revitalized Skin

Four LED wavelengths, each penetrating to different depths, enhance treatment effect and target various skin concerns.



	Blue	Skin soothing, Hydration
	Yellow	Brightening
	Red	Elasticity care
	IR	Texture care
	Pink	Firming, Facial contouring

Recommended Modes (Auto Care)

Effortlessly manage your morning, night, and scalp routines with this convenient auto care program.




Morning

Morning Mode

A quick program for busy mornings or before heading out. It helps smooth skin texture for better makeup application.

Functionality	LED	Time	Head Type
EP + Laser	Pink + IR	3min.	Face Head



[Pink + IR]



[Attach face head]

How to Use :
Apply an appropriate amount of skincare product to your skin. Gently press the head against the skin and roll outward from the center to aid absorption. Apply more product than you would during regular skincare. If your skin feels dry during use, apply additional product as needed.

Precautions Before Use :
Start at Level 1 when using the device for the first time. Gradually increase to a comfortable level over time.



Low intensity recommended for :

Sensitive areas around the eyes, nose, or mouth

Do not use on :

The throat area, wounds, or areas that have undergone procedures or surgery


*Avoid prolonged use on a single area and do not overuse, as it may cause irritation. Use at appropriate intervals.



Night
Night Mode

A focused care program ideal for evenings or when you have time to relax.

Functionality	LED	Time	Head Type
RF + Laser	Red + IR	7min.	Face Head



[Red + IR]



[Attach face head]

How to Use :
Apply a generous amount of skincare product to your skin.
Gently press the head against the skin and roll outward from the center and upward from the bottom, massaging in a lifting motion.
Apply more product than you would during regular skincare.
If your skin feels dry during use, apply additional product as needed.

Precautions Before Use :
Start at Level 1 when using the device for the first time. Gradually increase to a comfortable level over time.



Low intensity recommended for :

Sensitive areas around the eyes, nose, or mouth



Do not use on :

The throat area, wounds, or areas that have undergone procedures or surgery


*This program uses RF, which delivers heat to the skin. Avoid prolonged use on a single area and do not overuse, as it may cause irritation. Use at appropriate intervals.



Scalp
Scalp Mode

A program that helps massage the scalp and improve absorption of active ingredients.

Functionality	LED	Time	Head Type
MC + Vibration	Red + IR	5min.	Scalp Head



[Red + IR]



[Scalp care mode vibrations]



[Attach scalp head]

How to Use :
Apply an appropriate amount of scalp care product directly to the scalp. Gently press the head against the scalp and roll as you go.
Apply more product than you would during regular scalp care.
If your scalp feels dry during use, apply additional product as needed.

Precautions Before Use :
Start at Level 1 when using the device for the first time. Gradually increase to a comfortable level over time.



Do not use on :

Wounds or areas that have undergone procedures or surgery


*Avoid prolonged use on a single area and do not overuse, as it may cause irritation. Use at appropriate intervals.

**RF Shot****RF Shot Mode**

A mode designed for intensive elasticity care during the evening or when you have time to relax.

Functionality	LED	Time	Head Type
---------------	-----	------	-----------

RF	Red + IR	5min.	Face Head
----	----------	-------	-----------



[Red + IR]



[Attach face head]

How to Use :

Apply a generous amount of skincare product to your skin.
Gently press the head against the skin and roll outward from the center and upward from the bottom, massaging in a lifting motion.
Apply more product than you would during regular skincare.
If your skin feels dry during use, apply additional product as needed.

Precautions Before Use :

Start at Level 1 when using the device for the first time. Gradually increase to a comfortable level over time.



Low intensity recommended for :

Sensitive areas around the eyes, nose, or mouth



Do not use on :

The throat area, wounds, or areas that have undergone procedures or surgery


*This program uses RF, which delivers heat to the skin. Avoid prolonged use on a single area and do not overuse, as it may cause irritation. Use at appropriate intervals.

**Toning****Pigmentation & Dark Spot Laser Mode**

A mode that helps brighten and even out your skin tone.

Functionality	LED	Time	Head Type
---------------	-----	------	-----------

Laser + EP	Yellow	5min.	Face Head
------------	--------	-------	-----------



[Yellow]



[Attach face head]

How to Use :

Apply an appropriate amount of skincare product to your skin. Gently press the head against the skin and massage from the center outward.
Apply more product than you would during regular skincare.
If your skin feels dry during use, apply additional product as needed.

Precautions Before Use :

Start at Level 1 when using the device for the first time. Gradually increase to a comfortable level over time.



Low intensity recommended for :

Sensitive areas around the eyes, nose, or mouth



Do not use on :

The throat area, wounds, or areas that have undergone procedures or surgery


*Avoid prolonged use on a single area and do not overuse, as it may cause irritation. Use at appropriate intervals.



MT Shot
Booster Mode

A mode that helps improve absorption
when used with basic skincare product or sheet masks.

Functionality	LED	Time	Head Type
EP + Vibration	Blue	5min.	Face Head



[Blue]



[Booster mode vibrations]



[Attach face head]

How to Use :
Apply an appropriate amount of skincare product to your skin.
Gently press the head against the skin and roll outward from
the center to aid absorption.
Apply more product than you would during regular skincare.
If your skin feels dry during use, apply additional product as needed.

Precautions Before Use :
Start at Level 1 when using the device for the first time. Gradually increase
to a comfortable level over time.



Low intensity recommended for :
Sensitive areas around
the eyes, nose, or mouth



Do not use on :
The throat area, wounds,
or areas that have undergone
procedures or surgery


*Avoid prolonged use on a single area and do not overuse, as it may cause irritation.
Use at appropriate intervals.



V Shot
Microcurrent Mode

A mode that helps smooth and refine facial contours.

Functionality	LED	Time	Head Type
MC + Vibration	Pink	5min.	Face Head



[Pink]



[Microcurrent mode vibrations]



[Attach face head]

How to Use :
Apply an appropriate amount of skincare product to your skin.
Gently press the head against the skin and sweep upward,
focusing on target areas, such as the jawline.
Apply more product than you would during regular skincare.
If your skin feels dry during use, apply additional product as needed.

Precautions Before Use :
Start at Level 1 when using the device for the first time. Gradually increase
to a comfortable level over time.



Low intensity recommended for :
Sensitive areas around
the eyes, nose, or mouth



Do not use on :
The throat area, wounds,
or areas that have undergone
procedures or surgery

*Avoid prolonged use on a single area and do not overuse, as it may cause irritation.
Use at appropriate intervals.

Cautions / Warranty

Precautions During Use

1. Make sure to read the user manual carefully before use.
2. This product is intended for use on the face and designated areas of the neck only. The manufacturer is not responsible for any harmful results caused by use on other body parts or improper use.
3. Remove all metal accessories, such as necklaces, earrings, and rings, before use.
4. Do not use if you have a metal allergy or sensitive skin.
5. Do not drop the device or subject it to strong impacts during handling.
6. This product is not waterproof. You may clean the surface with a wet tissue, but do not use the device in water or rinse it under running water.
7. Do not use the product in humid environments or in water, such as in a bathtub. This may cause electric shock or malfunction.
8. Only use components provided by the manufacturer.
9. Do not modify, disassemble, or attempt to repair the product yourself.
10. You may experience slight tingling or pricking sensations during use. In this occurs, reduce the intensity. The sensations should ease or disappear shortly.
11. You may temporarily see flashing lights due to optic nerve stimulation during use. If this continues even when the device is not in use, discontinue use immediately.
12. If the product malfunctions or exhibits abnormal behavior, please contact the Service Center for inspection or repair.

Danger

Do not use the device together with the following types of medical electrical equipment under any circumstances.

Doing so may cause serious bodily harm due to malfunction or interference with the medical devices.

- Implantable medical devices (e.g., pacemakers)
- Life-support equipment (e.g., artificial heart-lung machines)
- Wearable medical equipment (e.g., electrocardiographs)
- External medical devices (e.g., insulin pumps)

Warning

Do not use the device on or near individuals with any of the following conditions :

1. Recent surgical scars or incisions
2. Vascular conditions (e.g., phlebitis, thrombophlebitis, capillary damage, varicose veins) or skin rashes, inflammation, infected, or swollen skin
3. Acute illnesses
4. Pregnancy or breastfeeding
5. Children or individuals under 18
6. Abnormal reactions to electrical stimulation
7. Metal allergies
8. Symptoms of capillary dilation
9. History of epilepsy or seizures
10. Bleeding from acute trauma or fractures
11. Heart disease or cardiovascular conditions
12. Multiple sclerosis
13. Cancer or cerebral infarction
14. Facial nerve damage
15. Thyroid disorders
16. Metal implants or internal medical devices
17. Active vasculitis

Precautions for Battery

The product contains a built-in battery. Before first use, ensure to fully charge the battery after purchase.

- 1. Batteries have a limited lifespan.
Total usage time may decrease with prolonged use.
- 2. Do not attempt to remove or replace the built-in battery.
- 3. If the battery needs to be removed or replaced, please contact the Service Center.

Precautions for Storage

- 1. Avoid storing the product in direct sunlight or in areas with high temperature or humidity.
- 2. After use, wipe the electrodes with a wet tissue followed by a dry cloth to remove any residue before storing.
(Do not clean with water or any liquid solution.)
- 3. Do not use strong chemicals (e.g., alcohol), as they may damage the surface of the product.
- 4. Store the device safely to prevent damage or breakage from external impact or pressure.
- 5. Keep out of reach of children.
- 6. Do not place volatile substances near the product.
- 7. Do not store or use the product in humid environments, under direct sunlight, or near heat sources.
- 8. To prevent excessive battery discharge, charge the battery at least once every three months. If the battery remains fully discharged for a long time, its lifespan may be reduced.
- 9. If the power does not turn on or shuts off immediately after turning on, the battery is likely discharged. Charge the device before use.
When fully discharged, charging may take longer than usual.

Product Disposal

When the device reaches the end of its lifespan, please dispose of it in accordance with local environmental regulations.

Product Warranty

Date of Purchase	
Buyer	
Seller	
Contact	

Product Warranty

- 1. The warranty period is one year from the date of purchase.
(The warranty period for the built-in battery is 6 months.)
- 2. If a manufacturing defect or natural malfunction occurs within the warranty period, the product will be repaired free of charge or compensated for in accordance with consumer compensation regulations.
- 3. The following cases will be subject to paid repairs, even during the warranty period:
 - Malfunctions caused by use outside the intended purpose or by improper handling
 - Damage or malfunctions resulting from repairs or modifications not performed by our company or an authorized Service Center
 - Malfunctions caused by using consumables or accessories not specified by our company (For designated consumables and accessories, please contact us.)
 - Malfunctions caused by natural disasters such as fire, salt damage, flooding, or other similar events
 - Cases involving user negligence, intentional damage, or failure to follow the instructions in the user manual

Product Specifications and Information

Product Name	I.st LUCENIA DUAL SHINE
Model Name	IST-D SHINE
Power Supply (Battery)	3.7 V, 3000 mAh
Power Consumption	Max. 15 W (Standby: Less than 1.5 W)
Dimensions and Weight	125(W) × 263(H) × 123(D) mm / 800 g (including cradle)
Manufacturer	Korea BTB Co., Ltd. 7 Dongwon-ro, Bundang-gu, Seongnam-si, Gyeonggi-do, 13547, Republic of Korea
Seller	Korea BNCCo., Ltd. 62 Seongseogongdan-ro 11-gil, Dalseo-gu, Daegu, 42713, Republic of Korea
Service Center Inquiry	070-7874-9399
Manufacturing Country	Made in Korea



This device is a cosmetic beauty device and is not classified as medical equipment.

FCC Information to User

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.