


mriYa®




Azure Cloud Mate06
Smart Scalp Serum Applicator

Product Introduction

This product is a scalp care device that can regularly dispense foamed liquid medicine and serum. It also has multiple functions such as EMS microcurrent, red light, blue light, red & blue light therapy, and micro-vibration massage, aiming to provide you with a comprehensive and convenient scalp care experience to care for your hair.

Product Structure



Function Introduction

1、Button Function Introduction

Button Name	Button Function
Power switch	Long press Power on/off Long press to enable the default factory mode (Red light + EMS level 1 + 1-minute liquid dispensing + low vibration)
	Short press Short press the power key once to retract the liquid compartment plunger to the initial position with one click, and reset the liquid dispensing output time. (Short press once when the plunger is in the middle position to pause the plunger without resetting the liquid dispensing time.)

2、Working Mode

Mode Name	Selectable Working Status
Light therapy mode	Red light / Blue light / Red & blue light / Off
EMS mode	Level 1-8 / Off
Vibration mode	Low / Medium / High / Off
Liquid dispensing mode	Dispensing time: 30-second dispensing / 1-minute dispensing / Off
Warm reminder: All the above working modes are set via the APP. Among them, the set working status of EMS, vibration, and liquid dispensing modes will only take effect when the device's touch sensing teeth are in contact with the scalp. The touch sensing status cannot be set.	

3、Indicator Light Introduction

(1) Power/charging indicator light:
After powering on, the power indicator light flashes for 15 seconds to prompt the user to connect via Bluetooth, and stays on after 15 seconds. The indicator light flashes when the device is charging and stays on when fully charged.
(2) Liquid dispensing mode status indicator light:
When the liquid dispensing mode is set in the APP: The indicator light stays on when the touch sensing teeth detect contact during a single dispensing period, and turns off when no contact is detected. The indicator light flashes while the liquid tank plunger is retracting, and turns off when retraction is complete. When no liquid dispensing mode is set: The indicator light is off.

4、Low battery reminder:

When the device is powered on and the battery is low, all lights will flash 5 times simultaneously. The device will automatically shut down after the reminder. Please charge it in time.

5、Smart power-off:

The device automatically shuts down after 5 minutes of standby. After turning on the device, if there is no APP control operation, no key operation, and no touch operation on the touch sensing teeth, it will shut down by default after 5 minutes.

6、Smart sensing:

(1) When the device is powered on, liquid dispensing, vibration, and EMS will only start when the touch sensing teeth are in contact with the scalp. To ensure the liquid reaches the scalp directly, liquid dispensing, vibration, and EMS will stop immediately if the touch sensing teeth leave the scalp or are blocked by hair.
(2) The effective sensing time for a single touch on the touch sensing teeth is about 15 seconds. After 15 seconds, you need to release the sensing teeth and touch again for the output to take effect again.

First-time use guide:

- 1.For Android or iOS phones, search for the "MY Smart Care" APP in the app market or App Store, then download and install it.
- 2.Before connecting, the APP will pop up a prompt to grant permissions. Click "OK" to grant the relevant permissions.
- 3.Turn on your phone's Bluetooth, power on the device, open the "MY Smart Care" APP, click "Add Smart Device", swipe the product list to find your device and click it. The APP will automatically connect to the selected care device.
- 4.After successful connection, the APP will automatically enter the product operation interface.

Usage Steps

- 1.Rotate the foam compartment cover counterclockwise and remove it.
- 2.Add an appropriate amount of foamed liquid medicine into the compartment, cover it, and tighten clockwise.
- 3.Turn on your phone's Bluetooth and open the "MY Smart Care" APP. Connect this device to your phone. Set your desired care mode via the APP interface.
- 4.Press the device head firmly against the skin. It is recommended to use continuous circular motions or slide at a constant speed.
- 5.After use, use the cleaning pin to remove residual liquid from the liquid outlet teeth, clean other parts, and store the device in a well-ventilated, dry place.

Precautions

- 1.When cleaning the device body and applicator head, do not use chemical detergents or cleaning products containing mercury, alkaline, or acidic ingredients.
- 2.Do not immerse the device in water or place it in humid environments (e.g., bathrooms) to prevent short circuits.
- 3.If liquid dispensing is not smooth during use, use the cleaning pin to unclog the liquid outlet teeth for smoother flow.
- 4.Do not use or store the device in high-temperature environments, as this may cause damage due to excessive heat.
- 5.If the device is not used for a long time, clean it regularly to maintain hygiene, performance, and lifespan. Keep the battery at around 50% charge and perform a full charge-discharge cycle periodically to maintain battery health.

- 6.Use only the original charging cable. Avoid using damaged cables and keep the device away from flammable and explosive materials while charging.
- 7.Do not disassemble, strike, or drop the device. Store it out of reach of children and pets.
- 8.When using serums of different specifications or brands, thoroughly clean the foam compartment to avoid interactions between different products.
- 9.This device is not a medical device and cannot be used to treat diseases. Consult a doctor before use if you have any medical conditions or health concerns.
- 10.Use only externally applied liquid medicines or serums. Do not use oral medications, corrosive liquids, flammable/explosive liquids, or highly oily liquids, as these may damage the device.
- 11.When using the EMS function for the first time, start at the lowest level and gradually adapt. Stop use immediately if you experience severe stinging or discomfort.
- 12.This device may receive irregular updates. Some functions may differ from this manual. Please refer to the actual functions of the product.

Warnings

The EMS electric pulse function is prohibited for the following groups or body parts:

- 1.People with heart disease, those wearing pacemakers or other implanted electronic devices.
- 2.Pregnant women, infants, young children, and people with epilepsy.
- 3.The neck, heart area, and skin areas that are broken, inflamed, or affected by eczema or dermatitis.

Package Contents

Intelligent scalp care device ×1, User manual ×1, Charging cable ×1, Cleaning pin ×1

Product Parameters

Product model	Xiaoqingyun mate06	Charging voltage	DC 5V	Operating voltage	DC 3.7V
Battery capacity	2000mAh	Foam compartment capacity	15ml	Product net weight	232g
Product dimensions	180×97.5×56mm				

FCC Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
-- Reorient or relocate the receiving antenna.
-- Increase the separation between the equipment and receiver.
-- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
-- Consult the dealer or an experienced radio/TV technician for help.
-- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
-- Children should be supervised not to play with the appliances.
-- This appliance contains batteries that are non-replaceable. When the battery is at end of life, the appliance shall be properly disposed of. This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions:
(1) this device may not cause harmful interference.
(2) this device must accept any interference received, including interference that may cause undesired operation.
FCC RF Exposure Warning Statement:
This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment can be used in portable exposure conditions.

Warranty card

Product model	
Product color	
Purchase address	
Customer name	
Purchase date	
Customer phone number	
Customer address	

Warm Tips

- 1.This product has a one-year warranty.
- 2.Use according to the precautions in the manual. It will be given a repair if the product fails within the warranty. (Please submit the product and warranty card when repairing.)
- 3.Within the warranty , a paid repair is required if:
(1)Failure and damage due to wrong use
(2)Failure and damage due to improper repair or modification by yourself.
(3)Failure and damage due to transportation, dropping, etc. after purchase.
(4)Failure and damage due to fire, earthquake, flood, electric shock and other weather abnormal voltage or caused by the specified external power source, etc.
- 4.This warranty card will not be reissued, please keep well.

brand owner
Mirya (Shenzhen) Technology co.,Ltd.
Manufacturer:
Shenzhen Kimairay Technology Co., Ltd.
Add.:Factory Building No. 5101, Zone 2, Fuqiao Industrial Zone,Heping Community, Fuhai Subdistrict, Bao'an District Shenzhen,China