

E013A-说明书

尺寸：152x214mm，材质：105g铜版纸

印刷：四色印刷+专色PANTONE Bright Orange C

成型：打钉，备注：先打样确认

此框内容不印刷

vitahunt

Professional 8 Electrodes Smart Body Composition Scale

Instruction Manual



Read this manual before using and save for future reference

THANK YOU FOR YOUR PURCHASE!

The Scale uses bio-electrical impedance analysis (BIA) technology to help you conveniently measure body weight, body fat rate, visceral fat, body water, skeletal muscle rate, muscle rate, muscle mass, bone mass, protein, BMR, body age and more.

Note: For more functional modules and personalized settings, please refer to the real-time interface of the APP.


WARNING

- DO NOT use the scale if you have implanted medical devices such as pacemakers.
- DO NOT stand on the edge of the scale or jump on it.
- DO NOT overload the scale (Max 396lb/180kg/28st).
- DO NOT drop the scale or drop objects on it as this may damage the sensors.
- DO NOT immerse the scale in water or use chemical cleaning agents. Clean the scale with a slightly damp cloth.
- Always place the scale on a hard, dry and flat surface before measurement.
- Make sure your feet are dry before stepping on the scale.
- DO NOT use the device if damaged. The continuous use of a damaged unit may cause injury or improper results.
- Please check the device before each use.
- Be careful when using on the wet and slippery surface.
- For people under 16 years old, the scale can only measure body weight.


PRODUCT SPECIFICATIONS

Product size:	300x328x36.5mm
TFT screen size:	78x53mm
Weight unit	kg/lb
Weight limit:	5kg-180kg /11lb-396lb
Weight division:	0.05kg/0.1lb
Power Supply:	4x AAA batteries
Operating temperature:	10-40°C
Operation humidity:	20%-90% RH

PRODUCT OVERVIEW



DISPLAY SCREEN



Note: The pictures shown in this manual are only for illustration. Please refer to the actual product for details.

BMITABLE

BMI					
Age	Underweight	Healthy weight	Overweight	Obesity	Severe obesity
18-100	<18.5	18.5-23.9	24-27.9	28-34.9	≥35

Body Mass Index (BMI)=Weight (kg)÷Height^2 (m)

For example: BMI=70kg÷ (1.75m*1.75m)=22.86 In the range of healthy weight

Note: The BMI index is a commonly used standard to measure body weight and health status.

TROUBLESHOOTING

Lo	Low battery - Please change new batteries
Err	Overload (Scale will show Err)
C	An error has occurred during measurement. Step off the scales and step onto it again to repeat the measuring process.


Use a DC power supply and adjust the voltage to below 4.0V. "Lo" will be displayed and the device will shut down after 3 seconds.



After powering on and returning to zero, load the weight to 183.0kg, "Err" will be displayed, and the machine will shut down after 3 seconds.




When the zero position of the weighing scale exceeds 0.5kg, "C" will be displayed as shown below, and it will return to zero after 3 seconds.



POWER ON

For first use

1. Press the battery cover with your thumb and push it up to remove the battery cover.
2. Insert the batteries. Be sure the polarity of the batteries is placed correctly. (If there are batteries in the the scale, please remove the insulating strip.)
3. Put the battery cover back and push it down to the bottom.



For subsequent use

Press to wake the body composition scale. The pressure should be bigger than the minimum range of the scale (5 kg).

Download and Install the App

Search Vitahunt in the App Store , or scan the QR code below to download and install the Vitahunt App.



Start the Measurement

- 1 Place the scale on a solid flat surface. Carpeted or uneven floors may affect the accuracy of the weight measurement .

Calibration: Gently step on the scale to wake it up and wait until the scale displays "0.00", which indicates that the calibration is completed (Remember to re-calibrate the scale if it has been moved).

- 2 Connecting to the Vitahunt App:
Turn on the Bluetooth of the mobile phone and open the app, Physical test score → enter the physical test page → measure weight

Note: Before using the body fat scale, you need to attach the Vitahunt ring first.

- 3 Correct posture for measurement

- ✓ Step on the scale with two bare feet. Ensure that the soles of the feet are in full contact with four metal electrodes and there is no contact between the left and right lower limbs (including thighs, calves and feet).





Your feet should not touch each other




Your feet should be in full contact with the electrodes



- ✓ Hold the handle with both hands. Ensure that the fingers evenly cover the four electrodes on the handle with no contact between the hands.



- ✓ Keep your arms straight at a 30-45 degree angle




Note: Please take off your shoes and socks, otherwise the scale can only measure your weight.

QUICK GUIDE

- 1 Install the battery and turn on the scale. After the full screen is displayed on the scale body and then the number 0.00 appears, you can place items on the scale for measurement.


Note: The Bluetooth symbol will not appear "Ⓔ" when the scale body is not connected to the mobile phone app.




- 2 After the weighing scale is connected to the "Vitahunt" App, the screen will display the current user information. When you step on the scale for measurement, after the weight value is locked, the trend indicator bar at the bottom will start to jump and enter the body fat measurement state. After the body fat measurement is completed, the values of weight, body fat, muscle mass, water content, and BMI will be displayed in sequence for 3 seconds each. After two rounds of display, the scale will automatically turn off.




1.The weighing scale recognizes the current user (P-1) in the App




2.Lock the user's weight




3.Display the user's height.
The jumping of the trend bar below indicates that the measurement is in progress




4.Show body fat rate



5.Show muscle rate




6.splay visceral fat



7.Show BMI

Offline operation

1. When the scale shows 0.00 after powering on, click the "set" button to enter the offline user setting mode. After entering, set the user information in the order of username, age, gender, and height.



1. Press the SET button lightly



2. Press and hold the SET button to set the age information.
[< Decrease, > Increase]



3. Press the SET button lightly to set the gender information.
[<> Select Male - Female]



4. Press the SET button lightly to set the height information.
[< Decrease, > Increase]



5. Press the SET button lightly.
Once the screen shows 0.00 after the setting is completed, you can step on the scale for measurement

2. In the user selection interface, use the left and right keys to cycle through and switch in the order of U-1, U-2, U-3, U-4, and P-0. Among them, U-1, U-2, U-3, and U-4 are fixed users, and P-0 is a quest user.



1. Press the weighing scale lightly, and the screen will display 0.00kg



2. Press the SET button lightly for the U-1 user



3. Press the > button lightly for U-1, U-2, U-3, U-4, and P-0 respectively



4. Press the > button lightly to cycle and switch in sequence

FCC STATEMENT

Operation is subject to the following two conditions.

(1) This device may not cause harmful interference.

(2) This device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

- Increase the separation between the equipment and receiver

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

- Consult the dealer or an experienced radio/TV technician for help

The users manual or instruction manual for an intentional or unintentional radiator shall caution the user that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. In cases where the manual is provided only in a form other than paper, such as on a computer disk or over the Internet, the information required by this section may be included in the manual in that alternative form, provided the user can reasonably be expected to have the capability to access information in that form.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

Email: service@vitahunt.com

EU	REP
----	-----

UK	REP
----	-----