## ***Intelligent Care Companion Device***

### **Quick Start Guide Guide de démarrage rapide Schnellstartanleitung**

For inquiries regarding Mbrainet, please visit <https://mbrainet.com/>.

@Shanghai Mbrainet Intelligent Technology Co., Ltd. reserves all rights.

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# ***MBR001***

<https://mbrainet.com/>

Thank you for purchasing this product. Please read this manual carefully before using it.

After reading, please keep this manual for future reference.

"Scan the code"  
to access more information.

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Love has never gone away; it's simply taken on a different form.

## 01. Charging requirements

- If the device is not used for a long period, please recharge it once every three months, fully charge it, then power it off and store it.
- For long-stored devices, the first use requires pre-charging for at least 30 minutes before the charging status appears.

### Names and Contents of Hazardous Substances in the Product

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This table is prepared in accordance with SJ/T 11364.  
 "○" indicates that the content of the hazardous substance in all materials of the part is below the limit requirements specified in GB/T 26572.  
 "×" indicates that the content of the hazardous substance in at least one homogeneous material of the part exceeds the limit requirements specified in GB/T 26572.

Shanghai Mabeiren Intelligent Technology Co., Ltd.

Mabrain official website: <https://mbrainet.com/>

Technical support e-mail: [support@mbrainet.com](mailto:support@mbrainet.com)

Made in China.

"Scan the code"  
to access more information.

## Precautions

Thank you for purchasing Mbrainet products! This user manual briefly introduces all the functions of the product, enabling you to operate it with ease and fully enjoy the great value brought by this product. Please read this user manual carefully before use.

- New batteries reach optimal performance only after 2–3 full charge and discharge cycles.
- This is a high-precision electronic product. Do not disassemble.
- Avoid strong impacts, drops, water, collisions, or splashes.
- Avoid use in high temperature, low temperature, high humidity, dusty, strong magnetic fields, or direct sunlight. Recommended use: altitude < 3000m, temperature 0°–35°C, humidity 5%–95%.

Use only manufacturer-approved batteries, and accessories. Use of other accessories violates warranty and may be dangerous.

## Friendly Reminder

- The nominal battery life of this product refers to the maximum duration under minimal power consumption. Actual battery life may vary due to usage and environmental factors, and there may be slight differences among individual products.
- The nominal battery capacity is calculated by measuring the charge during actual charging and discharging using third-party tools. Due to fluctuations in current, there can be significant deviations; this should not be used as a definitive standard.

## Button Definitions

- Power button: -----
- When the device is off, press and hold the power button for 3 seconds to power on.
- When the device is on, short press the power button to turn off the screen and enter sleep mode; to wake the screen, short press the power button again.
- When the device is on, press and hold the power button for 3 seconds to open the power-off/restart menu, or press and hold for 10 seconds to force a restart.

# 02. Welcome to your Alma health monitoring tablet


## Set up your Alma

When you turn on the device for the first time, select your preferred system language and time. If you need to change the language or adjust the time later, you can do so through the settings in Personal Center.


## Create a User

You can create a user profile for yourself, choose an avatar you like and set a name; family members or friends can also create new users. Each device supports up to five user profiles for household members.

## Delete a User


Press and hold the avatar and tap the red button **X** in the upper left corner of the avatar. Confirm whether you want to delete the user: tap Confirm to delete all data for that user, or tap Cancel to return to the user screen.

## Guest Mode



For temporary use, you can perform a single test through "Guest Mode".

**Temporary use mode:** In temporary use, the monitored items include only heart rate, blood pressure, heart rate variability, blood oxygen saturation, brain oxygen health, and stress levels. Data collected in these temporary measurements will not be saved to your health record.

## Prepare for a Test

Before starting, hold Alma or place it on a stand (recommended), adjust the position of your face and the tablet so that your entire face is within the frame. During testing, try to keep your body and head still and maintain about an arm's length between the tablet and your face. Refer to the diagram below.

Facial Positioning Diagram

To ensure accurate measurements, please use this product under the following lighting conditions:  
 Color temperature: between natural white light and daylight (about 4 000–6 500 K, neither yellowish nor bluish)  
 Color rendering: color rendering index (Ra) > 95 to accurately reproduce skin and object colors  
 Brightness: adequately bright (approximately equivalent to 500 lumens or more, such as a bright desk lamp or indoor daylight lamp)

Overexposed

Dim light

Yellow light

Incomplete face

Side face / Profile

Side face / Profile

Face too close


Face too far

Face out of frame

## Start Test

When sufficient ambient light is detected and your face is correctly positioned and stable, the "Start Test" button will turn green. Tap "Start Test" and wait for 20 seconds to obtain your current health metrics.

## Test Results

After completing the measurement, the results screen shows your current health status. You can scroll down to view additional metrics. If the results appear abnormal, please first check that the ambient lighting meets the usage requirements and that the device placement and the distance to your face are correct. Once verified, we recommend repeating the measurement.

Clicking each card rotates it to reveal the meaning of the metric on the back. Click again to return the card to its original state.


**Notes:**

**Depression risk:** An assessment result will be generated on the third day after you complete at least one health check each morning and evening.

**Sudden death risk:** An assessment result will be generated on the sixth day after you complete at least one health check each morning and evening.

**Physiological age and rate of ageing:** An assessment result will be generated on the ninth day after you complete at least one health check each morning and evening.

**Note:** The above risks and indicators are calculated based on continuous measurement data collected over multiple days; insufficient or intermittent measurements may affect the accuracy and stability of the results.

## Test Failed

If you do not follow the requirements, an abnormal detection prompt will appear. If the abnormal status persists for a long time, the system will indicate that the test has failed.

Make sure the lighting is sufficient, place your entire face inside the frame, and then tap to retest.

## Health Records

Each time you test, the results will be stored locally on Alma, forming long-term health data.


Tap Day, Week, Month or 6 Months to display the health trends for the day, week, month and past six months respectively.

You can use these long-term data to understand the changes in your health.

## Personal Center


Through the "Personal Center," you can manage your account, obtain help and perform various settings.

In Personal Center, tap the avatar to change your profile picture.


## Personal Center — Nickname

Through Personal Center—Nickname, you can change your display name.


## Personal Center — Gender

Through Personal Center—Gender, you can change your gender.


## Personal Center — Age

Through Personal Center—Age, you can change your date of birth.


## Personal Center — Time

Through Personal Center—Time, you can adjust the current date, time and time format.


## Personal Center — Help

Through Personal Center—Help, you can view the user guide.


## Personal Center — Language

Through Personal Center—Language, you can change the system language. Available languages include Chinese, English, Japanese, Korean, French and German.


## Personal Center — Feedback

Through Personal Center—Feedback, you can report any issues you encounter during use to customer service.


## Personal Center — About Almacarevue

Through Personal Center—About Almacarevue, you can view the product version number, privacy policy and product-related copyrights.

## Switch Account

If you want to switch users, you can tap the avatar on any screen and select your personal account.

## 03. Health Metric Descriptions

### AVI Score

The AVI score is a health quantification system developed by Mabrain and calculated precisely by intelligent algorithms.

You can find the meaning of each AVI score range in the table below. Enter your AVI score into the Health Card to easily compare your daily health status, and you can also view dietary recommendations via the Health Card.

| Level | Score range       | Level name       | Recommendation                                                                                                                                                                    |
|-------|-------------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| A 级   | $90 \leq x < 100$ | Excellent Health | Your health condition is excellent, keep up your good habits!                                                                                                                     |
| B 级   | $70 \leq x < 90$  | Good Health      | Your physical condition is good. Keep exercising and maintain a regular routine.                                                                                                  |
| C 级   | $50 \leq x < 70$  | Sub-health       | Some of your physical indicators are outside the ideal range. It is recommended to pay attention to your diet and rest.                                                           |
| D 级   | $30 \leq x < 50$  | Needs Attention  | Your physical condition may be in a potentially risky state. Regular monitoring is recommended.                                                                                   |
| E 级   | $0 \leq x < 30$   | High Risk        | Your physical condition has significantly deviated from the ideal range. It is recommended to improve your lifestyle as soon as possible and seek medical attention if necessary. |

### Biological Age

Biological age refers to the actual 'aging' state of an individual's physiological functions, organ systems, and physiological systems. It reflects the body's real health status and does not equal the age on your ID.

### Rate of Aging

The rate of aging is the speed at which physiological functions gradually decline from their optimal state, indicating how quickly a person transitions from being 'young and healthy' to functional decline.

0β-1β: Slow rate of aging

1β: Normal rate of aging

1β-2β: Fast rate of aging

## Heart Rate

Heart rate is the number of times the heart beats per minute and reflects the frequency of cardiac activity.

Enter your ideal heart rate value into the Health Card to easily compare your daily health status. You can also check dietary recommendations via the Health Card.

You can find the ideal heart rate values for each age group in Appendix 1.

## Heart Rate Variability

Heart rate variability (HRV) is the fluctuation of the intervals between heartbeats and reflects the autonomic nervous system's ability to regulate and the body's stress state.

A high HRV usually indicates strong recovery ability, flexible nervous system regulation and good physical and mental condition;

A low HRV may indicate high stress, fatigue, reduced immunity or a higher risk of chronic disease.

Enter your ideal HRV value into the Health Card to easily compare your daily health status, and you can also view dietary recommendations via the Health Card.

You can find the ideal HRV values for each age group in Appendix 2.

## Blood Pressure

Blood pressure is the pressure exerted by blood on vessel walls and is an important indicator of cardiovascular health.

Enter your ideal blood pressure values into the Health Card to easily compare your daily health status. You can also check dietary recommendations via the Health Card.

You can find the ideal blood pressure values for each age group in Appendices 3 and 4.

## Oxygen Saturation

Oxygen saturation is the degree of oxygen saturation in the blood and indicates whether the body is supplied with enough oxygen. Systolic pressure is the highest pressure experienced by the blood vessels when blood is pumped from the heart into the arteries during contraction.

Enter your ideal oxygen saturation value into the Health Card to easily compare your daily health status, and you can also view dietary recommendations via the Health Card.

## Brain Oxygen Health

Brain oxygen health indicates how well your brain is supplied with oxygen, which affects cognitive function and clarity of consciousness.

## Stress Level

Stress level indicates your current physiological or psychological stress state.

## Depression Risk

Depression risk reflects the likelihood of depressive symptoms or tendencies.

## Sudden Death Risk

Predicts the likelihood of experiencing a sudden cardiac or other unexpected fatal event.

## Cardiovascular and Cerebrovascular Health Risk and Blood Glucose Testing

Coming soon.

## Appendix 1

### Heart Rate

| Age group | Ideal value - male | Ideal value - female | Tolerance |
|-----------|--------------------|----------------------|-----------|
| 18-19     | 65                 | 68                   | ±5        |
| 20-24     | 66                 | 69                   | ±5        |
| 25-29     | 66                 | 69                   | ±5        |
| 30-34     | 67                 | 70                   | ±5        |
| 35-39     | 68                 | 71                   | ±6        |
| 40-44     | 69                 | 72                   | ±7        |
| 45-49     | 70                 | 73                   | ±8        |
| 50-54     | 71                 | 74                   | ±9        |
| 55-59     | 72                 | 75                   | ±10       |
| 60-64     | 73                 | 76                   | ±11       |
| 65-69     | 74                 | 77                   | ±12       |
| 70-74     | 75                 | 78                   | ±13       |
| 75-79     | 76                 | 79                   | ±14       |
| 80+       | 77                 | 80                   | ±15       |

#### Disclaimer

The ideal values are taken from publicly available scientific literature and authoritative databases and are intended to provide general health references for users.

Please note: these reference ranges do not constitute professional medical diagnostic standards and cannot replace a doctor's professional advice. There are significant physiological differences between individuals, influenced by factors such as age, gender, physical condition, lifestyle habits, and medical history.

## Appendix 2

### Heart Rate Variability

| Age group | Ideal value | Tolerance |
|-----------|-------------|-----------|
| 18-24 岁   | 160         | ±35       |
| 25-29 岁   | 155         | ±35       |
| 30-34 岁   | 150         | ±36       |
| 35-39 岁   | 145         | ±36       |
| 40-44 岁   | 140         | ±37       |
| 45-49 岁   | 135         | ±37       |
| 50-54 岁   | 130         | ±38       |
| 55-59 岁   | 125         | ±39       |
| 60-64 岁   | 120         | ±40       |
| 65-69 岁   | 115         | ±40       |
| 70-74 岁   | 110         | ±42       |
| 75-79 岁   | 105         | ±44       |
| 80+ 岁     | 100         | ±45       |

#### Disclaimer

The ideal values are taken from publicly available scientific literature and authoritative databases and are intended to provide general health references for users.

Please note: these reference ranges do not constitute professional medical diagnostic standards and cannot replace a doctor's professional advice. There are significant physiological differences between individuals, influenced by factors such as age, gender, physical condition, lifestyle habits, and medical history.

## Appendix 3

### Systolic Pressure (High Pressure)

| Age group | Ideal value | Tolerance |
|-----------|-------------|-----------|
| 18–19     | 108         | ±6        |
| 20–24     | 110         | ±7        |
| 25–29     | 111         | ±7        |
| 30–34     | 112         | ±7        |
| 35–39     | 115         | ±8        |
| 40–44     | 117         | ±8        |
| 45–49     | 120         | ±9        |
| 50–54     | 124         | ±9        |
| 55–59     | 128         | ±10       |
| 60–64     | 132         | ±10       |
| 65–69     | 136         | ±11       |
| 70–74     | 140         | ±11       |
| 75–79     | 143         | ±12       |
| 80+       | 145         | ±12       |

#### 1 Disclaimer

The ideal values are taken from publicly available scientific literature and authoritative databases and are intended to provide general health references for users.

Please note: these reference ranges do not constitute professional medical diagnostic standards and cannot replace a doctor's professional advice. There are significant physiological differences between individuals, influenced by factors such as age, gender, physical condition, lifestyle habits, and medical history.

## Appendix 4

### Diastolic Pressure (Low Pressure)

| Age group | Ideal value | Tolerance |
|-----------|-------------|-----------|
| 18–19     | 70          | ±6        |
| 20–24     | 71          | ±6        |
| 25–29     | 72          | ±6        |
| 30–34     | 72          | ±6        |
| 35–39     | 73          | ±6        |
| 40–44     | 74          | ±7        |
| 45–49     | 75          | ±8        |
| 50–54     | 76          | ±8        |
| 55–59     | 75          | ±9        |
| 60–64     | 74          | ±9        |
| 65–69     | 73          | ±9        |
| 70–74     | 72          | ±9        |
| 75–79     | 71          | ±9        |
| 80+       | 70          | ±9        |

#### 1 Disclaimer

The ideal values are taken from publicly available scientific literature and authoritative databases and are intended to provide general health references for users.

Please note: these reference ranges do not constitute professional medical diagnostic standards and cannot replace a doctor's professional advice. There are significant physiological differences between individuals, influenced by factors such as age, gender, physical condition, lifestyle habits, and medical history.

## FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harm full interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Specific Absorption Rate (SAR) information:

This Device meets the government's requirements for exposure to radio waves. The guidelines are based on standards that were developed by independent scientific organizations through periodic and thorough evaluation of scientific studies. The standards include a substantial safety margin designed to assure the safety of all persons regardless of age or health. FCC RF Exposure Information and Statement the SAR limit of USA (FCC) is 1.6 W/kg averaged over one gram of tissue. Device types: Device has also been tested against this SAR limit. This device was tested for typical body-worn operations with the back of the device kept 0mm from the body. To maintain compliance with FCC RF exposure requirements, use accessories that maintain an 0mm separation distance between the user's body and the back of the device. The use of belt clips, holsters and similar accessories should not contain metallic components in its assembly. The use of accessories that do not satisfy these requirements may not comply with FCC RF exposure requirements, and should be avoided.

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## ISED Statement

- English: This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

(1) This device may not cause interference.

(2) This device must accept any interference, including interference that may cause undesired operation of the device.

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This device meets the government's requirements for exposure to radio waves. The guidelines are based on standards that were developed by independent scientific organizations through periodic and thorough evaluation of scientific studies. The standards include a substantial safety margin designed to assure the safety of all persons regardless of age or health. ISED RF Exposure Information and Statement the SAR limit of and Canada (ISED) is 1.6 W/kg averaged over one gram of tissue. Device types: device has also been tested against this SAR limit. This device was tested for typical body-worn operations with the back of the device kept 0mm from the body. To maintain compliance with ISED RF exposure requirements, use accessories that maintain an 0mm separation distance between the user's body and the back of the device. The use of belt clips, holsters and similar accessories should not contain metallic components in its assembly. The use of accessories that do not satisfy these requirements may not comply with ISED RF exposure requirements, and should be avoided.

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The image shows two screenshots of a mobile application interface. The left screenshot, titled 'Persönliches Zentrum', displays a user profile card with a male avatar and the following details: Spitzname: test, Geschlecht: Frau, Alter: 1994-06-02, Zeit, Hilfe, Sprache, Feedback, and Über Almacarevue. Below the card is a '< Zurück' button. The right screenshot, titled 'Spitznamen bearbeiten', shows a text input field containing 'test' and a green 'Beenden' button. It also has a '< Zurück' button at the bottom.

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The image shows two screenshots of a mobile application interface. The left screenshot, titled 'Persönliches Zentrum', displays a user profile card with a male avatar and the following details: Spitzname: test, Geschlecht: Frau, Alter: 1994-06-02, Zeit, Hilfe, Sprache, Feedback, and Über Almacarevue. Below the card is a '< Zurück' button. The right screenshot, titled 'Geschlecht festlegen', shows two buttons: 'Männlich' (blue) and 'Weiblich' (red), with a green 'Beenden' button below them. It also has a '< Zurück' button at the bottom.

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The image shows two screenshots of a mobile application interface. The left screenshot, titled 'Persönliches Zentrum', displays a user profile card with a male avatar and the following details: Spitzname: test, Geschlecht: Frau, Alter: 1994-06-02, Zeit, Hilfe, Sprache, Feedback, and Über Almacarevue. Below the card is a '< Zurück' button. The right screenshot, titled 'Geburtsdatum bearbeiten', shows a date picker with the year 1994 selected, and the month 6 and day 2. Below the date picker is a green 'Beenden' button. It also has a '< Zurück' button at the bottom.

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