

# HEATING WAVE VIBRATION ELECTRIC DR.MEETI CERVICAL PILLOW

## ● Combination contents and names of each part


## ● Charging method

- 1 Insert the included USB cable into the USB Type-C port next to the device.
- 2 Connect the USB cable terminal to the USB adapter or mobile battery. \*Please prepare USB adapter and mobile battery separately. More than 5v 1a is recommended.
- 3 The charging signal red light flashes, and the green light is on when fully charged


## ● Turn off the power

- 1 Long press the power button (press for about 3 seconds).

After turning on, it can be used with a remote control, and M is the function button

## ● Instructions

### ◆ How to use

- 1 Press and hold the power for about 3 seconds to turn on the machine. The mode indicator light is on the first grid. You can press the increase button. When the number 1 is displayed in the gear area, the EMS function is turned on. At this time, you can press the increase or decrease button to select The required EMS gear (6 gears in total), when the gear area is number 0, the EMS function stops using.


- 2 Press the power button again. When the mode light is in the second grid, it is in heating mode. When the gear area is number 1, the heating mode is turned on. There are 3 levels in total. Warning: (Everyone chooses according to their own needs, beware of burns)
- 3 Press the power button again. When the mode light is in the third grid, the vibration mode is enabled. When the gear area is number 1, the vibration mode is turned on. There are 3 levels in total.
- 4 The power will be automatically turned off after 15 minutes. About using the heating function Do not use continuously for a long time. At lower temperatures (40-50°C), long term Even if there are no subjective symptoms such as pain, contact with skin contact parts may cause low-temperature burns. Please do not do anything in contact with water (cooking, washing, pets, etc.) while using the EMS.
- 5 EMS/heating/vibration functions can be used at the same time or separately.

**Attention: This product does not have medical functions for use.**

## ● Instructions

Please adjust the position so that the hot EMS plate touches the place you want to heat. As shown in the figure below, this product can be used regardless of the orientation. Actually try it on your neck and use it in the direction you like.  
This product is designed for the neck. Please do not use on other parts such as arms and legs.

**EMS is a structure that transmits stimuli through water and moisture.**  
**Wet the neck where the pad touches with water or lukewarm water before use.**


## ● Maintenance

- Please use and maintain the warm EMS board every time.
- Wring out the cloth soaked in water fully, warm it up
  - Wipe the EMS board well.
  - If you are concerned about dirt, soak a soft cloth with neutral detergent diluted in warm water, squeeze it dry, and then wipe off the dirt on the surface of the fuselage.
  - Please do not soak in water or wash all over. Do not use solvents such as alcohol, thinner, or benzene. Cause failure, damage to parts, discoloration, etc.

### FCC Caution:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.  
NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.  
If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:  
-- Reorient or relocate the receiving antenna.  
-- Increase the separation between the equipment and receiver.  
-- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.  
-- Consult the dealer or an experienced radio/TV technician for help.  
The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

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



## English Version

## How to use

1. Press and hold the power button for about 3 seconds to turn on the power, the mode light is in the first frame, you can press the increase button, when the number 1 is shown on the display, it indicates that the EMS function is turned on, at this time, you can press the increase or decrease button to select the desired EMS gear (6 speeds in total), the EMS function is stopped when the number 0 is shown on the gear area.
2. Press the power button again, the mode light is in the second frame for heating mode, when the block area is number 1, the heating mode will be turned on, a total of 3 speeds.
- Warning: (each person to choose according to their own needs, beware of scalding)
3. Press the power button again, the mode light is in the third grid for vibration mode, when the block area is number 1, turn on the vibration mode, a total of 3 gears.
4. Automatically turn off the power after 15 minutes. About the use of the heating function please do not use it continuously for a long time. At lower temperatures (40-50°C), prolonged use may result in low-temperature burns from contact with the skin contact point even if there are no subjective symptoms such as pain. Do not do anything that comes into contact with water (cooking, washing, pets, etc.) while the EMS is in use.
5. EMS/Heat/Vibration function can be used simultaneously or separately.

\*Note: This product is not for medical function use.