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95mm

95mm

95mm

95mm

95mm

95mm

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睡眠 B-Relaxer

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User's Manual

B-Relaxer

Please read this user's Manual carefully and retain it for future reference before operating the device.


1. Product Introduction

This product utilizes brain-computer intelligent chip technology, combining specific frequencies of audio, energy waves, and light waves to stabilize emotions, relieve stress, and improve sleep quality.

2. Package Contents

Main Unit x1	Type-C Charging Cable x1
Elastic Headband x1	Rear Headband Pad x1
Forehead Pad x1	Replacement Transducer Pad x2
Storage Case x1	Cleaning Cloth x1
Storage Bag x1	User's Manual x1

3. Product Structure Description

1. Power Touch Area	2. Light Effect Touch Area
3. Front Cover Light	4. Bottom Light Strip
5. Forehead Pad	6. Working Transducer
7. Replacement Transducer Pad	8. Type-C Charging Port
9. Elastic Headband	

4. Product Specifications

Version: X-N001	Frequency Response Range: 20Hz-10KHz
Playback Time: Approx. 3 hours	Standby Time: Approx. 90 minutes
Charging Time: Approx. 2 hours	Battery Specification: 1000 mAh
Operating Power: 10W (5W × 2)	Weight: 250 g

※The above data are laboratory-measured values. Actual performance may vary depending on usage conditions. Playback time is measured under 80% of device performance.

5. Operating Instructions

(1) **Power On**: Press and hold the power touch area for 1.5 seconds to turn on the device. The light strip will flash orange to indicate successful startup.

(2) **Connection**: After powering on, the light strip will blink slowly, indicating the device is in Bluetooth pairing mode. Use a Bluetooth-enabled phone to search for and connect to **B-Relaxer**. If no connection is made within three minutes, the device will automatically power off.

(3) **Power Off**: Press and hold the power touch area for 2 seconds or insert the charging cable. The light strip will flash quickly three times with a shutdown prompt tone to confirm power-off.

6. Light Effect Control

(1) After initial startup, lightly touch to select from three modes: light off, steady on, or breathing mode. The device will default to the last used mode in subsequent uses.

(2) **Breathing Mode**: Users can slowly adjust their breathing rhythm according to the specific flashing frequency of the light strip.

(3) **Low Battery Reminder**: The light strip will flash red rapidly accompanied by a voice prompt, indicating that the device needs to be charged.

7. Wearing Instructions

- Place the product on your forehead with the charging port facing downward. Adjust the elastic headband according to your head size.
- Once the headband is properly adjusted, position the main unit so that the working transducer fits snugly against your forehead.
- The rear pad on the elastic headband can be removed based on personal preference.

8. Playback Function Operation

Once the device is successfully connected to a phone via Bluetooth, download the **B-Relaxer APP** by following the instructions in the APP installation Guide to connect and use the device.

Control Functions in Connected Mode:

- Pause/Play**: Lightly touch the power touch area to pause playback. Touch again to resume.
- Session Duration**: The device will automatically power off after 30 minutes of continuous use, accompanied by a prompt tone. If paused before 30 minutes, the device will automatically shut down after a cumulative 90 minutes of inactivity.
- For APP installation, please refer to the attached **APP Installation Guide**.

9. Charging and Battery Precautions

(1) Charging the Device: Insert the Type-C charging cable. The charging port will glow red during charging, and the light will turn off when the battery is fully charged.

(2) Battery Precautions:

- This product contains a built-in rechargeable battery that cannot be removed.
- The charging current and voltage should be 5V/1000mA. Do not use a charger exceeding 5V/1000mA, or the battery might be damaged.

10. Maintenance and Care

- This is a wearable product. Before use, clean it with water or medical alcohol by the cleaning cloth included.
- The forehead-contact area of the product is fitted with shock-absorbing soft pads, which may wear out over time. Replacement pads are included in the package and additional ones can be purchased separately through our official customer service.
- The elastic headband is secured with buckles for easy installation and removal. It can be cleaned with a damp cloth and must be fully dried before reuse.

11. Safety Precautions

- To ensure personal safety, do not wear or place this product too close to a cardiac pacemaker (less than 10 cm).
- This product is not a toy and is intended for users aged 16 and above.
- Before turning on the device, set your phone's volume to the lowest level to avoid discomfort caused by high volume.
- For safety reasons, do not use the device while driving or cycling.
- Please do not disassemble, impact, or bend the device with excessive force.
- Please dispose of the product in an environmentally responsible manner. Please do not discard indiscriminately.

12. Names and Contents of Toxic or Hazardous Substances or Elements in the Product

This table is prepared in accordance with the provisions of SJ/T 11364-2014:

Component Name	Toxic or Hazardous Substances or Elements					
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent Chromium Cr (VI)	Polybrominated Biphenyls (PBBs)	Polybrominated Diphenyl Ethers (PBDEs)
Metal Parts	×	○	○	○	○	○
Plastic Parts	○	○	○	○	○	○
Wires	○	○	○	○	○	○
Circuit Board Assembly	×	○	○	○	○	○
Battery	×	○	○	○	○	○

※ Indicates that the content of the toxic or hazardous substance in all homogeneous materials of this part is below the limit specified in GB/T 26572.

※※ Indicates that the content of the toxic or hazardous substance in at least one homogeneous material of this part exceeds the limit specified in GB/T 26572.

The parts marked "x" in the table cannot yet be substituted with harmless alternatives due to global technological limitations.

13. Contact Information

Email: xncurs@outlook.com

Customer Service WeChat: B-Relaxer004

14. Product Certificate

Inspection Date:

Inspector:

Inspection Conclusion: ★ Qualified ○ Not Qualified

Company Name: Shanghai Guang'an Medical Technology Development Co., Ltd.

Address: Floor 26, No. 588-1 Longchang Road, Yangpu District, Shanghai, China

FCC Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

RF Exposure Information

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.