

Instructions(with remote control)

(Different styles/functions correspond to different buttons)

Single heating:



Red high /Blue medium/Green low



Long press (on/off)
Short press to shift

Massage heating:



Red high /Blue medium/Green low



Red high /Blue medium/Green low



Long press (on/off)
Short press to shift

Plug

Long press (on/off)
Short press to shift

Massage heating airbag:



Vibration→Charge and deflate→heat

Blue vibration/Green charge and deflate/Red heating



Long press (on/off)
Short press to shift

plug

Long press (On/Off)
Short press to switch



Precautions for use

1. This product is not suitable for children under 10 years old, pregnant women, lactating women, patients with hypertension and diabetes, people without self-care ability and people with skin allergies.
2. Please do not use this product after applying painkillers, laxatives, or other products.
3. Due to individual skin sensitivity differences, if you feel high temperature or other discomfort during use, please stop using immediately; If there are common symptoms such as rash, redness, and itching on the skin, please seek medical attention immediately.
4. The temperature of this product may vary depending on the usage environment.
5. Do not use this product together with other heating devices.
6. This product is a hot compress product and cannot replace medication.
7. It is recommended to use this product within 30 minutes to avoid skin burns/excessive massage caused by prolonged use.
8. This product only supports 5V voltage plugs with 3C certification for power on use. Without 3C certification or plugs with 12V or 24V voltage, there is a risk and their use is strictly prohibited.
9. Product details: Product+power cord (without charging head)
10. Cannot be washed.

Manufacturer: Guangzhou Shure Environmental Protection Technology Co., Ltd
Tel: 19898396007
Implementation standard: GB 18401-2010 Standard value
Add: No. 357, Shilou Road, Shilou Town, Panyu District, Guangzhou City (temporary house number), Building 22, No. 10 factory in the industrial zone

FCC Caution

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radiocommunications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

- Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for help.

*RF warning for Portable device: The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

