

LE-CARE
Smart Pill Organizer
PD01 Bluetooth Version
User Manual



Please read the instruction carefully before using this product and keep it properly.

Contents	01 Before Use	02 Package Contents	03 Product Overview	04 First-Time Setup	04 First-Time Setup	04 First-Time Setup	04 First-Time Setup
01 Before Use 02 Package Contents 03 Product Overview 04 First-Time Setup 05 Daily Use 06 Battery level 07 Care and Maintenance 08 Specifications 09 Compliance Statements 10 Support	Please read this manual before using the product. ■ This product is designed to assist with medication reminders and medication record management. ■ It does not replace medical diagnosis, prescriptions, or professional medical advice. ■ Always follow your physician's instructions and medication labels. ■ For children, older adults, or users with special medical conditions, setup and regular checking by a caregiver is recommended. ■ Before first use, charge the device, connect it to the applet, and set up medication schedules. *Records are based on compartment opening behavior and are not a medical confirmation of actual ingestion.	Smart Pill Organizer x 1 Data cable x 1 User Manual x 1	① Function button ② Indicator Light ③ Charging port ④ Speaker	Step 1: Download the APP Open the app store and search for "APP", or scan the QR code to download. 	Step 2: Power on Make sure the device has sufficient battery power. Press and hold the power button of the medicine box for 3 seconds. 	Step 3: Connect the device Method 1: Log in to the APP, click "Scan Code", and aim the scanner at the linear barcode on the back of the medicine box for scanning and recognition. *Before connecting: Turn on Bluetooth on your phone. Keep the phone close to the device. 	Step 3: Connect the device Method 2: Manually enter the product code. The product code is on the back of the medicine box. 

First-Time Setup

04 First-Time Setup

05 Daily Use

06 Battery level

07 Care and Maintenance

08 Specifications

09 Compliance Statements

09 Compliance Statements

10 Support

Step 4:
Set the plan

Set the drug name, dosage, usage, reminder time, etc.

1

2

3

4

5

6

7

8

9

10

11

12

Select Time

Range: 0:00-24:00

Select Compartment

Range: 1-12

Repeat Interval

Range: 1-120 min

Rem

Reminder cycle at 10 min

*For medications requiring moisture-proof, light-proof, refrigerated, or special storage conditions, follow the medication label or professional instructions.

Step 5:
Put in the medicine

According to the pre-set plan, place the medicine in the corresponding medicine compartment by category.

1

2

3

4

5

6

7

8

9

10

11

12

Select Time

Range: 0:00-24:00

Select Compartment

Range: 1-12

Repeat Interval

Range: 1-120 min

Rem

Reminder cycle at 10 min

*For medications requiring moisture-proof, light-proof, refrigerated, or special storage conditions, follow the medication label or professional instructions.

At the scheduled time, the medicine box will emit a sound and light up.

Open the blinking lighted medicine container and take the medicine as per the plan.

1

2

3

4

5

6

7

8

9

10

11

12

Select Time

Range: 0:00-24:00

Select Compartment

Range: 1-12

Repeat Interval

Range: 1-120 min

Rem

Reminder cycle at 10 min

*For medications requiring moisture-proof, light-proof, refrigerated, or special storage conditions, follow the medication label or professional instructions.

Battery is fully charged

The battery level is below 10%.

The medicine box contains lithium batteries. After charging once, the standby time is approximately 60 days.

When the battery power is low, please recharge promptly to avoid affecting normal use.

1

2

3

4

5

6

7

8

9

10

11

12

Select Time

Range: 0:00-24:00

Select Compartment

Range: 1-12

Repeat Interval

Range: 1-120 min

Rem

Reminder cycle at 10 min

*The daily usage battery life is affected by the reminder frequency and working mode. Please refer to actual usage for details.

■ The product should be stored in an environment that meets the storage requirements.
Storage conditions: Temperature -10°C to +45°C; Relative humidity < 80% R.H. Protected from light, well-ventilated, without corrosive gases.

■ Daily cleaning can be done by wiping the equipment with a dry soft cloth or alcohol cotton swabs, then drying it before use.

■ Do not wash, soak, or place the device in a dishwasher.

■ Do not use corrosive cleaners.

■ Keep it away from heating devices and open flames.

■ Avoid dropping, squeezing, or strong impact.

■ If the device is not used for a long time, charge it regularly and store it in a dry place.

Name

Smart Pill Organizer

Model

PD01 Bluetooth Edition

Dimension

248 x 160 x 47 mm

Weight

413g

Battery Capacity

2000 mAh

Power Input

5V

Rated Power Consumption

1 mW

Connection

Bluetooth

Supported OS

Android & iOS

■ **Child Safety Statement**
Keep this device and all medication out of reach of children. This product is not child-resistant.

■ **Hazardous Safety Statement**
Use only the specified charging input. Do not use if the device, cable, or battery is damaged, wet, overheated, swollen, or otherwise abnormal.

■ **Medical Disclaimer**
This product is a medication reminder tool and is not a medical device. It is intended to assist users in managing their medication schedules but does not guarantee that doses will not be missed. Users should always verify their medication intake. Consult your physician or pharmacist before making changes to any medication regimen.

■ **Environmental and Recycling Statement**
This product should not be disposed of along with household waste. Please follow the local environmental laws and guidelines for disposal. Protect the environment. Before discarding the product, please remove the battery and place it at a designated recycling point. Do not disassemble, discard or burn it at will.

■ **FCC Caution**
This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

■ **RFCA**
This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
-Reorient or relocate the receiving antenna.
-Increase the separation between the equipment and receiver.
-Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
-Consult the dealer or an experienced radio/TV technician for help. The device has been evaluated to meet general RF exposure requirements. The device can be used in portable exposure condition without restriction.

If you have any questions, please contact us.
Wenn Sie Fragen haben, kontaktieren Sie uns bitte.
Email: service@le-care.com