
DEX™

Accessories


EN
DE
FR
ES
DK
IT
JP
US

WIDEX®
HIGH DEFINITION HEARING


1

2

3

CONTENTS

DESCRIPTION OF DEVICE.....	8
Intended use.....	8
Accessories.....	9
CHARGING YOUR DEVICE.....	10
TURNING THE DEVICE ON AND OFF.....	11
PAIRING YOUR DEVICE.....	12
USING YOUR DEVICE.....	13
Sound adjustment.....	13
Visual indicators.....	13
Maintenance.....	14
TROUBLESHOOTING.....	15
Warnings.....	16
REGULATORY INFORMATION.....	17
Directive 1999/5/EC.....	17
SYMBOLS.....	18
GERÄTEBESCHREIBUNG.....	20
Bestimmungsgemäßer Gebrauch.....	20
Zubehör.....	20
AUFLADEN DES GERÄTS.....	21
EIN- UND AUSSCHALTEN DES GERÄTS.....	22

KOPPLUNG DES GERÄTS.....	23
BENUTZUNG DES GERÄTS.....	24
Lautstärkeanpassung.....	24
Visuelle Anzeigen.....	24
Aufbewahrung und allgemeine Handhabung.....	25
FEHLERSUCHE.....	26
Warnhinweise.....	27
INFORMATIONEN ZUR REGELKONFORMITÄT.....	28
Richtlinie 1999/5/EG.....	28
SYMBOLE.....	29
DESCRIPTION DE L'APPAREIL.....	31
Usage auquel l'appareil est destiné.....	31
Les accessoires.....	31
CHARGEMENT DE VOTRE APPAREIL.....	32
ALLUMER ET ÉTEINDRE L'APPAREIL.....	33
APPAIRAGE DE VOTRE APPAREIL.....	34
UTILISATION DE VOTRE APPAREIL.....	35
Le réglage du son.....	35
Indicateurs visuels.....	35
L'entretien.....	36
PROBLÈMES ET SOLUTIONS.....	37

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記号.....	86
DESCRIPTION OF DEVICE.....	88
Intended use.....	88
Accessories.....	89
CHARGING YOUR DEVICE.....	90
TURNING THE DEVICE ON AND OFF.....	91
PAIRING YOUR DEVICE.....	92
USING YOUR DEVICE.....	93
Sound adjustment.....	93
Visual indicators.....	93
Maintenance.....	94
TROUBLESHOOTING.....	95
Warnings.....	96
REGULATORY INFORMATION.....	97
Directive 1999/5/EC.....	111
FCC and IC statements.....	112
SYMBOLS.....	115

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Visual indicators

(See illustration 5)

Your device has two LEDs that flash with different colours. The purpose of these colours is to inform you of which mode your device is in at any given time.

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Maintenance

- Clean the device with a damp cloth. Do not use any liquids or chemicals.
- Do not immerse in liquids.
- When the device is not in use keep it in a dry location out of reach of children and pets.
- Never try to open or repair the device yourself. This may only be done by authorised technicians.

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SYMBOLS

Symbols commonly used by Widex A/S in medical device labelling
(labels/IFU/etc.)

Symbol	Title/Description
--------	-------------------

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WEEE mark

“Not for general waste”. When a product is to be discarded, it must be sent to a designated collection point for recycling and recovering to prevent the risk of harm to the environment or human health as a result of the presence of hazardous substances.

Symbol	Title/Description
	CE mark The product is in conformity with the requirements set out in European CE marking directives.
	RCM mark The product complies with electrical safety, EMC and radio spectrum regulatory requirements for products supplied to the Australian or New Zealand market.
	Interference Electromagnetic interference may occur in the vicinity of the product.

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1.6 Clip

Intended use

The COM-DEX Remote Mic is a microphone with an integrated radio intended for transmitting sound to a COM-DEX connected to the user's hearing aid(s). The aim is to improve the hearing aid user's ability to hear speech in difficult listening situations.

For optimal connectivity the device and the COM-DEX must be within sight of each other.

NOTE

Your Remote Mic is only intended to be used with a COM-DEX and it will not work with any other device.

Accessories

(See illustration 2)

USB cable

CHARGING YOUR DEVICE

(See illustration 3)

Connect the USB cable to your device and the power supply. It takes approximately 3 hours to charge your device.

NOTE

Before you can use the device, you must disconnect the USB cable.

Battery life: 8 hours streaming, 8 days standby.

TURNING THE DEVICE ON AND OFF

To turn on your device:

- Press the power button and keep it pressed for 3 seconds.

To turn off your device:

- Press the power button and keep it pressed for 3 seconds.

The device turns off automatically if it has not been streaming sound for more than 30 minutes.

PAIRING YOUR DEVICE

Pairing your device to your COM-DEX

To pair your device to your COM-DEX for the first time, follow these steps:

1. Make sure that your COM-DEX is in pairing mode. Consult the COM-DEX user guide.
2. Turn on your device. The first time you do this, the device will be set in pairing mode and it will automatically connect to your COM-DEX.

NOTE

Once you have paired your device to your COM-DEX, it will remain so and will connect automatically if your COM-DEX is on.

Pairing your device with another COM-DEX:

1. Turn on your device.
2. Reset your device to factory settings by pushing both buttons for 5 seconds. The device will be set to pairing mode.
3. Follow the pairing procedure described above.

NOTE

Your Remote Mic is only intended to be used with a COM-DEX and it will not work with any other device.

USING YOUR DEVICE

(See illustration 4)

Place the device on the shirt collar of the person you'd like to hear.

Turn on the device.

To start and stop streaming, simply push the streaming button on your Remote Mic. While streaming, the LED will show a green light.

Sound adjustment

The sound that is streamed from the device to COM-DEX can be adjusted in 3 levels. Simply push the power button on the side of the device and the sound level will be altered in a three step circle (highest level – lowest level – back to normal level).

Visual indicators

(See illustration 5)

Your device has two LEDs that flash with different colors. The purpose of these colors is to inform you of which mode your device is in at any given time.

- Constant red light on **both** LEDs: The device is charging.
- Constant green light on **both** LEDs: The device is fully charged.
- Flashing blue-red-blue-red light on **both** LEDs: The device is ready to be paired.
- Flashing blue light on the **top** LED: The device is on, but there is no connection to the paired COM-DEX.

- Flashing green light on the **top** LED: The device is on and connected to the COM-DEX, but it is not streaming.
- Constant green light on the **top** LED: The device is on and streaming.
- Occasionally flashing red light on **both** LEDs: The battery is low.

Maintenance

- Clean the device with a damp cloth. Do not use any liquids or chemicals.
- Do not immerse in liquids.
- When the device is not in use keep it in a dry location out of reach of children and pets.
- Never try to open or repair the device yourself. This may only be done by authorized technicians.

TROUBLESHOOTING

Problem	Possible cause	Solution
The device does not work	1. Remote Mic battery low	1. Charge Remote Mic battery
	2. COM-DEX battery low	2. Charge your COM-DEX
	3. Out of range	3. Move into range
No sound	1. Not paired to COM-DEX	1. Repeat the pairing process
	2. COM-DEX is out of the hearing aids' range	2. Move into range
	3. Hearing aid batteries are low	3. Change the hearing aids' batteries

If the problem persists, contact your hearing care professional.

Warnings

Do not use the device in mines or other areas with explosive gases.

Do not expose to extreme temperatures or high humidity.

When the device is connected via micro USB cable to external mains-operated equipment such as a charger, PC, audio source or similar, this equipment must comply with IEC 60065, IEC 60950-1, IEC 60601-1 or equivalent safety standards.

REGULATORY INFORMATION

The following Table summarizes the technical details of the Bluetooth and WidexLink technology as it is implemented in the COM-DEX Remote Mic.

Unit	COM-DEX Remote Mic	Hearing aid
Type	Bluetooth*	W-link
Antenna type	Embedded ceramic antenna	Inductive antenna
Antenna dimensions	NA	Ø1.8 mm L – 4.85mm
Modulation	FHSS/GFSK, $\pi/4$ DPSK, 8DPSK	FSK
Magnetic Field Strength	NA	-54dB μ A/m at 10m distance
Output power	2.5 mW EIRP +4dBm	29pW EIRP** -75 dBm
Range	<10m between COM-DEX Remote Mic and Bluetooth device	<1m remote unit to hearing aid <30cm between hearing aids
Center frequency	2.45GHz ISM	10.6MHz

Unit	COM-DEX Remote Mic	Hearing aid
Channels	80 channels	Single channel radio
Bandwidth	1MHz	660kHz (-15dB)
Duty Cycle ***	< 100 %	< 5 %
Data rate	3 Mbit/s	212kbit/s raw channel capacity
Data flow	Time division duplex (TDD)	Simplex or semiduplex
Protocol	Packet based protocol, time divided	Random access – no collision avoidance

* Bluetooth specification v 4.1 + EDR/LE published by the Bluetooth Special Interest Group (SIG)

** EIRP = Equivalent isotropically radiated power

*** Duty Cycle averaged over 1 hour of operation

Cables and transducers:

The COM-DEX Remote Mic is using a micro USB cable for charging.
Guidance and manufacturer's declaration

Electromagnetic emissions The COM-DEX Remote Mic is intended for use in the electromagnetic environment specified below. The customer or the user of a COM-DEX Remote Mic should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The COM-DEX Remote Mic uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The COM-DEX Remote Mic is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply
Harmonic emissions IEC 61000-3-2	Not applicable *)	

Emissions test	Compliance	Electromagnetic environment - guidance
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable *)	network that supplies buildings used for domestic purposes.

*) Battery powered equipment
Electromagnetic immunity

The COM-DEX Remote Mic is intended for use in the electromagnetic environment specified below. The customer or the user of a COM-DEX Remote Mic should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electro-static discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrical fast transients/burst IEC 61000-4-4	± 2 kV for power line supplies ± 1 kV for input/output lines	Not applicable *) Not applicable *)	Not applicable *)
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable *) Not applicable *)	Not applicable *)

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Not applicable *)	Not applicable *)

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

NOTE UT is the a.c. mains voltage prior to the application of the test level.

*) Battery powered equipment

Electromagnetic immunity – cont.

The COM-DEX Remote Mic is intended for use in the electromagnetic environment specified below. The customer or the user of a COM-DEX Remote Mic should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the COM-DEX Remote Mic, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Conducted RF			Recommended separation distance
IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	d = 1.2 ÖP
Radiated RF	3 V/m	3 V/m	d = 1.2 ÖP 80 MHz to 800 MHz
IEC 61000-4-3	80 MHz to 2.5 GHz		d = 2.3 ÖP 800 MHz to 2.5 GHz

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a, should be less than the

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
			compliance level in each frequency range b. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
---------------	-------------------------	------------------	--

-
- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the COM-DEX Remote Mic is used exceeds the applicable RF compliance level above, the COM-DEX Remote Mic should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or re-locating the COM-DEX Remote Mic.
- b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.
-

Recommended separation distances

Recommended separation distances between portable and mobile RF communication equipment and the COM-DEX Remote Mic.

The COM-DEX Remote Mic is intended for use in the electromagnetic environment in which RF disturbances are controlled. The customer or the user of the COM-DEX Remote Mic can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the COM-DEX Remote Mic as recommended below, according to the maximum output power of the communications equipment.

**Rated maximum
output power of
transmitter
W**

**Separation distance according
to frequency of transmitter
m**

	150 kHz to 80 MHz d = 1.2 ÖP	80 MHz to 800 MHz d = 1.2 ÖP	800 MHz to 2.5 GHz d = 2.3 ÖP
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

This COM-DEX Remote Mic may be interfered with by other equipment even if that other equipment complies with CISPR emission requirements.

EMC/EMI Compliance

The COM-DEX Remote Mic complies with the following EMC/EMI standards:

Standard	Test type	Note
47 CFR Part 15, subpart C	RF emissions	USA Federal Communications Commission (FCC) requirements to intentional radiators.
EN 300 328 V1.9.1	RF emissions incl. Spurious emission	EMC and radio spectrum matters for data transmission systems in the 2.4 GHz ISM band and using wide band modulation techniques
IEC 60601-1-2:2007	EMC emission Immunity, RF and ESD	Medical electrical equipment. General requirements for basic safety and essential performance. Electromagnetic compatibility.

Standard	Test type	Note
EN 301 489-17 V2.2.1	Immunity, RF and ESD	Standard for Broadband Data Transmission Systems in the frequency range 9 kHz – 40 GHz

Directive 1999/5/EC

Hereby, Widex A/S declares that this COM-DEX Remote Mic is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC.

A copy of the Declaration of Conformity can be found at:

www.widex.com/doc.

N26346

Information regarding disposal

Do not dispose of hearing aids, hearing aid accessories and batteries with ordinary household waste.

Hearing aids, batteries and hearing aid accessories should be disposed of at sites intended for waste electrical and electronic equipment, or given to your hearing care professional for safe disposal.

FCC and IC statements

FCC ID: TTY-CDRM

IC: 5676B-CDRM

Federal Communications Commission Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

NOTE:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

Changes or modifications to the equipment not expressly approved by Widex could void the user's authority to operate the equipment.

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Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada.

To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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SYMBOLS

Symbols commonly used by Widex A/S in medical device labelling (labels/IFU/ etc.)

Symbol	Title/Description
--------	-------------------

Manufacturer

The product is produced by the manufacturer whose name and address are stated next to the symbol. If appropriate, the date of manufacture may also be stated.

Catalog number

The product's catalog (item) number.

Consult instructions for use

The user instructions contain important cautionary information (warnings/precautions) and must be read before using the product.


Warning


Text marked with a warning symbol must be read before using the product.

Symbol	Title/Description
	WEEE mark “Not for general waste” When a product is to be discarded, it must be sent to a designated collection point for recycling and recovering to prevent the risk of harm to the environment or human health as a result of the presence of hazardous substances.
	CE mark The product is in conformity with the requirements set out in European CE marking directives.
	RCM mark The product complies with electrical safety, EMC and radio spectrum regulatory requirements for products supplied to the Australian or New Zealand market.
	Interference Electromagnetic interference may occur in the vicinity of the product.

4

5

WIDEX A/S

Nymoellevvej 6, DK-3540 Lynge, Denmark
www.widex.com

Manual no.:
 9 514 0380 080 #01
 Issue:
 2016-07