UPPER ARM

Series 800 BLOOD PRESSURE MONITOR

INSTRUCTION BOOKLET

Model #BP3MW1-4YCVS
CVS Item #800229

Please Note:

This medical instrument must be used according to instructions to ensure accurate readings

Important Product and Safety Information

Follow Instructions for Use. This document provides important product operation and safety information regarding this Blood Pressure Monitor. Please read this document thoroughly before using the device and keep for future reference.

This blood pressure monitor is an automatic digital blood pressure measuring device for use by adults, on the upper arm at home or in the doctor's office. It enables a very fast and reliable measurement of the systolic and diastolic blood pressure, as well as the pulse rate, by way of the oscillometric method during cuff inflation. This device detects the appearance of irregular heartbeats during measurement and provides a warning signal when the irregular heartbeat is detected. The device can be connected via Bluetooth to a smart mobile device running the App.

This blood pressure monitor IS intended to be used:

- For self-measurement/monitoring of blood pressure and pulse in adults; whereas the person being measured may be the user/operator of the device.
- Within a home healthcare environment.
- With a cuff located upon the user's upper arm; ½ inch above the elbow over the artery as indicated on the cuff.

This blood pressure monitor IS NOT intended for use with:

- Pregnancy
- Preeclampsia
- Children under 12
- Neonatal patients

WARNING – Potentially hazardous situation that if not avoided may result in serious injury or death.

- Self-measuring means monitoring, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you independently alter the dosages of any drugs prescribed by your physician.
- Consult your physician before using this device if any of the following or similar conditions are present: arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, age, pregnancy, preeclampsia, renal diseases. Motion during measurement, including trembling or shivering may affect the measurement.

- Though not for use with children under 12, ensure that any children around this device are supervised; some parts are small enough to be swallowed and any present tubes or cables may provide a risk of strangulation.
- Ensure the cuff tubing is not kinked during use as harmful injury may occur due to the effect of blood flow interference caused by high pressure in the cuff not releasing.

DO NOT:

- ✗ Use this device if you think it is damaged or anything appears unusual (e.g. sporadic operation, open packaging upon purchase).
- ✗ Use the displayed pulse for checking the frequency of heart pacemakers as this device is not suitable for this action.
- ✗ Open/modify this device; inaccuracy and/or harmful injury may result.
- ✗ Conduct many frequent measurements as this may result in harmful injury due to blood flow interference.
- ✗ Place the Cuff over a wound as this may cause further injury.
- ✗ Place and pressurize the Cuff over/near any present intravascular access or therapy, or arteriovenous shunt, as this may cause blood flow interference and result in harmful injury.
- ✗ Place and pressurize the Cuff over a limb near the side of a mastectomy as this may cause harmful injury.
- ✗ Use the Blood Pressure Monitor on a limb simultaneously with other medical equipment on the same limb due to possible interference with such medical equipment.
- ✗ Maintain pressure in the Cuff applied to the limb for a prolonged amount of time. Ensure that circulation in the limb is not impaired by checking circulation if prolonged/repeated exposure to pressure occurs.
- ✗ Service the device when being used or when power is supplied. When the device is not in use and power is removed, cuff assembly, batteries and ac adapter (if applicable) may be replaced by the user with Microlife supplied replacements. No other parts/components are accessible.

CAUTION – Potentially hazardous situation that if not avoided may result in minor/moderate injury, property damage, and/or damage to the device

- This device contains sensitive electronics components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g., mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring

accuracy. Move the device to another location if interference is determined.

- When not using the Blood Pressure Monitor for extended lengths of time, remove the batteries to avoid potential battery leakage and damage to the monitor.
- When replacing the batteries, ensure all the batteries are replaced simultaneously to avoid battery damage and potential damage to the monitor. Microlife USA does not recommend using rechargeable batteries.
- To avoid inaccurate measurements and to lessen any discomfort from Cuff pressure, ensure the Cuff is placed correctly on the limb and fits correctly when snug (not tight), as indicated by markings with the Cuff.
- Consult your physician in cases of frequent irregular heartbeat detections.
- This Risk Indicator feature is provided in order to help you understand your potential blood pressure risk. However, this feature is neither a diagnosis nor a substitute for a medical examination. It is important to consult with your physician to determine your risk.

DO NOT:

- ✗ Drop this device or expose it to strong vibrations; sensitive components may be affected resulting in inaccuracies and/or operational issues.
- ✗ Use the Blood Pressure Monitor outside of its specified operation temperature and humidity rating, or if stored outside of its specified storage temperature and humidity rating. Avoid storage in direct sunlight.
- ✗ Use this device in a moving vehicle; inaccurate measurements may result.
- ✗ Use third party accessories. Only use Microlife authorized accessories, such as cuffs or AC adapters, as those not approved for use with the device may provide inaccurate measurements, injury, and/or damage the device.

STANDARDS

In addition to the standards stated in the Instruction Manual:

- This medical device is compliant with medical device and non-invasive blood pressure monitor standards IEC 60601-1, IEC 60601-1-2, IEC 60601-1-11, and AAMI/ANSI/IEC 80601-2-30, and
- Electromagnetic standards IEC 60601—1-2 along with FCC Part 15, and
- Clinical Testing per standard ISO 81060-2:2013 was conducted on blood pressure device using the same measurement technology.

Please note: According to international standards, your monitor should be checked for accuracy every 2 years.

TYPE

IP20 Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water. Keep Dry. 

Type BF Applied Part

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Expected Life

Monitor 5 Years

Cuff 2 Years

FCC

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Changes or modifications to the product are not approved by Microlife USA and could void the user's authority to operate the equipment under FCC jurisdiction.

This equipment has been verified to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment on and off, the user is encouraged to try to correct the interference by increasing the distance between the product and the affected device; or if applicable, reorientation the receiving antenna, or moving the product's power plug to another receptacle.

Trademark Usage:

Apple, the Apple logo, iPad, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

Android and Google Play are both trademarks of Google Inc.

The Bluetooth® Smart word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks in this Blood Pressure Monitor is under license. Other trademarks and trade names are those of their respective owners.

Automatic Blood Pressure Monitor

Instruction Manual

Table of contents

1. Introduction

- 1.1. Your automatic blood pressure monitor
- 1.2. Gentle+™ Measurement System

2. How is blood pressure measured?

- 2.1. What is blood pressure, and why does it fluctuate?
- 2.2. Measurement guidelines
- 2.3. What can I do to change my blood pressure?

3. Components of your blood pressure monitor

4. How do I get started?

- 4.1. Inserting the batteries
- 4.2. Using the touch pad
- 4.3. Setting the date and time
- 4.4. Using the AC power adapter
- 4.5. Tube connection
- 4.6. Select the user
- 4.7. Select the measuring mode: standard or measurement averaging mode
- 4.8. Measurement averaging mode technology

5. How do I take a measurement?

- 5.1. Preparing to take a measurement
- 5.2. Common errors
- 5.3. Using the cuff
- 5.4. Taking a measurement
- 5.5. Memory: reviewing readings
- 5.6. Stopping a measurement
- 5.7. Setting the medication reminder
- 5.8. Hypertension risk indicator
- 5.9. Irregular heartbeat detector
- 5.10. Battery change indicator
- 5.11. AM/PM Averaging

6. Bluetooth functions & application (App)

6.1. Activating your Bluetooth®

6.2. Pairing your device

7. Error messages/troubleshooting

8. Care and maintenance

9. Limited lifetime warranty

10. Technical specifications

11. How to contact us

1. Introduction

1.1. Your automatic blood pressure monitor

Thank you for purchasing a fully automatic blood pressure monitor. Your monitor is designed to provide fast and reliable digital readings of your pulse, and systolic / diastolic blood pressure using the oscillometric method on your upper arm. It offers clinically proven accuracy and has been designed to be user friendly.

Before using your blood pressure monitor, please read this instruction manual carefully to ensure correct use. If you have additional questions regarding blood pressure measurements please contact your doctor.

1.2. Gentle+™ Measurement System

Gentle+™, our most advanced system to date, combines clinical accuracy and premium features including:

- A quieter, more comfortable measurement experience that reduces the amount of time your arm is under pressure. Less pressure = more comfort!
- An ergonomic cuff design that is easier to use and fits a wider range of arm sizes (for arms 8.7"-16.5").

2. How is blood pressure measured?

2.1. What is blood pressure, and why does it fluctuate?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is affected by fine muscles in the blood vessel walls.

Your level of arterial blood pressure changes periodically during heart activity. During the "blood ejection" (Systole), the value is highest (systolic blood pressure value). At the end of the heart's "rest period" (Diastole), pressure is lowest (diastolic blood pressure value).

Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2. Measurement guidelines

Blood pressure is very high if your systolic blood pressure is over 160 mmHg and/or your diastolic pressure is above 100 mmHg, while at rest. In this case, please consult your physician immediately. Long-term values at this level endanger your health due to continual damage to the blood vessels in your body.

If your systolic blood pressure values are between 140 mmHg and 159 mmHg and/or the diastolic blood pressure values are between 90 mmHg and 99 mmHg, consult your physician. Regular self-checks are necessary.

If you have blood pressure values that are too low (i.e., systolic values under 105 mmHg and/or diastolic values under 60 mmHg), consult your physician.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. **Never use the results of your measurements to independently alter the medication doses prescribed by your physician.**

Which values are normal?

The following standards for assessing high blood pressure (in adults) have been established by the National Institutes of Health JNC7, 2003.

Category	Systolic (mmHg)	Diastolic (mmHg)
Normal	<120	and <80
Pre-Hypertension	120-139	or 80-89
Hypertension		
Stage 1 Hypertension	140-159	or 90-99
Stage 2 Hypertension	≥160	or ≥100

Additional information

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called "labile hypertension." Consult your doctor.

2.3. What can I do to change my blood pressure?

- Consult your doctor.
- Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result from arteriosclerosis. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension and secondary hypertension. The latter group can be ascribed to specific organ malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- There are measures which you can take to reduce and even prevent high blood pressure. These measures must be permanent lifestyle changes.

1) Eating habits

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

2) Previous illnesses

- Consistently follow all medical instructions for treating illness such as:
 - Diabetes (diabetes mellitus)
 - Fat metabolism disorder
 - Gout

3) Habits




- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption (e.g., coffee).

4) Physical constitution

- **After a preliminary medical examination, do regular exercise.**
 - Choose sports which require stamina and avoid those which require strength
 - Avoid reaching the limit of your performance.
 - With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your exercise routine. You must receive advice regarding the type and extent of exercise that is appropriate for you.

3. Components of your blood pressure monitor


a) Measuring unit

b) Wide range soft cuff:

For arm circumference 22 cm - 42 cm (8.7" - 16.5")

If you ever need to buy a replacement cuff, call CVS/pharmacy® Blood Pressure Support toll-free at 1-866-464-6184.


**Please Note:**

Arm circumference should be measured with a measuring tape in the middle of the relaxed upper arm. Do not force cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port. If the cuff is too small, call 1-866-464-6184 for more information.

4. How do I get started?**4.1. Inserting the batteries**

After you have unpacked your device, insert the batteries. The battery compartment is located on the bottom of the device.

- a) Remove the battery cover.
- b) Insert the batteries (4 x AA 1.5 V), observing the indicated polarity.
- c) If a battery warning appears in the display, the batteries are discharged and must be replaced.

Low Battery Indicator

Attention!

- After the low battery indicator appears, the device won't function until the batteries have been replaced.
- Please use AA Long-Life or Alkaline 1.5 V batteries and replace them all at once.
- Rechargeable batteries are not recommended for use.
- If the blood pressure monitor is not used for long periods, remove the batteries from the device.

4.2. Using the touch pad

This monitor features three touch pad "buttons" which operate similarly to traditional buttons, but require only a light touch of the finger to operate.

- The buttons are MEMORY, START/STOP, and TIME.

4.3. Setting the date and time

1. After the new batteries are activated, the year number flashes in the display. You can advance the year by pressing the "M" button. To confirm and then set the month, press the TIME button.
2. You can now set the month using the "M" button. To confirm and then set the day, press the TIME button.
3. Please follow the instructions above to set the day, hour and minutes.
4. Once you have set the last minute and pressed the TIME button, the date and time are set and the time is displayed.
5. To change the date and time, press and hold TIME for approximately 3 seconds until the MAM symbol and ON/OFF flash onscreen. Release and press TIME again and the year number will flash. Enter the new values as described above.

Memory Button


Time/Date Button

4.4. Using the AC power adapter

You may also operate this monitor using the included AC adapter. Use only the included AC adapter to avoid damaging the unit.

- a) Ensure that the AC adapter and cable are not damaged.
- b) Plug the adapter cable into the AC adapter port on the right side of monitor.
- c) Plug the AC adapter into a 110 V power socket (U.S. or Canada).
- d) Test that power is available by pressing the START/STOP button.

Note:

- No power is taken from the batteries while the AC adapter is connected to the instrument.
- If the power is interrupted during a measurement (e.g., by removal of the adapter from the wall socket), the instrument must be reset by removing the plug from the instrument.

If you have any questions regarding the AC adapter, call CVS/pharmacy® Blood Pressure Support at 1-866-464-6184.

4.5. Tube connection

Insert the cuff tube into the opening on the left side of the monitor.

4.6. Select the user

This blood pressure monitor is designed to store 120 measurements for each of two users. In addition, there is a guest mode in which results are not stored.

Before taking a measurement, be certain that the correct user has been selected.

- With the unit off, press and release the TIME button to cycle through users ("1", "2" or guest mode (OFF)).
- Press the START/STOP button to make your selection.

4.7. Select the measuring mode: standard or measurement averaging mode (MAM)

This instrument enables you to select either standard (single measurement) or measurement averaging mode (automatic triple measurement).

- With the unit off, press the TIME button and hold for 3 seconds, the screen will show the current mode.
- Press the "M" button to switch between standard and averaging mode.
- To select Standard mode, select off.
- To select Averaging mode, select on.

If you select off, then only one measurement will be taken.
If you select on, the unit will inflate and deflate three times resulting in one final average.

4.8. Measurement averaging mode technology (MAM)

- In Measurement Averaging Mode (MAM), three measurements are automatically taken in succession and the result is then automatically analyzed and displayed. Because your blood pressure constantly fluctuates, a result determined in this way is more reliable than one produced by a single measurement.
- After pressing the START/STOP button the selected Measurement Averaging Mode appears in the display as the MAM symbol.
- The bottom, right-hand section of the display shows a 1, 2 or 3 to indicate which of the 3 measurements is currently being taken.

5. How do I take a measurement?

Please note: You should always be seated before and during measurement.

5.1. Preparing to take a measurement

- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors influence the measurement result. Find time to relax by sitting in an armchair that gives back and arm support in a quiet atmosphere for about 5 minutes before your measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left) and do not cross your legs during the measurement.
- Compare measurements at the same time of day, since blood pressure changes during the course of the day (as much as 20–40 mmHg).

5.2. Common errors

Note: Comparable blood pressure measurements always require the same conditions. Conditions should always be quiet.


- All efforts by the user to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not flex any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower or higher than the heart, an erroneously high or low blood pressure will be measured. Each 15 cm (6") difference in height between your heart and the cuff results in a measurement error of 10 mmHg.
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is extremely important. Cuff size is dependent upon the circumference of the arm (measured in the center). The permissible range is printed on the cuff. If this is not suitable for your use, please call 1-866-464-6184.
- A loose cuff or a sideways protruding air pocket causes false measurement values.
- With repeated measurements, blood accumulates in the arm, which can lead to false results. Consecutive blood pressure measurements should be repeated after a 1 minute pause in order to allow the accumulated blood to flow away.

5.3. Using the cuff


- a) Pass the end of the cuff through the flat metal ring so that a loop is formed. The hook and loop material must be facing outward. (Ignore this step if the cuff has already been prepared).

- b) Slide the cuff over the left upper arm so that the tube is closer to your lower arm.

- c) Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff, with the red edge, lies approximately $\frac{1}{2}$ " (1 to 2 cm) above the elbow and that the tube is closer to the inner side of the arm.

Important! The red strip on the edge of the cuff (Artery Mark) must lie over the artery which runs down the inner side of the arm.

TIP: Align red artery mark to pinky finger.

- d) Tighten the cuff by pulling the end and close the cuff by affixing the hook and loop material.

- e) There should be little free space between the arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Clothing must not restrict the arm. Any piece of clothing which does must be removed. Cuffs that don't fit properly result in false measurement values. Measure your arm circumference if you are not sure of proper fit.

- f) Lay your arm on a table so the cuff is at the same height as your heart. Make sure the tube is not kinked.


- g) Remain seated quietly for 5 minutes before you begin the measurement.

Comment

If it is not possible to fit the cuff to your left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.

Comparable blood pressure measurements always require the same conditions (relax for several minutes before a reading).

Cuff on
right arm

5.4. Taking a measurement

After the cuff has been appropriately positioned the measurement can begin. Remain still during your measurement, do not flex muscles.

- Remain seated quietly for 5 minutes before you begin the measurement.
- Press the POWER button to begin. After pressing the POWER button you will hear a series of short beeps and the pump begins to inflate the cuff. The measurement is performed while cuff inflates. The pumping speed can vary, which is completely normal. The increasing cuff pressure will show on the display.
- When the device has detected the pulse, the heart symbol in the display begins to blink, and a cuff signal is shown.
- When the measurement has been concluded, the air will automatically release from the cuff. The measured systolic and diastolic blood pressure values, as well as the pulse, are now displayed.
- The measurement results are displayed until you switch the device off. If no button is pressed for 1 minute, the device switches off automatically.
- When the unit is set to the MAM (Measurement Averaging Mode) setting, 3 separate measurements will take place in succession, after which your result is calculated and displayed as a single, averaged measurement. There is a 15 second resting time between each measurement. A count-down indicates the remaining time and a beep will sound for 5 seconds before the 2nd and 3rd readings begin.

Pumping
Pressure

162

Measuring

146
♥

Systolic

129

Diastolic

78
75

Pulse

Measurement
Complete

If one of the measurements causes an error message, it will be repeated one more time. If any additional error occurs, the measurement will be discontinued and error code displayed.

5.5. Memory: reviewing readings

At the end of a measurement, this monitor automatically stores each result with date and time. This unit stores 120 memories for each of 2 users.

Viewing the stored values

With the unit off, press the "M" button. The display first shows "A," then shows an average of all measurements stored in the unit. Please note: Measurements for each user are averaged and stored separately. Be certain that you are viewing the measurements for the correct user.

Pressing the "M" button again displays the previous value. To view a particular stored memory, press and hold the "M" button to scroll to that stored reading.


Memory full

When the memory has stored 120 results, a new, measured value is stored by overwriting the oldest value.

Clear all values

If you are sure that you want to permanently remove all stored values, hold down the "M" button (the instrument must have been switched off beforehand) until "CL" appears and then release the button. If you do not want to clear the values, press the START/STOP button. To permanently clear the memory, press the "M" button while "CL" is flashing.

Individual values cannot be cleared.


5.6. Stopping a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g., the patient feels unwell), the START/STOP button can be pressed at any time. The device then immediately lowers the cuff pressure automatically and enters sleep mode.

5.7. Setting the medication reminder

This instrument allows you to set two alarm times at which an alarm signal will then be triggered. This can be a useful aid, for instance as a reminder to take medication or to remind you to take your blood pressure at the same time each day.

1. To set an alarm time, press the TIME button (the instrument must have been switched off beforehand), and immediately afterwards the "M" button, and hold both down until the bell symbol appears in the bottom left of the display. Then release both buttons. The flashing "1" in the display indicates that the first alarm time can now be set.
2. Press the TIME button to set the hours – the hours display flashes and pressing the "M" button allows you to set the alarm hour. To confirm, press the TIME button.
3. The minute display will now flash. The minutes can be set using the "M" button. To confirm, press the TIME button again.
4. The bell symbol will now flash. Use the "M" button to select whether the alarm time is to be active (bell) or inactive (crossed-out bell). To confirm, press the TIME button.

TIME PM 3:16

TIME PM 3:16

To set a second alarm time, proceed as above but if the "1" flashes, press the "M" button to select "2" and confirm with the TIME button.


- An active alarm time is indicated by the bell symbol in the display.
- The alarm sounds at the set time every day.
- To switch off the alarm when it is sounding, press the START/STOP button.
- To permanently switch off the alarm, proceed as above (steps 1-4) and select the crossed-out bell symbol. This will then disappear from the display.
- The alarm times must be re-entered each time the batteries are replaced.

5.8. Hypertension risk indicator

The bars on the left-hand edge of the display show you the range within which the indicated blood pressure value lies. Depending on the height of the bar, the readout value is either within the normal (green), borderline (yellow) or danger (orange, red) range.

The classification is based on standards established by the National Institutes of Health JNC7, 2003.

Refer to the chart in section 2.2 of this instruction manual for details of the classifications.

The traffic light bar illuminates according to your measurement.

- If the bar is in the green zone, your measurement is "Normal," according to NIH standards.
- If the bar is in the yellow zone, your measurement is "Pre-Hypertension."
- If the bar is in the orange zone, it is "Stage 1 Hypertension."
- If the bar is in the red zone, it is "Stage 2 Hypertension."

5.9. Irregular heartbeat detector

The appearance of this symbol  indicates that certain pulse irregularities were detected during the measurement. In this case, the result may deviate from your normal basal blood pressure – repeat the measurement. In most cases, this is no cause for concern. However, if the symbol appears on a regular basis (e.g., several times a week with measurements taken daily), we advise you to tell your doctor.

Please show your doctor the following explanation:

Information on frequent appearance of the irregular heartbeat symbol

This instrument is an oscillometric blood pressure monitor device that also analyzes pulse frequency during measurement. The instrument is clinically tested.


If pulse irregularities occur during the measurement, the irregular heartbeat symbol is displayed with the measurement.

If the symbol appears frequently or if it suddenly appears more often than usual, we recommend the patient seek medical advice. The instrument does not replace a cardiac examination but serves to detect pulse irregularities at an early stage.

5.10. Battery charge indicator

Batteries almost discharged

When the batteries are approximately 75% used, the battery symbol will flash a few times as soon as the instrument is switched on. Although the instrument will continue to measure reliably, you should obtain replacement batteries.

Batteries discharged – replacements required

When the batteries are discharged, the battery symbol will appear, unblinking, as soon as the instrument is switched on. You cannot take any more measurements and must replace the batteries.

1. Open the battery compartment on the bottom of the instrument.
2. Replace all of the batteries at the same time – ensure they are correctly connected, as shown on the symbols in the compartment.

3. The memory retains all values although date and time (and possibly also set alarm times) must be reset – the year number will flash automatically after the batteries are replaced.
4. To set date and time, follow the procedure described in Section 4.3.

Note:

Use four new, Long-Life 1.5 V AA batteries. Do not use batteries beyond their expiration date. If the monitor is not going to be used for a prolonged period, the batteries should be removed.

Rechargeable batteries are not recommended for use with this device.

5.11. AM/PM Averaging

Doctors recommend that you compare blood pressure measurements taken at the same time each day to help ensure consistent conditions during measurements. To assist you, this blood pressure monitor is equipped with AM/PM averaging capabilities.

To view AM/PM readings:

- With the monitor off, press “M” (Memory). The average reading for the current user will appear.
- Press “M” again. The average of all AM readings (4 AM-12 PM), designated by the sun icon, will be displayed on the screen.
- Press “M” again. The average of all PM readings (6 PM-12 AM), designated by the moon icon, will be displayed on the screen.
- After 30 seconds, the display will show the last reading taken.
- To scroll through all individual readings, press “M” repeatedly.

6. Bluetooth functions and application (App)

Please download the Microlife Connected Health App

(referred to as App in this document) from Apple's App Store® or Google Play® before connecting your devices.

Questions?


- To access the tutorial, go to the App home screen (Microlife Connected Health) and press the "Settings & Help" button.
- Press "Tutorial" and select your topic.
- If you have any issues with your setup, please contact customer service toll free: 1-866-464-6184 for additional support

6.1. Activating your Bluetooth®

The bluetooth® icon indicator on your blood pressure monitor, located in the upper left portion of the screen, is designed to provide information about the connection between your device and the monitor. Please note the following:

- Locate the bluetooth® switch on the left side of the monitor, and slide it to the "ON" position
- The bluetooth® icon on your blood pressure monitor will flash initially when ready to connect with a device, and continue flashing for 2 minutes until a connection occurs.
- If a connection does not occur within 2 minutes, the Bluetooth® function on the monitor will automatically turn off.
- To manually activate the Bluetooth® on your monitor, press and hold the power button for approximately 8 seconds or until the Bluetooth® icon starts blinking.
- Note: At least one blood pressure reading must be taken and recorded on the monitor before the data can be transferred


Bluetooth switch

Bluetooth® icon

6.2. Pairing your device:

- Download the App to your phone.
- Create User Profile on App.
- Confirm phone's Bluetooth® is turned "ON".
- Open the App on your phone.
- Select "Dashboard" from the App's home screen.
- Confirm your monitor's Bluetooth® is turned "ON"
- Confirm at least one blood pressure reading has been taken and recorded on the monitor. This needs to occur before the data can be transferred.
- Press and hold the POWER/START button on the monitor for approx. 8 seconds or until the Bluetooth® icon starts blinking.
- On the App press the "Sync" button located on the "Dashboard" screen.
- Data in monitor should successfully transfer over to the App on your phone.

Bluetooth® Switch

Setting up your blood pressure monitor:

Confirm the correct user is selected on your blood pressure monitor:

- Confirm the user icon on the top left hand side of the monitor screen is set to user 1 or 2.
- To change users on your monitor, see section 4.6.

7. Error messages/troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (example: Error no. 2).

**Err
2**

Error No.	Possible cause(s)/Solutions
ERR 1	The tube may have loosened, or no pulse was detected.* Ensure cuff connections are tight with proper cuff placement. See section 5.3.
ERR 2	Unnatural pressure impulses influenced the measurement result. Reason: The arm was moved during the measurement (artefact). Repeat measurement, keeping still and quiet.
ERR 3	Inflation of the cuff takes too long. The cuff is not correctly seated or the hose connection is not tight. Re-position cuff and repeat the measurement.
ERR 5	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.
HI	The cuff pressure is too high. Relax for 5 minutes and repeat the measurement.*
LO	The pulse is too low (less than 40). Repeat the measurement.*

***If this or any other problem occurs repeatedly, please consult your doctor.**

Other possible errors and their solutions

If problems occur when using the device, the following points should be checked and, if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains blank when the instrument is switched ON although the batteries are in place.	• Check battery installation. • If the display is unusual, remove the batteries and then exchange them for new ones • Check polarity
The pressure does not rise although the pump is running.	• Check the connection of the cuff tube and connect properly.
The device frequently fails to measure the blood pressure values, or the values measured are too low or too high.	1. Fit the cuff correctly on the arm. (see 5.3) 2. Before starting the measurement make sure that the cuff is not too tight and that clothing is not exerting pressure on the arm. Take off articles of clothing if necessary. 3. Measure blood pressure again in complete peace and quiet.
Every measurement results in different values, although the device functions normally and normal values displayed.	• Please read the following information and the points listed under "Common sources of error".(5.2) Repeat the measurement. Blood Pressure changes constantly. The observed readings may accurately reflect your pressure.
The blood pressure values differ from those measured by my doctor.	• Record daily values, consult your doctor. • Pressure readings in your doctor's office may be higher due to anxiety. • Your blood pressure fluctuates throughout the day so record the time along with your reading.
After the instrument has inflated the cuff the pressure falls very slowly or not at all. (No reasonable measurement possible).	• Check cuff connections. Ensure the unit has not been tampered with.

8. Care and maintenance

- a) Do not store the unit under direct sunlight, at a high temperature, or in high humidity or dust. Performance may be degraded.
- b) The cuff contains a sensitive air-tight bubble. Handle this cuff carefully and avoid all types of stress through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and mild liquid dish soap. The cuff with bladder must not be washed in a dishwasher, clothes washer, or submerged in water.
- d) Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- e) Never open the monitor. This voids the manufacturer's warranty.

9. Limited lifetime warranty

Your Automatic Blood Pressure Monitor is warrantied for the estimated life of the Monitor by Microlife USA Inc, against manufacturer defects for the original purchaser only, from date of purchase.

The lifetime warranty applies to the monitor only. The following accessories are warranted for 1 year: cuff and adapter. Batteries are not covered by this warranty. There are no user serviceable parts inside.

The warranty does not apply to consequential and incidental damages, or damage caused by batteries, improper handling, and accidents. Professional use, not following the operating instructions, and alterations made to the monitor or accessory by third parties, are also not included in this warranty. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

Microlife USA Inc will investigate your concern. A monitor or accessory as defined by this warranty, determined to be out of specification, will be replaced and shipped to you at no cost. A monitor or accessory as defined by this warranty, determined to be within specification, will be returned to you with a report of findings, at no cost.

Please use the customer service contact information to reach Microlife USA Inc. regarding any warranty concerns. We ask that you please contact us before sending any product back in order to better identify, and more quickly process, your concern.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

10. Technical specifications:

Weight:	.468 g (with batteries)
Size:	.110 (W) x 120 (L) x 85 (H) mm
Storage temperature:	..-20 to +55°C (-4° to +131°F)
Humidity:	.15 to 90% relative humidity maximum
Operation temperature:	.10 to 40°C (50° to 104°F)
Display:	.LCD (Liquid Crystal Display)
Measuring method:	.Oscillometric
Pressure sensor:	.Capacitive
Measuring range:	
SYS	.60 to 255 mmHg
DIA:	.40 to 200 mmHg
Pulse:	.40 to 199 per minute
Cuff pressure display range:	.0-299 mmHg
Memory:	.Automatically stores the last 120 measurements for 2 users (total 240)
Measuring resolution:	.1 mmHg
Accuracy:	.Pressure within ± 3 mmHg or 2% of reading >200 mmHg .Pulse $\pm 5\%$ of the reading
Power source:	a) 4 AA batteries, 1.5 V b) AC adapter 6 V DC 600 mA (voltage 4.5 V DC to 6 V DC)
Battery life:	.Approximately 500 measurements
IP Classification:	.IP20
Accessories:	.Cuff type: Softfit™ cuff for arm circumference 22-42 cm (8.7"-16.5") Storage case

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© 2021 CVS/pharmacy
Distributed by: CVS Pharmacy, Inc.