



Vscan Air Dual Charger

User Manual



Model
5978755

GE HealthCare is a trademark of GE HealthCare Technologies Inc.
Manufactured under trademark license
Vscan Air is a trademark of GE HealthCare



Overview

The Vscan Air Dual Charger is specifically designed to charge the Vscan Air probe.

Signal Indicator

		Status
	BLUE — 5 sec	Charger initially plugged in - no probe on charger.
	BLINKING GREEN — — —	Probe battery charging.
	ORANGE — — —	Charging failure. Probe battery not full, charging paused.
100%	STATIC BLUE — — —	Probe battery fully charged.

Cautions

1. Do not place any foreign materials, such as metal objects, magnets, and magnetic stripe cards between the Vscan Air probe and the charger.
2. Use this charger in an environment temperature between 3°C ~ 40°C.
3. Do not disassemble the charger.
4. Use original or certified cable and charger for safe charging.
5. Do not defeat the safety purpose of the polarized or grounding-type plug.

Product Information

Product name: Vscan Air Dual Charger

Model: 5978755

Authorized by: GE HealthCare

Address:

BizLink (Kunshan) Co., Ltd.

No. 168, Nanhe Rd., Kunshan Economic & Technology Development Zone, Kunshan City, Jiangsu 215300, China

Manufacturer: BizLink (Kunshan) Co., Ltd.

Specifications

Input	100-240V~, 50/60Hz, 1.0A Max
Output	5Vdc, 5W Max x2
Size	Width: 190±3mm Depth: 155±3mm Height: 39.6±2mm
Weight	≤500g

Product Labels

Please refer to Appendix A. The label is located at the bottom side.

Service Life

5 years

List of Suggested Disinfecting Wipes

- | | |
|---|---|
| 1. Wet Wipe | 11. Mikrobac Tissues |
| 2. Soap Water | 12. Mikrozid Sensitive Wipes |
| 3. ALCOHOL (ISOPROPANOL) 70% | 13. Oxivir TB Wipes |
| 4. CaviWipes | 14. Sani-Cloth AF3 Germicidal Disposable Wipe |
| 5. Clinell Universal wipes | 15. Sani-Cloth Bleach Germicidal Disposable Cloth |
| 6. Clorox Healthcare Bleach Germicidal Wipes | 16. Sani-Cloth Prime Germicidal Disposable Cloth |
| 7. Sani-Cloth Plus Germicidal Disposable Wipe | 17. Protex |
| 8. SONO ULTRASOUND WIPES | 18. Easy Screen Cleaning Wipes |
| 9. Super Sani-Cloth Germicidal Disposable Wipes | 19. Tristel Trio |
| 10. Cleanisept Wipes | |

Note

Always consult the cleaner or disinfectant manufacturer's instructions for proper use of their product. Wear appropriate PPE as indicated by the manufacturer. Always follow the manufacturer's instructions concerning the use of the disinfectant and follow the contact time to be sure the disinfectant's kill claims are accomplished.

Compliance with 2014/53/EU Radio Equipment Directive

The following table provides information on the frequency bands used and the maximum RF transmit power of the product for sale in the EU and UK:

Operating Frequency: 127.7 kHz

Output Power (Field Strength): -2.9 dBμA/m at 10 meters

Service Contact

GEHC - Americas	GEHC - Europe	GEHC - Asia
Online Services Ultrasound Phone: 1-800-321-7937 (USA/Canada) 1-262-524-5300 (Latin America)	EMEA Online Center Phone: +49 (0)212-2802-652 (-OLC) (English/German all segments incl. training): Fax: +49(0)212-2802-431 EAGM Online Center Egypt Service Center: 00202 2322 1252 UAE Service Center: 00971 8003646	Online Services Ultrasound Asia China +(86) 800-810-8188 India +(91) 1800-425-8025 +(91) 1800-102-7750 +(91) 1800-425-7255 Japan +(81) 42-648-2940 Korea +(82) 2-1544-6119

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Caution:

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This device and its antenna(s) must not be co-located or operating in conjunction with any other antenna or transmitter.

FCC Radiation Exposure Statement

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 20cm between the radiator & your body.

For Canada

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

L'émetteur/récepteur exempt de licence contenu dans le présent appareil est conforme aux CNR d'Innovation, Sciences et Développement économique Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes:

1. L'appareil ne doit pas produire de brouillage.
2. L'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

IC Radiation Exposure Statement

This equipment complies with IC RSS-102 radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 20cm between the radiator & your body. Cet équipement est conforme aux limites d'exposition aux rayonnements IC établies pour un environnement non contrôlé. Cet équipement doit être installé et utilisé avec un minimum de 20cm de distance entre la source de rayonnement et votre corps.

EU/UK Declaration of Conformity

For the following equipment:

Vscan Air Dual Charger

(Product Name)

5978755

(Model Designation)

is herewith confirmed to comply with the requirements set out in the Council Directive on the Approximation of the Laws of the Member States relating to Electromagnetic Compatibility of Radio Equipment Directive 2014/53/EU and UK Radio Equipment Regulations 2017 SI 2017/1206. For the evaluation regarding the Directives, the following standards were applied:

EN 300 330 V2.1.1 (2017-02)

EN 50665:2017

EN 62368-1:2014+ A11:2017

EN 303 417 V1.1.1 (2017-09)

EN 62311:2008

BS EN 62368-1:2014+

EN 301 489-1 V2.2.3 (2019-11)

A11:2017

EN 301 489-3 V2.3.2 (2023-01)

RED Notified Body:

The notified body, Kiwa Netherlands (0063), performed an assessment of Article 3.1a, Article 3.1b and Article 3.2 and issued the EU-Type Examination certificate

UK Approval Body:

The Notified Body, Kiwa Netherlands (0558), performed the conformity assessment according to Schedule 6.1a, 6.1b, and 6.2 of the Radio Equipment Regulations 2017/UK and issued the UK-type examination certificate 3

Object of the declaration (further identification of the radio equipment allowing traceability; it may include a color image for the identification of the radio equipment



The following importer/manufacturer is responsible for this declaration:

BizLink (Kunshan) Co., Ltd.

(Company Name, Manufacturer)

No.168, Nanhe Road Economic & Technological Development Zone, Kunshan City, Jiangsu Province 215300, CHINA

(Company Address, Manufacturer)

Person responsible for this declaration:

Amiee Zhu

(Name, Surname, Manufacturer)

Engineer

(Position/Title)



(Legal Signature)

China

(Place)

2025/11/7

(Date)