SKG Health Technologies Co., Ltd.

Address: 23A Floor, Building 3, Zhongke R&D Park, No. 009, Gaoxin South 1st Road, High-tech Zone Community, Yuehai Street, Nanshan District, Shenzhen, Guangdong, China.

Website : https://global.skg.com/

If you want to know more about product specification informaation, please visit: https://www.skg.com/

Made in China

SKG


Eye Massager
KS700-E


USER MANUAL

英语


Enjoy the power-on experience

EN





Diagonally above



*The illustrations of the products and accessories are for illustration purposes only. The appearance shall be subject to the actual product in the packaging.

* Before using this product, please adjust the elastic band to a suitable length, then wear it around the eyes, adjust the position and tightness of the massage head according to the wearing situation, and then turn on the main unit to experience the massage effect. (See the diagram on the right for how to wear it)

Operation steps

1. Press and hold the main button for 1.5 seconds to turn on the device (accompanied by voice prompts and the indicator light turns on).
2. Short press the main button to switch modes; (To ensure the cooling effect, please avoid frequent mode switching or power on/off)
3. Double-click the main button to switch between mute and voice mode; (press twice quickly in succession)
4. Press and hold the main button for 1.5 seconds to shut down the device. A voice prompt will appear and the indicator light will go out (no voice prompt in silent mode).
   Note: After shutting down, in order to extend the service life of the product, the built-in fan of the product continues to run for about 2 minutes. This is a normal phenomenon and not a product abnormality.
5. Can connect to Bluetooth and play online music.
   Turn on the Bluetooth of the mobile phone system, search for "SKG_K7_XXXX", and after pairing and saving, it will automatically connect every time you turn on Bluetooth. When the mobile phone Bluetooth connects to the device and plays mobile music, the music on the host will stop and play mobile music.

Before use:

1. Please remove contact lenses, glasses, false eyelashes, hairpins and hair accessories first to avoid discomfort during massage.
2. For sensitive skin and facial allergic dermatitis, please use this product under the guidance of a doctor.
3. Do not use if you have symptoms of eye diseases such as glaucoma, cataracts, retinal detachment, or conjunctivitis.
4. After turning on the device, the cumulative time from using the cold compress mode to automatic shutdown is 6 minutes, and the cumulative time for other modes is 12 minutes. If the mode is switched during the use of the cold compress mode, the cumulative time will become 12 minutes before automatic shutdown.

After use:

1. After use before going to bed, be sure to remove it and avoid using it while sleeping.
2. The massage area may feel red and hot while feeling comfortable, which is normal.
3. After use, slight marks may be left on the face, which will return to their original state after a while.
4. Due to long-term eye-closing massage, your eyes may become blurred for a short time. Close your eyes and rest for 1 to 3 minutes, and your eyes will gradually become clear.

Product basic parameters

EN

Product Name: Eye Massager	**Model:** KS700-E
Rated charging voltage: 5V⎓	**Rated input power:** 10W
Charging port: Type-C	**Battery capacity:** 2000mAh/7.4Wh

Battery indicator light description

EN

product status	light status	Battery description
while working	Blue light always on	Electricity
	Orange light always on	Weak current (remaining time must be >5min)
	Orange light flashes 5 times, power off	Low battery (about to shut down)
	blue breathing light	Product cold compress status
	orange breathing light	Product hot compress status
While charging	Orange light flashes slowly	Charging
	Blue light always on	The battery level is higher than 98%

Daily maintenance of the product

EN

- When the surface of the massager is dirty, wipe it with a clean, soft cloth dampened with a little water. Do not wipe with volatile liquids like alcohol. The product cannot be disassembled or washed in water.
- Please cover the product when not in use for an extended period of time to prevent dust accumulation.
- Prevent water or other liquids from entering the machine during cleaning to avoid short circuits during use.

Troubleshooting Instructions

EN

If the product does not work properly, please check if the following situations have occurred. If our suggestions do not solve the problem, please return the product to the dealer for repair.

Condition	Possible Cause	Solution
Unable to Power On	Internal Battery Low	Charge the Product
	Product is Charging	Use After Fully Charged
Unable to Charge /No Charging Indicator	Charging Cable Not Properly Connected to Main Unit and Power Adapter	Reconnect Charging After Unplugging

⚠ Tips

- Please use under normal operating conditions: Temperature: 0℃ to 40℃; Relative Humidity: <80% RH, no condensation; Atmospheric Pressure: 50Kpa to 106Kpa.
- Please select normal transport and storage conditions: Environmental Temperature Range: -20℃ to +55℃; Relative Humidity Range: ≤95% RH, no condensation; Atmospheric Pressure: 50Kpa to 106Kpa.
- It is recommended to use the product no more than twice per day for no more than a half-hour at a time, based on the user's specific requirements.

Product Disposal Precautions

EN

Please support the recycling of rechargeable batteries (lithium-ion batteries)

⚠ **Danger**	This product's battery is a dedicated rechargeable battery that should not be used in other products. Do not charge the battery once it has been removed. • Do not throw into the fire or expose to excessive heat. • Do not puncture with nails, subject to heavy impact, or disassemble • Do not directly connect the battery's positive and negative terminals with metal • Do not store it with necklaces, hairpins, etc. • Do not charge, use, or store in high-temperature locations such as near a stove or under direct sunlight • All of the abovementioned scenarios can result in overheating, fire, or bursting.
⚠ **Warning**	Place the removed rechargeable battery out of the reach of infants to prevent accidental ingestion • If swallowed, seek medical attention immediately. • If battery leakage occurs, do not touch the liquid with bare hands. • Follow the safety instructions below: • There is the risk of blindness if battery fluid contacts the eyes. Do not rub the eye with your hands. Rinse with clean water immediately and seek medical attention. If battery liquid gets on the body or clothing, it may cause skin inflammation or injury. Rinse with clean water, and if discomfort persists, seek medical attention mediately. • Disassembly is only permitted when the product is being discarded. Failure to comply may result in fire, electric shock, injury, or other accidents. Please contact the retailer where you purchased the product for repairs (except battery replacement).

The lithium-ion battery used in this product is a valuable recyclable resource. Please remove the lithium-ion battery prior to disposal and hand over it to a professional organization capable of recycling and processing rechargeable batteries.

- The battery must be removed from the device before discarding the appliance.
- The device must be turned off before removing the battery.
- The battery should be disposed of safely.

Removal Steps:

1. Pry open the main unit's inner and outer shells using professional tools.
2. Unscrew the fixed screws of the internal bracket using a screwdriver.
3. Unscrew the screws at the circuit board.
4. Cut the battery leads using scissors.
5. Remove the battery.

Precautions

EN

⚠ Note

- Please read and comply with the guidelines in the "User Manual" before using this product. The company is not responsible for consequences resulting from improper usage or non-compliance with the product instructions.
- Technical specifications are subject to change without notice as we continually improve our products.
- Do not use this product after applying makeup to your face since it will destroy your makeup and may contaminate the skin's surface. If skin contamination occurs, please wipe with clean water.
- Do not use an eye massager while driving.
- Do not use an eye massager in excessively hot or humid areas.
  This product should not be used by individuals with stroke, high or low blood pressure.
- Do not use this product if you have glaucoma, cataracts, retinal detachment, conjunctivitis, or other eye disorders.
- Please avoid using an eye massager after consuming alcohol.
- Patients who have undergone foot surgery or have heart diseases are prohibited from using this product.
- Before using, ensure that there are no contact lenses in the eyes and that any hair clips or jewelry are removed.
- This is a personal care massage product and cannot be used to substitute medications or medical devices.
- Children, pregnant women, and women on their menstrual period are prohibited from using this product.
- Contains lithium battery. Do not place in a high-temperature environment.
- Use this product under the supervision of a doctor if you have sensitive skin, facial allergic dermatitis, or face areas with cosmetic materials.
- This product has a heating function; please be cautious if you insensitive to heat.
- If you feel discomfort during use, please stop immediately.
- This product should not be used to massage parts of the body other than the eyes.
- Do not disassemble this product yourself.
- This product is not waterproof; please do not wash with water or use in high-humidity environments. If the lining becomes dirty, you can gently wipe it with a moist cloth.
- The appliance it must only be supplied at safety extra low voltage corresponding to the marking on the appliance.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance.
- Cleaning and user maintenance shall not be made by children without supervision.
- If you experience skin irritation when using the body massager, please refrain from using it, and wait two to three days to see if your symptoms ease. If symptoms persist or worsen, please consult a physician.
- If you have eczema, allergies or asthma, you may be more likely to experience a skin irritation or allergy from using the massager.

This device complies with part 15 of the FCC Rules. Operation is subject to the condition that this device does not cause harmful interference (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
Any changes or modifications not expressly approved by the party responsible for compliance
could void the user's authority to operate the equipment.
NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.
If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
-- Reorient or relocate the receiving antenna.
-- Increase the separation between the equipment and receiver.
-- Connect the equipment into an outlet on a circuit different
from that to which the receiver is connected.
-- Consult the dealer or an experienced radio/TV technician for help.
The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

This device was tested for typical body-worn operations. To comply with RF exposure requirements, a minimum separation distance of 0cm must be maintained between the user's body and the handset, including the antenna. Third-party belt-clips, holsters, and similar accessories used by this device should not contain any metallic components. Body-worn accessories that do not meet these requirements may not comply with RF exposure requirements and should be avoided. Use only the supplied or an approved antenna. This device in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU. All essential radio test suites have been carried out.
1. CAUTION : RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE. DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS
2. The device complies with RF specifications when the device used at 0mm form your body EU Compliance Statement:
SKG Health technology Co.,Ltd. Here by declares that this device is in compliance with essential requirements and other relevant provisions of the Directive 2014/53/EU. A copy of the EU Declaration of Conformity is available online at https://www.skg.com/pages/skg-declaration-of-conformity.

Old electrical appliances must not be disposed of together with the residual waste, but have to be disposed of separately. The disposal at the communa collecting point via private persons is for free. The owner of old appliances is responsible to bring the appliances to these collecting points or to similar collection points. With this little personal effort, you contribute to recycle valuable raw materials and the treatment of toxic substances.

FR
Cet appareil et sa batterie se recyclent
À DÉPOSER EN MAGASIN
OU
À DÉPOSER EN DÉCHÈTERIE
Points de collecte sur www.quefairedemesdechets.fr
Privilégiez la réparation ou le don de votre appareil !

FR
DÉPLIANT
+BARQUETTE
+ÉTUI
BAC DE TRI

EU REP
OST_EU_20230214000009

SUCCESS COURIER SL
CALLE RIO TORMES NUM.1, PLANTA 1, DERECHA, OFICINA 3, Fuenlabrada, Madrid, 28947 Spain
successservice2@hotmail.com