

FCC ID Label Location

The product does not have enough space to place the label, and FCC ID is reflected in the usermanual and package.

User manual as below:

- The plug of the adapter is considered a disconnect device. The socket-outlet shall be installed near the equipment and shall be easily accessible.
- Rechargeable batteries are only to be charged by adults.
- The power delivered by the charger must be between min [1.5] Watts and max [5] Watts required by the radio equipment in order to achieve the maximum.

Eliminating the Risk, keep away from the flammable materials, heat sources and ignition sources while a battery is charged.
Charging should only take place in a well ventilated area and when someone is about to be aware and react to a potential issue. To reduce risk of fire, a battery should not be charged for more than 4 hours and disconnect charger from line voltage (AC outlet) by removing main plug even if the battery is not completely charged. Don't leave a device unattended while charging. Disconnect the charger when the device is fully charged: do not overcharge the device.



WARNING: CHOKING HAZARD—Small Parts.
Not suitable for children under 3 years old.

Authorized Accessories List

- Contact Reteless for assistance regarding repairs and service.
- For a list of Reteless-approved accessories for your radio model, visit the website <http://www.reteless.com>

FCC ID: 2A3NOTE103

31

32

Package as below:

RETEKESS

US **Hearing Aids**

Model: TE103 **Input:** 5V $\pm$ 1A **RF Power:** <10mW
Frequency: 2400MHz~2483.5MHz
Battery:
Earphone-Button Lithium-Ion Cell: LIR1040, 3.6V, 35mAh, Weight 0.98g
Base-Li-Ion Polymer Battery Pcd07: 3.7V, 400mAh/1.48Wh, Weight 4.5g

FCC ID: 2A3NOTE103
SRN: 23071211031100001
FDA Registration Type: TE103

WARNING: If you are younger than 18, do not use this.
You should go to a doctor, preferably an ear-nose-throat doctor (an ENT), because your condition needs specialized care. Over-the-counter hearing aids are only for users who are age 18 or older.

This information and other labeling, including the user instructional brochure, are available on the internet at: www.reteless.com

You may also write to support@reteless.com to request a copy of this information and other labeling.

Manufacturer's return policy:
To return the item by email to support@reteless.com. We'll work to ensure that both your return and credit refund are processed accurately.

This hearing aid is for adults with signs of mild to moderate hearing loss.
How do you know if you have this?

- You have trouble hearing speech in noisy places
- You find it hard to follow speech in groups
- You have trouble hearing on the phone
- Listening makes you tired
- You need to turn up the volume on the TV or radio, and other people complain it's too loud

Some people with hearing loss may need help from a hearing healthcare professional. How do you know if you need to see one?

- You can't hear speech even if the room is quiet.
- You don't hear loud sounds well, for example, you don't hear loud music, power tools, engines, or other very noisy things.

If your hearing loss makes it hard to hear loud noises, this hearing aid may not be your best choice without help from a professional. If this hearing aid does not help you enough, ask for help from a hearing healthcare professional.

WARNING: When to see a doctor
If you have any of the problems listed below, please see a doctor, preferably an ear-nose-throat doctor (an ENT).

- Your ear has a birth defect or an unusual shape. Your ear was injured or deformed in an accident.
- You saw blood, pus, or fluid coming out of your ear in the past 6 months.
- Your ear feels painful or uncomfortable.
- You have a lot of ear wax, or you think something could be in your ear.
- You get really dizzy or have a feeling of spinning or swaying (called vertigo).
- Your hearing changed suddenly in the past 6 months.
- Your hearing changes: it gets worse then gets better again.
- You have worse hearing in one ear.
- You hear ringing or buzzing in only one ear.

HONG KONG RETELESS CO., LIMITED
ADD: FLAT/RM 1019B 10/F LIVEN HOUSE NO.61-63
KING YIP STREET KWUN TONG KL

FC **CE** **UK** **CA**

MADE IN CHINA