


Labels and S/N Medimi™ Pro CDMA C200W

1. The Product and FCC label


1.1 Product label on the product



Color: **White**
Material: **MT300 Series, UL approved Polyester**
Size: **40 x 20-22 mm**
Print: **Black/white/gray**
Protection: **The label is placed under a transparent plastic parts and the user cannot get in direct contact with the label.**

W: **Production year and week: yyww (1619 = year 2016 week 19)**

1.2 Placement of the label on the product



Label
protected
under the
transparent
plastic parts

The label is protected under the transparent plastic tablet frame.
So the user will not be to touch the label.


Document No: **Labels and SN Spec Medimi CDMA**

Prepared and Approved by: 2016-02-16 / HH

Updated: 2016-07-14 / Håkan Harvigsson

2. The carton Product and FCC label

2.1 Product label on the carton



Material: **White, plastic**
Size: **125 x 52 mm (L x H)**
Print: **Black/white/gray**

W: **Production year and week: yyww (1619 = year 2016 week 19)**

2.2 Placement of the label on the product



Document No: **Labels and SN Spec Medimi CDMA**

Prepared and Approved by: 2016-02-16 / HH

Updated: 2016-07-14 / Håkan Harvigsson

3. S/N

3.1 Creating the S/N


S/N 1ddddywwnnnn (12 digits)

1 = Medimi
ddd = Producer ID no or batch/production line no (producers own choice)
yyww = Production year and week, 2016 w.24 -> yy=16 ww=24
nnnn = running production no, start with 0001 every new production week

The S/N are also added in the memory of the board, it's done during downloading the software into the board.


3.2 The S/N label on product

Placement of the label (small) on the product



Color: White
Material: MT300 Series, UL approved Polyester
Size: 45 x11 mm (L x H)
Print: Black and white
Barcode: 128

3.3 Placement of the product S/N label



Placed on the battery pack.

Document No: **Labels and SN Spec Medimi CDMA**

Prepared and Approved by: 2016-02-16 / HH


Updated: 2016-07-14 / Håkan Harvigsson

3.4 The S/N label on the carton



Material: **White, paper**
Size: **70 x 36 mm (L x H)**
Print: **Black/white/gray**
Barcode: **128**

3.5 Placement of the S/N label on the carton



Document No: **Labels and SN Spec Medimi CDMA**


Prepared and Approved by: 2016-02-16 / HH

Updated: 2016-07-14 / Håkan Harvigsson

4. MEID


4.1 The MEID label on product

Placement of the label (small) on the product



Color: **White**
Material: **MT300 Series, UL approved Polyester**
Size: **45 x11 mm (L x H)**
Print: **Black and white**
Barcode: **128**

4.2 Placement of the MEID label




Placed on the battery pack.

Document No: **Labels and SN Spec Medimi CDMA**

Prepared and Approved by: 2016-02-16 / HH

Updated: 2016-07-14 / Håkan Harvigsson

4.3 The MEID label on the carton



Material: **White, paper**
Size: **70 x 36 mm (LxW)**
Print: **Black/white/gray**
Barcode: **128**

4.4 Placement of the MEID label on carton

