


Label & Location

Back Side



Label Location

MEDICAL-GENERAL MEDICAL EQUIPMENT AS TO ELECTRICAL SHOCK,
FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH
ANSI/AAMI ES 60601-1(2005)+AMD(2012) and CAN/CSA-C22.2

No. 60601-1 (2008) + (2014)

Product Name: Tablet PC

Model : CyberMed Rx

Power : 19Vdc $\approx$ 3.78A

FCC ID: 2AHZW-CYBRXX10
IC: 24472-CYBRXX10



This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

USE EDAC POWER ELECTRONICS CO., LTD

Adapter model EM10681V, AC Input 100-240 Vac, 2.0-1.0 A, 50-60 Hz

Tested To comply With FCC class B Standards

www.cybernet.us

RoHS/WEEE



Product Name: Tablet PC

Power : 19Vdc $\approx$ 3.42A

Model : Rugged X10

USE FSP GROUP INC.

Adapter Model : FPS065-REBN2 , AC Input 100-240 Vac, 1.5A, 50/60Hz

Tested To comply With FCC class B Standards

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.

RoHS/WEEE



www.cybernet.us

FCC ID: 2AHZW-CYBRXX10
IC: 24472-CYBRXX10

